IKA

designed for scientists

VC 10 lite VC 10 pro

Betriebsanleitung	DE	5
Ursprungssprache		
Operating instructions	EN	29
Mode d'emploi	FR	53
Руководство пользователя	RU	77
使用说明	ZH	101

Breve guía	ES	124	Bezpečnostné pokyny	SK	141
Brevi istruzioni	IT	125	Ohutusjuhised	ET	142
Veiligheidsaanwijzingen	NL	127	Drošības norādījumi	LV	144
Säkerhetsanvisningar	SV	128	Saugos nurodymai	LT	145
Sikkerhedshenvisninger	DA	130	Инструкции за безопасност	BG	146
Sikkerhetsinformasjon	NO	131	Indicații de siguranță	RO	148
Turvallisuusohjeet	FI	132	Υποδείξεις ασφαλείας	EL	150
Instruções de segurança	PT	134	Güvenlik talimatları	TR	151
Wskazówki dotyczące bezpieczeństwa	PL	135	안전 지침	KO	153
Bezpečnostní pokyny a informace	CS	137	安全上の注意	JA	154
Biztonsági utasítások	HU	138	คำแนะนำเกี่ยวกับความปลอดภัย	TH	155
Varnostna navodila	SL	139	Petunjuk keselamatan	ID	157

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
Fig. 1

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Contents

	Page
Device setup	2
Declaration of conformity	29
Explication of warning symbols	29
Warranty	29
Safety instructions	30
Unpacking	31
Intended use	31
Useful information	31
Setting up	33
Commissioning	36
Interfaces and outputs	46
Maintenance and cleaning	47
Error codes	48
Accessories	49
Product contact parts	49
Technical data	50

Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/35/EG, 2014/30/EG and 2011/65/EG and conforms with the standards or standardized documents EN 61010-1, EN 60529, EN 61326-1 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Explication of warning symbols

 DANGER	Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
 WARNING	Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.
 CAUTION	Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
 NOTICE	Indicates practices which, if not avoided, can result in equipment damage.

Warranty

In accordance with warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Safety instructions

For your protection

NOTICE

Read the operating instructions in full before starting up and follow the safety instructions.

- Keep the operating instructions in a place where it can be accessed by everyone
- Ensure that only trained staff work with the instrument.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

DANGER

Always wear personal protective equipment in accordance with the hazard class of the media being worked with, otherwise there is a risk of splashing liquids, release of toxic or inflammable gases!

- The device and device parts must not be used on humans or animals.
- Do not expose parts of the human body or animals to the vacuum.
- Position the device as specified in the "Commissioning" section and connect the connection cables and interfaces as described.
- Never work with an incorrectly connected or defective device.
- Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapours, dusts or biological and microbiological materials can be hazardous to user. Make sure that all connections are tight and leak-free when working with such media.
- Prevent the release of the above mentioned substances. The use of suitable extraction systems e.g. fume cupboards is recommended in the place of installation. Take protective measures for personnel and the environment.
- Due to the residual leakage rate of the device, this can result in the release of medium.
- The VC 10 lite/pro vacuum controller is not designed to be installed in potentially explosive atmospheres
- The VC 10 lite/pro vacuum controller is not designed for use with self inflammable substances, substances which are inflammable without air or explosive substances.
- Take suitable precautions to prevent any formation of explosive mixtures. Use inert gas for venting and/or dilution, if necessary.
- Take into consideration interactions or possible chemical or physical reactions when working with media under reduced pressure and elevated temperature.
- Never work with the device under water or underground.
- The device must only be operated with supervision.
- Electrostatic events can occur between the medium and device and lead to a direct hazard.
- Safe operation is only guaranteed with the accessories described in the "Accessories" section.
- The instrument can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- Do not use the device if the AC adapter is defective.
- Follow the operating instructions of the additional equipment carefully (e.g. a rotary evaporator, vacuum pump) with which the VC 10 lite/pro vacuum controller is operated.

- The VC 10 lite/pro vacuum controller may only be operated under the conditions described in the "Technical data" chapter.
- The pressure at the gas inlet and gas outlet must not exceed 1100 mbar. For pressure over 1100 mbar, the total allowable pressure is exceeded and the device no longer displays the pressure accurately.
- Only use flexible hose lines.
- Flexible elements can be compressed under vacuum.
- In the event of a power failure, the integrated venting valve automatically aerates the connected recipient.
- Be aware of your emergency measures in the event of a power failure and ensure that the system is put into a safe state.
- The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly, contrary to the manufacturer's specifications or if the instrument or the printed circuit board are modified by third parties.

For protection of the equipment

CAUTION

The voltage stated on the type plate must correspond to the mains-voltage.

- Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation, so long as no other connection is made at this point. This will prevent the infiltration of foreign objects, liquids etc.
- The device must only be operated with the original plug-in power supply unit supplied.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts (authorized staff).
- To ensure sufficient cooling of the VC 10 lite/pro vacuum controller, the ventilation slots on the housing must not be covered.
- Use only original spare parts for repair and maintenance in order to ensure the reliable operating condition of the device.
- Pay attention to water condensation inside and outside of the device. If the device is brought from a cold environment, allow the device to warm up.
- Never attach the vacuum controller above a heating bath.
- Ensure that no solids and/or liquids enter the VC 10 lite/pro vacuum controller via the suction power of the pump. This will result in damage to the pressure sensor and the valves.

WARNING

Check the compatibility of the substances used with the device with the materials listed under "Product contact parts"; see chapter "Product contact parts".



- Place the VC 10 lite/pro vacuum controller in an elevated position if possible, so that in the event of a fault, no condensate can collect on the pressure sensor.

Unpacking

Unpacking

- Please unpack the device carefully.
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder)

Delivery scope

Vacuum controller VC 10 lite/pro		Vacuum hose 1 m		Fig. 4
Universal plug-in power supply unit		Y connector		Fig. 5
USB A - micro B cable		Operating instructions		
		Warranty card		

Intended use

Use

Together with the accessories recommended by , the VC 10 lite/pro vacuum controller is suitable for the controlled evacuation of air (gas) from laboratory devices and equipment (e.g. rotary evaporator or laboratory reactors, and also for classical separation or filtration or drying tasks in the laboratory) in conjunction with a suitable vacuum source (e.g. MVP 10 basic / VACSTAR). The VACSTAR can be operated automatically with the VC 10 pro

vacuum controller.

Therefore, it can be used for automatic boiling point recognition, programmed pressure-time curves or programs from the solvent library, for example.

Mode of operation: Stand device

Area of use

Indoor environments similar to that a laboratory of research, teaching, trade or industry area.

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- If the device is operated with accessories that are not supplied or recommended by the manufacturer;
- If the device is operated improperly or contrary to the manufacturer's specifications;
- If the device or the printed circuit board are modified by the third parties.

WARNING

The device must not be used for:

- **vacuation of biological biotopes,**
- **evacuation of explosive, corrosive or gas volumes,**
- **throughput/use of fluids!**

Useful information

In addition to a high-resolution Al_2O_3 pressure sensor, the VC 10 lite/pro vacuum controller has an integrated vacuum valve and a venting valve. As a result, the VC 10 lite/pro vacuum controller can be used without any additional peripheral devices.

The VC 10 pro provides two basic modes of operation: Two-position control and analogue speed control.

The change between the two modes is performed automatically by plugging or unplugging the connector cable at the "VACSTAR" interface. Depending on the selection made, some submenu items may not be available (greyed out on the display).

With two-position control, parameters from the solvent library can be applied or target values specified in manual mode.

With analogue speed control, an automatic mode as well as the setting of the percentage pump speed can also be preset. Space requirement: To connect all peripheral devices you will need approximately WxDxH 150 x 200 x 155 mm³ space for the device with connected plugs.

Two-position control VC 10 pro (only)

A vacuum is made in the glassware with the help of a vacuum pump. The vacuum pump operates at constant speed which generally cannot be adjusted. Once the target value setting is reached, the integrated vacuum valve interrupts the suction line.

The device will fall slightly short of the target value setting due to the time difference in the millisecond range with the vacuum valve, the target value comparison and the switching of the vacuum valve. Natural leakage from the apparatus now causes another increase in the pressure value of the system, which in turn is detected by the integrated pressure sensor.

If the target value is exceeded, the vacuum valve is reopened and the vacuum pump in operation lowers the pressure. The pressure curve shown on the display fluctuates around the set target value. The pressure difference between activating/deactivating the valve can be set using the vacuum hysteresis value.

The quality of the control depends on the suction power of the pump (rotation speed), the tightness of the system and signal delay time of the electronics used.

Analogue speed control

Imprecision of the two-position control is avoided by using the analogue speed control.

With this kind of control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value. Once the target value is reached, the pump only operates according to the leakage rate of the system. This type of control offers quieter operation, more exact vacuum control and a long service life of movable pump parts.

Automatic boiling point recognition is possible with this type of control, i.e. the system reaches and keeps the solvent at its boiling point in automatic mode.


Analogue speed control with the VACSTAR:

The VACSTAR is automatically detected as soon as the connection cable is plugged into both devices. The analogue speed control is enabled, the two-position control is disabled and cannot be selected.

Menu items such as automatic boiling point recognition, programmed pressure-time curves or programs from the solvent library can be used.


Refer to the operating instructions for the VACSTAR vacuum pump.

Schematic view of speed-vacuum control

Automatic boiling point recognition is not possible when using two-position control.

Schematic view of a two-position vacuum control

*An external temperature sensor can be connected to the "PT1000" temperature sensor interface. When the menu view is enabled, the temperature value is shown in the display.

Time sequences can be controlled using a timer.

The operating modes A, B, C and D have preset access options for adjusting and saving target values.

Display settings can be changed using the "Display" menu item.

Access to the device can be restricted by setting up a password in the "Security" menu item.

Common settings such as language, pressure and temperature units can be set via the "Settings" menu item. The display background can also be changed from black to white as well as brightness settings adjusted.

The integrated solvent library calculates the theoretical boiling point at a specified steam temperature value, or vice versa.

The upper and lower limits or the switching points of the vacuum control can be set in the submenu items "Limit values" and "Hysteresis".

**only VC 10 pro*

Setting up

NOTICE

Observe the general information, always connect the recipients (load/vacuum vessel/glass cooler) at their highest points to the suction line. This will prevent the risk of fluid entering into the vacuum controller or pump.

Install a separator (e.g. Woulff bottle) in front of the intake mounting to protect against the ingress of liquid.

In the case of induced solvent vapours, a vapour condenser (accessory of pump) installed downstream of the pump helps to condense the vapours and prevents release into the atmosphere.

Liquid in the pump chambers deteriorates the pump characteristics.

Attach the VC 10 lite/pro vacuum controller to a stand (d=16 mm) close to the pump.

Secure the controller with the stand locking screw (N).

Fig. 5

To do so, use a stand with diameter 16 mm, for example, stand VC 10.400, VSS 1 vacuum safety set, the stand on VACSTAR vacuum pump or attach the controller to the RV 8 rotary evaporator, or at LR 1000 laboratory reactor.

Fig. 6

VC 10 lite/
pro at stand
VC 10.400

Fig. 7

VC 10 lite/pro to VSS 1 with MVP
10 basic

Fig. 8

VC 10 lite/pro to VACSTAR

Fig. 9

VC 10 lite/pro
to RV 8

Connecting the interfaces

Connect the vacuum hoses, pay attention to the following connections (also see **Fig. 1**):

J: Load

Hose connection d=10 mm to the system to be evacuated (recipient).

K: Pump

Hose connection d=10 mm to the vacuum source (e.g. pump).

L: Vent

Ventilation connection d=10 mm, also suitable for the connection of inert gas.

Secure the hoses with hose clamps, if necessary.

Lay the vacuum hoses according to your device configuration.

Now connect the connecting cables according to your system configuration:

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Connect the VC 10 lite/pro vacuum controller to a PC using the USB 2.0 cable (A/Micro-B).

Software updates to any existing equipment can be installed using the FUT software tool.

M: RS 232 Interface

You can connect the VC 10 lite/pro vacuum controller to a PC using an RS 232 interface cable. The pump can be operated in conjunction with other devices using labworldsoft® laboratory software. For more information see Chapter "Interfaces and outputs".

P: PT 1000 temperature sensor (VC 10 pro)

Connect the optional PT 1000.60, PT 1000.61 or PT 1000.70 temperature sensor (accessory). The measured temperature (e.g. tempering bath) is shown on the display.

The value displayed is for information purposes only and is not processed logically.

Note that the value displayed is dependent on the positioning of the sensor in the medium, its mixing, as well as the calibration of the probe and can deviate from other measured temperature values.

Incorrect or fluctuating measured values can be displayed by using a heating bath, e.g. by laying the sensor on the heating bath floor, instead of positioning in the medium.

Q: VACSTAR membrane vacuum pump (Mini DIN)

You can connect the vacuum controller and the VACSTAR vacuum pump to the analogue connection cable (accessory) for exact speed-controlled vacuum control. The vacuum controller detects the pump and switches to speed-vacuum control mode. The two-position control is disabled. The pump speed is regulated depending on the measured pressure.

R: Connector for 24 W AC adapter

Connect the universal AC adapter (included) to the socket on the rear. Observe the country-specific variants of the AC adapter!

Explanation of abbreviations used in the following schematic drawings and configurations:

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Separator (Woulff bottle)
- 3 Venting valve
- 4 Pressure sensor
- 5 Vacuum valve/ball valve
- 7 Vapour condenser
- 8 Connection cable (analogue)
- 9 Vacuum source (pump, in-house vacuum)
- 10 VC 10 lite/pro.300 Check valve

Two-position control

Vacuum pump/vacuum source with VC 10 pro vacuum controller.

Once the target value is reached, the suction line is closed by the vacuum valve integrated in the VC 10 pro vacuum controller. The pump speed can be changed by manually adjusting the speed setting.

For pumps which have the speed setting option, it is recommended to work at low speed.


NOTICE

When using the VACSTAR vacuum pump: The pump must be operated in mode "B" (restart after voltage drop).

The vacuum source can also be a central in-house vacuum network.

Illustration two-position control - vacuum pump/vacuum source with vacuum controller VC 10 pro

Vacuum source with two VC 10 pro vacuum controllers.


Once the target value is reached in a recipient, the suction line is closed by the vacuum valve integrated in the VC 10 pro vacuum controller.

A sufficiently stable vacuum source (in-house vacuum, pump) can serve multiple consumers.

For pumps which have the speed setting option, it is recommended to work at low speed.


A second consumer can be controlled independent of the first using a second vacuum controller.

Illustration two-position control - vacuum source with two vacuum controllers VC 10 pro


The following working screen appears automatically in the display.
Display of two-position control, e.g. with MVP 10 basic.

Hysteresis display for two-position control

Display speed control: VACSTAR to analogue connection

No hysteresis display for speed control

Explanation of symbols on the working screen

The symbols displayed change depending on the status and settings of the vacuum controller:

Lock key

This symbol means that the function of the key and of the rotating/pressing knob for controlling the vacuum controller are disabled. The symbol no longer appears if the functions are enabled once again by pressing the lock key a second time (min. 1 sec.).

Operating mode

This symbol indicates the operating mode currently selected (A, B, C, D).

USB

This symbol means the vacuum controller is communicating via a USB cable. The symbol no longer appears if no USB cable is being used for communicating with the station.

Temperature sensor (VC 10 pro)

This symbol appears when the temperature display is active in the display.

Timer

This symbol appears when the timer display is active in the display.

PC

PC control

This symbol means that the vacuum controller is connected to a computer and the device is being controlled from the computer.

PR

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This symbol means that the vacuum controller is connected to a computer and the device is being controlled from the selected program.

Continuous mode

This symbol indicates that the vacuum controller is in continuous mode.

Vent

This symbol indicates that venting is enabled.

Vacuum

This symbol indicates that the current pressure is less than the existing atmospheric pressure.

Detection

This symbol is displayed during the automatic boiling point search, the system pressure is reduced further.

Boiling

This symbol indicates that a boiling point has been found automatically, the system pressure is kept constant.

Refrigerated circulator

This symbol indicates that refrigerated circulator is operating.

Navigation menu

Fig. 9


- ☞ Press the **"Menu"** key (C).
- ☞ Select the menu by turning the rotating/pressing knob (D) to the right or left to select the desired menu or sub-menu, which can then be selected by pressing the rotating/pressing knob (D).
- ☞ Press or turn the rotating/pressing knob (D) again to select the desired menu option and edit the values or settings, or activate/deactivate a function.
- ☞ Turn the rotating/pressing knob (D) to **"OK"** and press the **"Back"** (B) or **"Menu"** key (C) to end the procedure and return to the previous menu or working screen.
- ☞ Navigating in Menu "Graph": Press the **"Back"** (B) at main screen to navigate to graph menu, showing real-time pressure change on graph.

NOTICE

The menu option activated is highlighted yellow in the display. Menus which are grayed out are inactive.

Navigation menu: ▼ Press the **"Menu"** key (C) and turn the rotating/pressing knob (D)
 ▲ Press the **"Back"** key (B) or the **"Menu"** (C)

NOTICE


If you press the **"Menu"** key (C), the system skips directly back to the working screen.
 If you press the **"Back"** key (B), the system skips back to the previous display.

Menu structure


Menu	Sub menu	Option	Action	Factory setting	
Pump	Modes	Automatic (only VC 10 pro).....		Deactivated	
		Manual.....		Activated	
		Pump %.....		Deactivated	
	Refrigerated circulator	Program.....		Deactivated	
		Start After.....		-	
		Temperature.....		-	
		Pump speed.....		-	
		Cool now.....		-	
		Stop now.....		-	
	Cleaning	Start after.....		05:00 [mm:ss]	
		Duration.....		06:00 [mm:ss]	
		Pump speed.....		20 %	
		Clean now.....		Deactivated	
	Limits	Manual	Minimum.....		1 mbar
			Maximum.....		1100 mbar
		Pump %	Minimum.....		0 %
			Maximum.....		100 %
	Hysteresis (VC 10 pro)	Automatic.....		Activated	
	Solvent library (VC 10 pro)	Manual.....		Deactivated	
		Acetonitril.....		-	
		N-Amylalkohol.....		-	
		n-Pentalon.....		-	
		Benzol.....		-	
	Actions after stop	Open venting valve.....		Activated	
		Cleaning.....		Deactivated	
		Stop refrigerated circulator.....		Deactivated	
	Actions after start	Turn on refrigerated circulator.....		Deactivated	
	Leakage Test	Leakage rate.....		3 mbar	
		Tolerance.....		1mbar	
		Test point.....		500 mbar	
Leakage test duration.....			1 minute		
Temperature (VC 10 pro)	Sensor temperature.....		-		
	Display.....		Deactivated		
	Calibration.....		20.0 °C		
Timer	Timer function	Time.....		00:00:01 [hh:mm:ss]	
		Beep after stop.....		Activated	
		Activate.....		Deactivated	
	Display.....		Deactivated		
Operating mode	A.....		Deactivated		
	B.....		Activated		
	C.....		Deactivated		
	D.....		Deactivated		
Display	Timer.....		Deactivated		
	Temperature (VC 10 pro).....		Deactivated		
Programs	Program 01.....		-		
	...		-		
	Program 10	Select.....		-	
		Edit.....	Edit, insert, delete	-	
		Delete.....		-	
	Last measurement	Save as.....		-	
	Edit.....	Edit, insert, delete	-		
	Delete.....		-		
Safety	Password.....		0 0 0		
Service	Valves	Venting valve.....	Close, open	Open	
		Vacuum valve (VC 10 pro).....	Close, open	Close	
	Pump.....		Deactivated		
Refrigerated circulator.....		Deactivated			
Settings	Languages	English.....		Activated	
		German.....		Deactivated	
	Units	Pressure	mbar.....		Activated
			hPa.....		Deactivated
			mmHg.....		Deactivated
		Temperature (VC 10 pro)	Torr.....		Deactivated
			°C.....		Activated
			°F (not available for Japan).....		Deactivated
	Display	Background	Black.....		Activated
			White.....		Deactivated
	Brightness.....		100 %		
	Firmware update info.....		Activated		
	Sound	Volume.....		10 %	
		Key tone.....		Deactivated	
	Factory setting	Menu values.....	Reset to factory setting	-	
Programs.....			-		
All.....			-		
Communication	Device name.....		VC 10 lite/pro		
Information	Firmware update info.....		-		
	Display version.....		-		
	Logic version.....		-		
	Operating mode.....		B		
	Pressure max.....		1100 mbar		
	Pressure min.....		1 mbar		
	Pump % max.....		100 %		
	Pump % min.....		0 %		

Menu details

Pump

Modes

Automatic (VC 10 pro):


In the "Modes" menu, the user can enable automatic boiling point recognition by selecting the "Automatic" menu item. No other parameters must be set.

The boiling point is detected automatically. For rotary evaporators with heating bath, it must be ensured that the tempering medium and solvents have a constant temperature (e.g. 60 °C).

Manual:

In the "Manual" menu item, the target value can be specified (e.g. in "mbar"). The system is evacuated until the target value.

Refrigerated circulator


Start after...:

Specifies a delay before the refrigerated circulator start operating.

Temperature:

Determines the temperature, that the refrigerated circulator is operated.

Cleaning

The cleaning of the pump ensures that solvent residues are removed in the pump after the measurement. For this purpose, the pump is operated at a low level while the system is vented. The air flowing through the pump ensures that the solvent residues are discharged from the pump. If the cleaning option is selected, this starts once the measurement has finished (see "Actions after stop").

Pump %

By selecting the "Pump %" menu item, the pump can be operated continuously with a running performance of between 100 % and 1 %.

Program

Under the "Programs" menu, 10 user-defined pressure-time profiles can be created. The last measurement taken can also be viewed in this menu and saved as a program.

The target value and the running performance can only be entered within the set limits (see menu item "Limits").

Changing to the working screen and pressing the rotating/pressing knob starts the process. Pressing the knob again stops the process.

NOTICE

If no speed-controlled pump with analogue connection is connected to the "VACSTAR" interface, the system can only work in manual mode and in "two-position control" mode. The menu items "Auto" and "Pump %" are greyed out and cannot be activated.

Pump speed:

Determines the speed in rpm, that the refrigerated circulator is operated

Cool now:

If you want to perform a cooling immediately, you can select it with this setting.

Stop now:

If you want to stop the cooling immediately, you can select it with this setting.

NOTICE

The "Cleaning" menu item increases the service life of the pump, since solvent residues do not remain inside the pump.

You can also clean the pump without using the corresponding menu item. To do this, remove the hoses and operate the pump at the end of a working cycle (test, test series, end of the day ...) at idling speed.

Start after...:

Specifies after which minimum duration of a measurement the cleaning should take place. The default time is 5 minutes. This means that a measurement must be run for at least 5 minutes, so that cleaning is performed at the end of the measurement.

Duration:

Specifies how long the cleaning process should be performed.


Pump speed:

Determines at what power (in %) the pump is operated during cleaning. If the power is set too high this can lead to a pressure drop in the system.

Clean now:


If you want to perform a cleaning, you can select it with this setting.

Limits

In this menu, limits for the target value can be set.

Hysteresis (VC 10 pro)

With the "Hysteresis" option, the switching frequency and control accuracy can be influenced when in a steady state. The hysteresis value determines the upper and lower limits of the target value and therefore the closing and opening of the pump and valve (only for two-position control).

Automatic: Hysteresis is always 10 % of the actual pressure.

Manual: Manual value specification of hysteresis.

NOTICE


Limits which are too narrow shorten the service life of the devices.

Solvent Library (VC 10 pro)


In the "solvent library", the most common solvents are assigned a function to determine the steam temperature at the set pressure and vice versa. This provides assistance when setting these parameters in manual mode, for example, evaporation processes with a rotary evaporator.

Actions after stop


This menu allows the user to define actions that should be performed when a test has ended.

Actions after start

This menu allows the user to define actions that should be performed when a test has started.

Leakage test

Leakage test: User to confirm the vacuum connection. Press ok if the connection is good.

Set Value: User to input the required leakage rate with tolerance. Default leakage rate is 3mbar with tolerance 1mbar.


Test point: User to input the pressure which the leakage test will

be executed at. Default value is 500 mbar.

Duration: User to input duration of the leakage test. Default value is 1 minute.

In progress: The pump will start running to reach the test point. (Graph below: Phase 1 and 2)

The leakage measurement will start once the pressure has reached the interested test point. Meanwhile, the timer will start counting down from the pre-set duration. (Graph below: Phase 3)

Temperature (VC 10 pro)

Temperature sensor

In the "Temperature" menu the user can specify that the temperature sensor is displayed on the display/working screen. A tick (✓) shows that the option is activated.

The precondition for this is that a temperature sensor is connected to the vacuum controller. If no temperature sensor is connected, or an error is present or the temperature rises above 350 °C, three dashes are shown as the temperature value.


Please comply with the temperature measurement range of the external temperature sensor specified in the "Technical data" section. The temperature display can be in either °C or °F (see the "Settings" menu).


Calibration

In order to calibrate the external temperature sensor, it must first be connected to the vacuum controller. The temperature sensor is then immersed in a vessel with water, for example. The water temperature is determined using a second, already calibrated temperature measuring device. The measured temperature must then be entered in the menu under "Temperature calibration". Once the entry has been confirmed with "OK", the external temperature sensor is calibrated.

Timer

Timer functions

In the "Timer" menu the user can specify that the timer is displayed on the display/working screen. A tick (✓) shows that the

option is activated. This setting allows the user to specify the actual time for the evacuation procedure.

A default time can also be set for the timer. This setting allows the user to start evacuation for a standard time. The device can emit a signal (beep) after the set target time has elapsed.

The user can stop the evacuation procedure before expiry of the set time. In this case the countdown of the timer is interrupted.

Operating mode

Operating mode A:

In this operating mode, the set target value is not saved when the current run comes to an end or the device is switched off.

Operating mode B:

In this operating mode, the set target value is saved when the

current run comes to an end or the device is switched off, and the value can be changed.

Operating mode C:


In this operating mode, the set target value is saved when the current run comes to an end or the device is switched off, and the value cannot be changed.

Operating mode D:

In this operating mode, the target value is applied if it does not change after 3 seconds, or when the yellow background is displayed behind the target value. In operating modes A and B, the target value is applied immediately when the target value is changed by turning the knob.

In operating mode C, the target value cannot be changed.

Display

In the "Display" menu the user can specify what information will be displayed on the main screen.

Graph

In the "Graph" menu, user can observe the real-time pressure change on the shown graph.

Programs

Select:

To select a program, use the rotating/pressing knob and press "Select". Only programs that have at least one program segment can be selected.

When a program has been successfully selected, this is indicated by a tick mark (✓). The vacuum controller is then in program mode. In the main screen, the program view can then be selected by turning the knob to the right. Turning the knob to the left returns the user to the main screen.

Start:

To start a selected program, the push knob must be pressed in the main screen.

Edit:

Edit the selected program parameters. Start to edit the selected program parameters by pressing on menu option "Edit" with rotating/pressing knob. The user can edit, delete or insert one selected program segment in the program.

Delete:

Deletes the selected program. If a selected program is deleted by pressing on menu option "Delete" with rotating/pressing knob, all the program parameters will be emptied. The tick (✓) disappears.

Details for editing the program

Program 1		
No.	Pressure	hh:mm:ss
01	600 mbar	00:00:10
02	550 mbar	00:00:20
03	500 mbar	00:00:30
04	450 mbar	00:00:40
05	400 mbar	00:00:50
06	350 mbar	00:01:00
07	300 mbar	00:01:10
08	250 mbar	00:01:20
edit insert delete		

When a program is being edited, following screen appears.

In this program, user can define up to 10 segments. The selected segment is highlighted. Then, the user can edit, delete or insert a segment in this program. The program is save automatically.




Edit:

If the background of a selected value is yellow, the user can change the setting of the pressure value or the time value.

Delete:

When a highlighted segment is deleted, the subsequent segments move upward, closing the gap created by the deleted segment.


Sample for editing the program

Example of saving the last measurement


Safety


In the "Password" menu, the user can protect the vacuum controller settings using a password.
(factory setting: 000)

Service

In the "Service" menu, the valves, the pump or the refrigerated circulator can be operated individually and also checked to ensure that they are operating correctly. If no VACSTAR is connected, the menu item "Pump" is greyed out and not active.

Settings

Languages:

The "Languages" option allows the user to select the desired language by turning and pressing the rotating/pressing knob (D). A tick (✓) indicates the language that is set for the system.

Units:

The "Units" option allows the user to select the desired unit for displaying the temperature and pressure by turning and pressing the rotating/pressing knob. For the temperature, the choice is between "°C" or "°F". For the pressure, the choice is "mbar", "hPa", "mmHg" or "Torr". A tick (✓) indicates the unit that is set for the system.

Display:

The "Display" option allows the user to change the background color and brightness of the working screen.

Sound:

The "Sound" option allows the user to activate/deactivate the key-press sound and to set the volume.

Factory settings:

Select the "Factory settings" option by turning and pressing the rotating/pressing knob. The system requests confirmation to re-create the factory settings. You can choose if you want to reset only the "Menu values", or only the "Programs". Alternatively, you can reset all by selecting "All". Pressing the "OK" key resets the system settings to the original standard values set at dispatch from the factory (see "Menu structure" illustration).

Communication:

The "Device name" option allows the user to input a device name. The device name appears on the opening screen when the device is switched on.

This can be useful when using multiple devices with different settings. This also facilitates the identification of each device when communicating with a PC via USB, RS232 and Bluetooth, for example.

Information:

The "Information" option offers the user an overview of the most important system settings of the vacuum controller VC 10 lite/pro.

Interfaces and outputs

The device can be operated in "Remote" mode via an RS 232 or USB interface using the laboratory software *labworldsoft*[®].

The RS 232 interface at the back of the device is fitted with a 9-pole SUB-D port which can be connected to a PC. The pins are assigned serial signals.

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Use the USB interface in conjunction with *labworldsoft*[®] for operation in "Remote" mode and also to update the firmware.

Installation

First, download the latest driver for devices with USB in-terface from <http://www.com/lws/download/usb-driver.zip> and install the driver by running the setup file. Then connect the device through the USB data cable to the PC.

Serial interface RS 232 (V 24)

Configuration:

- The functions of the interface connections between the device and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- Transmission procedure: asynchronous character transmission in start-stop mode.
- Type of transmission: full duplex.
- Character format: character representation in accordance with data format in DIN 66 022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- Transmission speed: 9600 bit/s
- Data flow control: none
- Access procedure: data transfer from the device to the computer takes place only at the computer's request

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- Commands are generally sent from the computer (Master) to the stirrer machine (Slave).
- The stirrer machine sends only at the computer's request. Even fault indications cannot be sent spontaneously from the device to the computer (automation system).
- Commands are transmitted in capital letters.
- Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- Each individual command (incl. parameters and data) and each response are terminated with Blank CR LF (Code: hex 0x20 hex 0x0d hex 0x20 hex 0x0A) and have a maximum length of 50 characters.
- The decimal separator in a number is a dot (Code: hex 0x2E).

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The NAMUR commands and the additional specific commands serve only as low level commands for communication between the stirrer machine and the PC. With a suitable terminal or communications programme these commands can be transmitted directly to the stirrer equipment.

The software package, *labworldsoft*[®], provides a convenient tool for controlling the equipment and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

The following table summarises the (NAMUR) commands understood by the control equipment.

NAMUR Commands	Function
IN_NAME	Read device name
IN_PV_3	Read PT 1000 value
IN_PV_66	Read current pressure value
IN_SP_66	Read pressure target value
IN_SP_70	Read hysteresis target value
IN_MODE_66	Read current evacuation mode
IN_ERROR	Read error status
OUT_SP_66	Set pressure target value
OUT_SP_70	Set hysteresis value
OUT_MODE_66	Set evacuation mode
START_66	Start evacuation
STOP_66	Stop evacuation
RESET	Switch to normal operating mode

PC 1.1 Cable (Device to PC)

Required for connecting the 9-pin socket to a PC.


Fig. 10

Connection VC 10 lite/pro - PC





Fig. 11

Connection Refrigerated circulator

PC1.3 Cable (Device - Refrigerated circulator)

Connection (Device - Refrigerated circulator)

Connection cable

This cable is required to connect the VC 10 lite/pro to the VACSTAR (speed mode).

Fig. 12

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

- For cleaning disconnect the mains plug!
- Use only cleaning agents which have been approved by to clean devices. These are water (with tenside) and isopropanol.
- Wear protective gloves during cleaning the devices.
- Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the device when cleaning.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with that this method does not destroy the device.

Spare parts order

When ordering spare parts, please give:

- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.com.
- software version

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For repair, please request the "Safety Declaration (Decontamination Certificate)" from or use the download printout of it from website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Error codes

Error:

The fault is shown by an error message in the display as following if the error occurs, e.g. Error 4.
Proceed as follows in such cases:

- Switch the device switch off,
 - carry out corrective measures,
 - switch the device switch on again.
- ⇒ Measurement is stopped.

Warning:

- Measurement is stopped.
- Warning is displayed.
- Warning can be confirmed.

Information:

- Measurement continues.
- Information is displayed.
- Information can be confirmed.

Description	Detection	Reason	Action
Information No Pressure Change	Pressure gradient analysis does not show any deviation.	Pump is not running. Vacuum hose isn't connected. Failure in hose system. Recipient isn't closed.	Check pump control cable connection. Check pump power supply. Switch on pump. the position. Check correct vacuum hose connection (IN/ OUT/VENT). Close recipient.
Information System Not Tight	Pressure gradient analysis shows pressure deviation, but set value cannot be reached.	Vacuum hose connection isn't tight. Recipient isn't tight. Vacuum pump power isn't sufficient. Venting valve leakage. Boiling point of solvent is reached before.	Check vacuum hose connection. Check recipient. Check technical data of the pump. Call service department. Check set value.
Information Venting Error	Pressure gradient analysis does not show any deviation after pressing the "Venting" button.	Recipient wasn't evacuated at all. Failure in hose system. Venting valve leakage. Front foil button is damaged.	Evacuated recipient. Check correct vacuum hose connection (IN/ OUT/VENT). Call service department.
Pump not connected	Connection to Vacstar(speed mode) not present.	Damaged or no cable connected	Check connection cable is present and well connected.
Warning Boiling Point Error	Boiling point cannot be detected. Pressure gradient analysis does not show boiling point.	Heating bath does not heat up. Solvent with extreme low boiling point is used.	Check heating bath. Solvent distillation only manually.
Warning External Temperature Is Too High	Measured temperature is too high. External temperature is higher than maximum value	Medium temperature is too high. External sensor is broken.	Check the medium temperature. Call service department.
Warning External Temperature Is Too Low	Measured temperature is too low. External temperature is lower than minimum value.	Medium temperature is too low. External sensor is short-connected.	Check the medium temperature. Call service department.
Warning Control Pump State Error	Pump is disconnected during analog speed control mode. Pump is connected during 2-point control mode.	Pump is disconnected during analog speed control mode. Pump is connected during 2-point control mode	Connect the pump. Disconnect the analog pump.
Warning RS232 Communication Error USB Communication Error	Watchdog time elapsed.	labworldsoft® uses wrong device for VC 10 lite/pro. Wrong tools are used to communicate with VC 10 lite/pro. Unstable connection. Watchdog was set in device or via PC command.	Check the labworldsoft® setting. Check the PC tools. Check the connection. Disable watchdog when not necessary. Increase sample rate.

Description	Detection	Reason	Action
Warning Temperature Calibration Error	Temperature calibration value is out of range.	Wrong simulator resistance is chosen. Set wrong temperature calibration value.	Check the simulator resistance. Call service department.
Warning Pressure Calibration Error	Pressure calibration value is out of range.	Input wrong vacuum actual value when calibration. Set wrong vacuum calibration value. Sensor is broken.	Calibrate again. Check the set calibration value. Call service department.
Error 3 Device Temperature Error.	Analysis of internal PCB temp. Sensor. Inside temperature reaches limit value.	Room temperature > 40 °C. PCB or valve is broken.	Switch off device for cool down. Call service department.
Error 9 Logic Storage Error. Display Storage Error.	Read or write internal memory error.	Read or write internal memory error.	Call service department.
Error 68 Pressure Out Of Range.	Pressure is out of range. Pressure is bigger than room pressure.	Pressure in the recipient is too high. Sensor is broken.	Check the air flow of pump and vent recipient. Call service department.
Error 69 Pressure Sensor Error.	Pressure sensor output is too low.	Sensor is not connected. Sensor is broken.	Check the sensor connection. Call service department.
Error 71 Internal Communication Error	Internal Watchdog elapsed.	Communication lost between logic board and display board.	Restart device (unplug mains cable, plug in mains cable). Check communication between logic board and display board. Call service department.

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- contact the service department,
- send the instrument for repair, including a short description of the fault.

Accessories

Accessories see www.ika.com.

Product contact parts

Designation	Material
Connection spigots	PP
Distributor	PPS
Pressure sensor	FPM/AL2O3
Vacuum valve	PEEK; EPDM; FKM
Venting valve	PEEK; EPDM; FKM

Technical data

	Unit	VC 10 lite	VC 10 pro
Connection diameter suction side	mm	8	
Connection diameter pressure side	mm	8	
Connection diameter venting side	mm	8	
Input pressure min.	mbar	1	
Input pressure max.	mbar	1050	
Boiling point detection		-	yes
Solvent library			yes
Two-point control		-	yes
Analog speed vacuum control			yes
Display			TFT
Pressure unit / scale			mbar, hPa, mmHg, Torr
Vacuum sensor			yes
Vacuum sensor type			ceramic Al ₂ O ₃
Pressure max. for pressure sensor	bar	1.6	
Measurement range (absolute) min.	mbar	1	
Measurement range (absolute) max.	mbar	1100	
Control range min.	mbar	1	
Control range max.	mbar	1100	
Resolution pressure	mbar	1	
Measurement uncertainty	mbar	1	
Medium temperature (gas) min.	°C	5	
Medium temperature (gas) max.	°C	40	
Vacuum valve		-	yes
Venting valve			yes
Connection for ext. temperature sensor		-	PT 1000
Temperature unit		-	°C/°F
Temperature measuring range min.	°C	-	-10
Temperature measuring range max.	°C	-	200
Temperature measurement resolution	K	-	1
Accuracy of temperature measurement	K	-	±1
Timer			yes
Time setting min.	s	1	
Time setting max.	min	6000	
Vacuum speed control interface			VACSTAR
Material in contact with medium			Al ₂ O ₃ , PTFE, FPM, PPS
Housing material			PBT
Fastening			stand/clamp
Fastening diameter	mm	16	

	Unit	VC 10 lite	VC 10 pro
Dimensions (B x H x T)	mm	95 x 150 x 110	
Weight	kg	1.5	
Permissible ambient temperature min.	°C	5	
Permissible ambient temperature max.	°C	40	
Permissible relative humidity	%	80	
Protection class according to DIN EN 60529		IP 20	
RS 232 interface		yes	
USB interface		yes	
Voltage	V	100 - 240	
Frequency	Hz	50/60	
Power input	W	24	
Power input standby	W	2	
DC Voltage	V	24	
Current consumption	mA	1000	

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Phone: +90 216 394 43 43

eMail: sales.turkey@ika.com

Discover and order the fascinating products of IKA online:
www.ika.com

IKAworldwide

IKAworldwide /// #lookattheblue

@IKAworldwide

IKA

designed for scientists

VACSTAR lite

ENGLISH


Fig. 1

	EU Declaration of conformity	6
	Explication of warning symbols	6
	Safety instructions.....	7
	Intended use.....	10
	Unpacking	11
	Useful information	12
	Operator panel and display	14
	Installation	15
	Operation	16
	Interfaces and outputs	18
	Maintenance and cleaning.....	20
	Accessories	20
	Material of parts in contact with medium.....	20
	Error codes	21
	Technical Data	22
	Warranty	23

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Explication of warning symbols

/// Warning symbols

-  **Danger!** Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
-  **Warning!** Indicates a hazardous situation, which, if not avoided, can result in death, serious injury.
-  **Caution!** Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
-  **Notice!** Indicates practices which, if not avoided, can result in equipment damage.
-  **Caution!** Indicates crushing risk to fingers / hands.

/// General Symbols

-  **A** — Position number
Indicates device components relevant to actions.
-  **Correct / result**
Indicates the correct execution or the result of an action step.
-  **Wrong**
Indicates the incorrect execution of an action step.
-  **Note**
Indicates steps of actions that require particular attention.

Safety instructions

/// General information

- › **Read the operating instructions completely before starting up and follow the safety instructions. If this device is used with another device, observe also its operating instructions.**
- › Keep the operating instructions in a place where it can be accessed by everyone.
- › Ensure that only trained staff work with the device.
- › Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- › The device must only be used in a technically perfect condition.

Notice!

- › Pay attention to the marked sites in **Fig. 1**.

/// Device design

Caution!

- › Ensure that your hands / fingers do not get crushed when setting down the device.

Notice!

- › Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- › Set up the device in accordance with the chapter "Installation" and connect the connection lines and interfaces as described.
- › Never work with a faulty or incorrectly connected device.
- › Set up the vacuum pump in appropriate and functioning extractor hoods, or set up an exhaust gas line to them. Ensure that the exhaust gas line cannot be kinked! The maximum permitted exhaust gas line length is 2 metres.
- › Prevent occurrence of explosive mixtures; if appropriate connect inert gas for ventilation and / or thinning.

/// Working with the device

Danger!

- › Do not use the device in explosive atmospheres, it is not EX-protected.
- › With substances capable of forming an explosive mixture, appropriate safety measures must be applied, e.g. working under a fume hood.
- › To avoid body injury and property damage, observe the relevant safety and accident prevention measures when processing hazardous materials.
- › The device and all parts of the device must not be used on people or animals.
- › Do not expose human or animal body parts to vacuum.
- › You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.
- › The device is not suitable for use with self-igniting materials, materials that are flammable without an air supply, or explosive materials.
- › Do not operate the vacuum pump while it is open. Otherwise there is a danger of severe injuries as a hand may be pulled into the moving parts.

Warning!

- › Inhaling or coming into contact with media such as poisonous liquids, gases, spray mists, vapours, dusts or biological materials can endanger the health of the user. Ensure that all connections are well sealed and free of leaks if you are working with such media.
- › Prevent release of the materials listed above. Take measures to protect staff and the environment.

Notice!

- › The vacuum pump must only be operated under the conditions escribed in the chapter "Technical data".
- › Pay attention to possible interactions or chemical or physical reactions when working with media at reduced pressure and increased temperature.
- › There can be electrostatic processes between the medium and the device which can lead to direct danger.
- › Some medium may be released due to the residual leakage rate of the device.
- › Only switch the pump on if the pump is standing vertically.
- › Connect the hose connections (INLET-OUTLET) and interfaces in accordance with the labelling on the device and the operating instructions.
- › Ensure that the temperature of the evacuated medium is below its ignition temperature. The pumping process (compression) increases the temperature of the medium additionally.
- › Please note that vapours containing solvent can be sucked into the pump.
- › Do not use the pump to create pressure.
- › Observe the permitted pressure at the inlet and outlet side; see the chapter "Technical data".
- › The gas flow must only be regulated / throttled in the suction-side line.
- › Use a solenoid valve or a check valve in the suction line if there are several load units.
- › When using an emission condenser ensure that the coolant is flowing freely.
- › Ensure that the substances used with the device are compatible with the materials of the device parts that come into contact with the product; see the chapter "Material of parts in contact with medium".
- › Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation, unless anything else is connected at this point. This will prevent the infiltration of foreign objects, liquids and other contaminants.
- › Unwanted liquids that have entered can be removed by extracting air during no-load operation.
- › To ensure sufficient cooling for the vacuum pump, the ventilation slits on the housing must not be covered.
- › Only use original replacement parts for maintenance to ensure that the device works reliably.
- › Watch out for water condensation inside and outside the device. Warm up the device first, if it has been brought in from a cold environment.
- › Never fasten the vacuum pump over a heating bath.
- › Ensure that no solids and/or liquids can enter the diaphragm of the vacuum pump through the suction line of the pump. This damages the diaphragm and other internal parts of the pump. This reduces delivery capacity and the final pressure can no longer be achieved. Deposits may be produced in the interior, reducing service life and causing leaks.

/// Accessories

- › Protect the device and accessories from bumps and impacts.
- › Check the device and accessories for damage before each use. Do not use damaged components.
- › Safe operation is guaranteed only with the use of original IKA accessories.
- › Always disconnect the plug before attaching accessories.
- › Observe the operating instructions of the accessories.
- › Elastic elements may be pressed together under vacuum.
- › Only use flexible hose lines.

/// Power supply / Switching off the device

- › The voltage stated on the type plate must correspond to the power voltage.
- › The device can only be disconnected from the power supply by pulling out the power plug or the connector plug.
- › The device must only be operated with the original power cord set.
- › The socket for the power cord must be easily accessible.
- › Socket must be earthed (protective ground contact).

/// Maintenance

- › The device must only be opened by trained specialists, even during repair. The device must be unplugged from the power supply before opening. Live parts inside the device may still be live for some time after unplugging from the power supply.

/// Disposal instructions

- › The device, accessories and packaging must be disposed of in accordance with local and national regulations.

Intended use

/// Use

⚠ Warning!

- › The device must not be used for:
 - evacuating biological habitats,
 - evacuating explosive, corrosive or similar gases,
 - throughput / use of liquids.
- › The device is not suitable for:
 - pumping up objects,
 - pumping liquids,
 - use in potentially explosive atmospheres,
 - underwater,
 - underground.

The laboratory diaphragm vacuum pump IKA VACSTAR lite is used together with other laboratory devices to create a vacuum.

For specific vacuum control, an additional vacuum controller IKA VC 10 lite / pro or IKA rotary evaporator RV 10 auto is required.

Standalone operation

The vacuum pump can be operated without a controller.

The pump can only run at a fix speed if it is not connected to the controller (e.g. VC 10 lite / pro or RV 10 auto). The pump will go into the energy-saving mode after running continuously for 1 hour. The final pressure is not controlled.

Operation with accessories

Along with the accessories recommended by IKA, the vacuum pump is suitable for controlled evacuation to a specified final pressure e.g. for operating an IKA rotary evaporator.

It can also be used for classical separation, filtration or drying in the laboratory.

With an appropriate vacuum controller (e.g. VC 10 lite / pro) or with a rotary evaporator (e.g. RV 10 auto), the VACSTAR lite can be operated automatically, for automatic boiling point detection, scheduled pressure-time curves or programs from the solvent library, for example.

Intended use: Tabletop device

/// Area of use

- › Indoor environments similar to that a laboratory of research, teaching, trade or industry area.
- › The safety of the user cannot be guaranteed:
 - if the device is operated with accessories that are not supplied or recommended by the manufacturer,
 - if the device is operated improperly or contrary to the manufacture's specifications,
 - if the device or the printed circuit board are modified by third parties.


Unpacking

/// Unpacking

- › Unpack the device carefully. Any damage should immediately be reported to the carrier (mail, rail or freight forwarding company).

/// Scope of delivery

1	VACSTAR lite	5	Silencer
2	Power cord set	6	User guide
3	Connecting cable	7	Warranty card
4	Vacuum hose, 1 m		

Useful information

⊗ Danger!

You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.

⊗ Caution!

Special care must be taken when working with gas mixtures that include condensable gases (e.g. solvents). If diaphragm pumps are used with such gases they condense in the pump's air flow. This leads to a build-up of pressure in the vacuum chamber, which in turn influences the effectiveness and service life of the diaphragms and valves.

If liquid or condensate enters the vacuum chambers, the specified minimum working pressure is no longer achieved.

To protect the inner valves and diaphragms against liquid, use an upstream condensate separator, such as a Woulff bottle.

It is generally recommended that you run the pump with no load at the end of the process (no-load operation, approx. 3 up to max. 5 minutes) to dispel any residues in the pump.

If applicable, vent the recipient.

To do this, remove the suction line during manual operation.

⚠ Notice!

Observe the operating instructions of the relevant devices.

Diaphragm vacuum pumps are oscillating displacement pumps. The diaphragms expand the volume of the vacuum chamber to produce an under pressure, which sucks air into the chamber. The air is then pressed out of the vacuum chamber by subsequent compression. Changes between the inlet and outlet are controlled by diaphragm valves.

The physical unit for pressure is 1 Pascal [Pa].

However people often still use the unit millibar [mbar] colloquially.

1 mbar = 100 Pa

1 bar = 10^5 Pa

1 Pa = 10^{-5} bar

1 Pa = 0.01 mbar

A distinction is made between different ranges in vacuum technology.

Rough vacuum: 10^5 to 10^2 Pa (1000 to 1 mbar)

Fine vacuum: 10^2 to 10^{-1} Pa (1 to 10^{-3} mbar)

High vacuum: 10^{-1} to 10^{-5} Pa (10^{-3} to 10^{-7} mbar)

Ultra high vacuum $< 10^{-5}$ Pa ($< 10^{-7}$ mbar)

The periodic movements of the diaphragms create a vacuum at the pump inlet. Diaphragm vacuum pumps work with a rough vacuum of up to 10^2 Pa. Lower pressures cannot be achieved with these pumps.

To create fine, high and ultra high vacuums, oil-sealed vacuum pumps, screw pumps, diffusion pumps or turbomolecular pumps must be used, for example. Diaphragm vacuum pumps are often used as pumps for the preliminary stage in these contexts.

Compressing and releasing the pressure on air creates a noise. The pump is quieter the less air is transported. This means that there is less operating noise during normal operation if the pressure is reduced.

A fitted exhaust hose with a sound absorber can also be used to reduce noise.


In automatic operation with controlled speed the pump reduces its speed when it reaches the target value. If the leakage rate in the entire system is low enough, the speed may even be set to "zero".

When combined with the vacuum controller VC 10 lite / pro or the rotary evaporator with integrated vacuum controller RV 10 auto, the pump can also be used for "two-position control" (valve-controlled) or "analogue speed control" of the vacuum.

Two-position control

When the target value is reached, the suction line is interrupted. Due to the small delay between the times when the pressure value is detected, when comparison is made with the target value, when a vacuum valve is switched off and when the suction line of the vacuum-producing system are switched off, it is easy to go below the set target value. If the pressure in the system goes up again due to a natural leakage rate, the suction line valve opens.


Schematic view of two-position vacuum control


Analogue speed control

With this control, the speed of the pump, and therefore its suction power, is reduced the closer the measured pressure reaches the target value. Once the target value is reached, the pump only operates according to the leakage rate of the system. This enables quiet functioning, precise vacuum control and a long service life for the moving pump parts. Exact automatic boiling point recognition is possible with this type of control, meaning that the system reaches and keeps the solvent at its boiling point in automatic mode.

Schematic view of speed vacuum control

Operator panel and display

A	Power switch ("O" « off, "I" « on)	G	USB Interface
B	Start / stop button	H	Hose connection for suction line d = 8 mm (INLET)
C	LED bar	I	Ventilation slit
D	Mounting thread for VC 10.600 Stand pole	J	Hose connection for pressure line d = 8 mm (OUTLET)
E	Connection for VSE 1	K	Power socket
F	Connection for vacuum controller VC 10 lite / pro or rotary evaporator RV 10 auto (Mini DIN)	L	Fuse

Installation

Notice!

- › Ensure that you follow the general rule, that the recipient (load / vacuum container/glass condenser) should always be connected to the suction line at its highest point. This reduces the risk of liquids entering the vacuum pump.
 - › Liquid in the pump chambers makes the pump properties worse.
 - › Connect the hose connections (INLET / OUTLET) and interfaces in accordance with the labelling on the pump and the operating instructions.
- › The connections are barbed hose connectors. Ensure that the hoses have hose clips if applicable.
 - › Lay the hoses in such a way that no condensate can flow into the pump.
 - › Ensure that the air can enter the ventilation slit unhindered.
 - › Install a condensate separator (e.g. a Woulff bottle) in the suction line before the intake spigot of the pump to protect against solvent coming in.
 - › If solvent vapours are sucked in, an emission condenser (accessory) after the pump helps to condense it, and prevents it being released into the atmosphere.

/// Connection of interfaces

Notice!

- › Observe the relevant connections (see "Operator panel and display").

OUTLET: Hose connection d = 8 mm

You can link this connection to the emission condenser (accessory) with a vacuum hose, or fit a sound absorber to the end.

Notice!

- › Put the end of the hose in the extractor hood!
- › Check there is a free outlet on the pressurised side!
- › Do not use a throttle on the pressurised side and do not close the outlet! Connect the exhaust line to this connection.

Connection of vacuum controller VC 10 lite / pro or rotary evaporator RV 10 auto via Mini DIN

You can connect the vacuum pump and the vacuum controller VC 10 lite / pro or the rotary evaporator RV 10 auto with the analogue connection cable for precise speed-controlled vacuum control.

The vacuum controller detects the pump and switches to speed-vacuum control mode. Two-position control is deactivated.

USB interface

Connect the vacuum pump to a PC with a USB A - USB C cable. Any device software updates can be loaded using the IKA Firmware-Update-Tool.

INLET: Hose connection for suction line d = 8 mm

Connect the suction line to this connection.

Link this connection to the recipient (rotary evaporator condenser, laboratory reactor etc.) with a vacuum hose.

Power supply cable connection

Check that the voltage information on the rating plate matches your mains supply. Connect the power supply cable to supply electricity.


⚠ Notice!

Observe the operating instructions for the devices.

Automatic setting of the operating mode “Speed regulation” as soon as the analogue connecting cable connects the pump with the vacuum controller or rotary evaporator.

When the target value is reached, the pump speed will then be regulated automatically to maintain at desired pressure. Once the target value is reached, the pump applies suction according to the leakage rate of the vacuum system. Automatic boiling point recognition is possible.

1	Recipient (load, e.g. rotary evaporator, reactor)
2	Vacuum separator (e.g. Woufff bottle)
3	Vent valve
4	Pressure sensor
5	Vacuum valve / ball valve
6	Emission condenser
7	Analogue connecting cable

Operation

/// Switching on


- › Switch the device on using the power switch. The LED bar lights up static green, indicating that the pump is in standby mode.

- › To start the pump, press the start / stop button. The pump starts to run and the LED bar pulsates in green.
- › You can stop the pump by pressing the start / stop button again. In this case, the pump returns to standby mode.
- › If the pump is operated with the vacuum controller VC 10 lite / pro or the rotary evaporator via the analog connection cable or is connected to labworldsoft® via the USB cable, the LED bar will light up static yellow.
- › During the first two minutes of operation, the pump speed is limited to a maximum of 900 rpm in order to reduce the mechanical stress on the pump and thus extend its service life.
- › If the pump is operated without a control system, the pump speed is automatically reduced to 900 rpm after one hour of continuous operation in order to save energy.

/// Switching off

- › Switch off the device using the power switch.
- › After switching off, the LED bar will remain illuminated in green or yellow for around 10 seconds. This indicates that the internal discharge process is underway to safely dissipate the energy stored internally.
- › This indicates that the device is currently discharging its internal power supply.

Note: Please wait until the discharge is complete before opening the device or performing any further steps.

Interfaces and outputs

The device can be operated by computer via an USB interface using the laboratory software labworldsoft®.

The device software can also be updated with a PC via the USB port.

Notice!

Please comply with the system requirements together with the operating instructions and help section included with the software.

/// USB interface:

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized.

/// USB device drivers:

Connect the IKA device through the USB data cable to the PC. The data communication is via a virtual COM port.

From Windows 10 and onwards the standard Windows USB driver is automatically loaded and a COM port number is assigned (find details in Windows Device Manager: "USB Serial Port (COMxx)"). If you have problems with USB communication, first ask your IT system administrator whether access to the USB interface is restricted for data security reasons.

/// Command syntax and format:

The following applies to the command set:

- › Commands are generally sent from the computer (Leader) to the device (Follower).
- › The device sends only at the computer's request. Even fault indications cannot be sent spontaneously from the device to the computer (automation system).
- › Commands are transmitted in capital letters.
- › Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- › Each individual command (incl. parameters and data) and each response are terminated with Blank CR LF (Code: hex 0x20 hex 0x0d hex 0x0A) and have a maximum length of 50 characters.
- › The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).

The NAMUR commands and the additional specific IKA commands commissioning serve only as low level commands for communication between the device and the PC. With a suitable terminal or communications program these commands can be transmitted directly to the device. The IKA software package, Labworldsoft®, provides a convenient tool for controlling device and collecting data under MS Windows, and includes graphical entry features, for motor speed ramps for example.

NAMUR Commands	Function
IN_NAME	Read the device name
IN_PV_4	Read the actual speed value
IN_SP_4	Read the set speed value
OUT_SP_4	Adjust the set speed value
START_4	Start the evacuation
STOP_4	Stop the evacuation
RESET	Switch to normal operating mode
OUT_NAME	Set device name
OUT_SP_42@m	Setting WD safety limit speed with set value echo
OUT_WD1@m	Watchdog mode 1
OUT_WD2@m	Watchdog mode 2
IN_ERR	Error Status

"Watchdog" functions, monitoring of the serial data flow:

If, once this function has been activated (see NAMUR commands), there is no retransmission of the command from the computer within the set time ("watchdog time"), the evacuation function is switched off in accordance with the set "watchdog" function or are changed to the set target values.

The data transmission may be interrupted by, for example, a crash in the operating system, a power failure in the PC or an issue with the connection cable between the computer and the device.

"Watchdog" – mode 1:

If event WD1 should occur, the evacuation function is switched off and PC1 will be returned when IN_ERR is sent. Set watchdog time to m (20 - 1500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.

"Watchdog" – mode 2:


If there is an interruption in data communications (longer than the set watchdog time), the speed target value is changed to the set WD safety speed limit. PC2 will be returned when IN_ERR is sent. The WD2 event can be reset with OUT_WD2@0 - this also stops the watchdog function.

Set watchdog time to m (20 - 1500) seconds, with watchdog time echo. This command launches the watchdog function and must be transmitted within the set watchdog time.

/// Connections between device and external devices


USB cable A – C:

This cable is required to connect USB port to a PC.

Connecting cable:

This cable is required to connect VC 10 lite / pro to RV 10 auto with the diaphragm vacuum pump VACSTAR lite.

Maintenance and cleaning

- › The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

/// Cleaning

- › For cleaning disconnect the mains plug!
- › Use only cleaning agents which have been approved by IKA to clean the devices: Water containing surfactant / isopropyl alcohol.
- › Wear protective gloves during cleaning the devices.
- › Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- › Do not allow moisture to get into the device when cleaning.
- › Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

/// Ordering spare parts

- › When ordering spare parts, please give:
 - device type.
 - serial number, see type plate.
 - position number and description of spare part, see www.ika.com.

/// Repairs

- › Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.
- › For repair, please request the “**Safety Declaration (Decontamination Certificate)**” from IKA or use the downloaded printout of it from IKA website at www.ika.com.
- › If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Accessories

- › For accessories see www.ika.com.

Material of parts in contact with medium

Connection spigot	PPS	Diaphragm valve	FFPM
Connection piece	PPS	Diaphragm	NBR / PTFE
Distributor	PPS	Clamp	PPS
Head piece	PPS		

Error codes

- › If an error occurs, the display light bar will be shown in red. The error code can be read via USB using IN_ERR command or labworldsoft®. Proceed as follows in such cases:
 - Switch off the device using the power switch. After switching off, the red LED bar will remain lit for about 10 seconds during the internal power supply discharge.
 - Important:** Always wait at least 10 seconds before attempting to restart the device. Restarting too soon may prevent a proper reset.
 - Carry out corrective measures.
 - Restart the device. Once the red LED bar has gone out, you can safely restart the device by switching it back on.

Error code | Causes | Effect | Solutions

E22273824

Causes	› Overheat
Effect	› The pump does not run
Solutions	› Switch off the pump, disconnect it from the mains, let it cool down.

E22240702

Causes	› Motor blocked / overload
Effect	› The pump does not run
Solutions	› Switch off the pump, disconnect it from the mains, let it cool down. › Ventilate the vacuum system to normal pressure.

E22243020

Causes	› Hall sensor defective
Effect	› The pump does not run
Solutions	› Service

E22546920

Causes	› BLP self-test / EEPROM error
Effect	› The pump does not run
Solutions	› Service

E22034115

Causes	› Connection problem
Effect	› The pump does not run
Solutions	› Service: check the motor cable connection

- › If the actions described fails to resolve the fault or another error code is displayed then:
 - contact the service department.
 - send the device for repair, including a short description of the fault.

Technical data

Pump rate max. (50/60 Hz)	m³/h	1.0
Pump rate max. (50/60Hz)	l/min	16.67
Final pressure	mbar	9
Suction levels		2
Cylinder		2
Connection diameter suction side	mm	8
Connection diameter pressure side	mm	8
Input pressure	mbar	10 ... 1030
Analog speed-vacuum control		yes (with controller)
Speed range	rpm	285 ... 1100 (with controller)
Noise at low pressure	dB(A)	54
Material in contact with medium		PTFE, FFPM, PPS, NBR
Housing material		alu-cast coating / thermoplastic polymer
Dimensions (W x D x H)	mm	115 x 280 x 292
Weight	kg	8.5
Permissible ambient temperature	°C	+ 5 ... + 40
Permissible relative humidity		80% (up to 31 °C), decreasing linearly to max. 50% (@40 °C)
Protection class according to DIN EN 60529		IP 20
USB interface		yes
Voltage	V	100 ... 240
Frequency	Hz	50 / 60
Power input	W	120
Power input standby	W	1.5
Fuse		2x T4A 250 V
Protection class		I
Overvoltage category		II
Contamination level		2
Operation at a terrestrial altitude	m	max. 2000
Firmware update		yes

Subject to technical changes!

Warranty

- › In accordance with IKA Terms and Conditions of Sale, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the device direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.
- › The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating instructions.

designed for scientists

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10,
79219 Staufen, Germany
Phone: +49 7633 831-0
eMail: sales@ika.de

USA

IKA Works, Inc.
Phone: +1 910 452-7059
eMail: sales@ika.net

KOREA

IKA Korea Ltd.
Phone: +82 2 2136 6800
eMail: sales-lab@ika.kr

BRAZIL

IKA Brasil
Phone: +55 19 3772 9600
eMail: sales@ika.net.br

MALAYSIA

IKA Works (Asia) Sdn Bhd
Phone: +60 3 6099-5666
eMail: sales.lab@ika.my

CHINA

IKA Works Guangzhou
Phone: +86 20 8222 6771
eMail: info@ika.cn

POLAND

IKA Poland Sp. z o.o.
Phone: +48 22 201 99 79
eMail: sales.poland@ika.com

JAPAN

IKA Japan K.K.
Phone: +81 6 6730 6781
eMail: info_japan@ika.ne.jp

INDIA

IKA India Private Limited
Phone: +91 80 26253 900
eMail: info@ika.in

UNITED KINGDOM

IKA England LTD.
Phone: +44 1865 986 162
eMail: sales.England@ika.com

VIETNAM

IKA Vietnam Company Limited
Phone: +84 28 38202142
eMail: sales.lab-vietnam@ika.com

THAILAND

IKA Works (Thailand) Co. Ltd.
Phone: +66 2059 4693
eMail: sales-lab.thailand@ika.com

TURKEY

IKA Turkey A.Ş.
Phone: +90 216 394 43 43
eMail: sales.turkey@ika.com

KENYA

IKA Works Kenya Ltd.
Phone: +254 112 323 745
eMail: sales.kenya@ika.com

UGANDA

IKA Works Kampala Limited
Phone: +254 112 323 745
eMail: sales.uganda@ika.com

SPAIN

IKA Works Spain, S. L.
Barcelona
eMail: sales.spain@ika.com

Register now: www.ika.com/register

www.ika.com

IKAworlwide // #lookattheblue

Technical specifications may be changed without prior notice.

IKA

designed for scientists

RC 2 lite

ENGLISH


Fig. 1

	EU Declaration of conformity	6
	Explication of warning symbols	6
	Safety instructions.....	7
	Intended use.....	10
	Unpacking	11
	Operator panel and display	12
	Installation.....	13
	Operation	16
	Menu navigation and structure	19
	Transportation and storage	23
	Fluids (Standard information for IKA fluid)	24
	Interfaces and outputs	25
	Maintenance and cleaning	28
	Error codes	29
	Accessories	30
	Technical Data	31
	Warranty	32
	Pump performance curve	32

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the standards or other normative documents: EN 61010-1, EN 61010-2-011, EN 61326-1, EN 60529, EN ISO 12100 and DIN 12876-1, -2, -3.

A copy of the complete declaration of conformity or further declarations of conformity can be requested at sales@ika.com.

Explication of warning symbols

/// Warning symbols

-  **Danger!** Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
-  **Warning!** Indicates a hazardous situation, which, if not avoided, can result in death, serious injury.
-  **Caution!** Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
-  **Notice!** Indicates practices which, if not avoided, can result in equipment damage.
-  **Warning!** Indicates the risk of fire or explosion!

/// General Symbols

- A** — Position number
-  Correct / Result
Shows the correct execution or the result of an action step.
-  Wrong
Shows the incorrect execution of an action step.
-  Note
Displays action steps that require particular attention to detail.

Safety instructions

/// General information

- > **Read the operating instructions completely before starting up and follow the safety instructions.**
- > Keep the operating instructions in a place where it can be accessed by everyone.
- > Ensure that only trained staff work with the device.
- > Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- > The device must only be used in a technically perfect condition.

Notice!

- > Pay attention to the marked sites in Fig. 1.

/// Device setup

Caution!

- > The power switch of the device must be accessed immediately, directly and without risk at any time. If access to the power switch cannot be ensured, an additional emergency stop switch that can be easily accessed must be installed in the work area.

Notice!

- > Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- > Insufficient ventilation may result in the formation of explosive mixtures. Only use the device in well ventilated areas.

/// Working with the device

Danger!

- > Do not use the device in explosive atmospheres, it is not EX-protected.
- > With substances capable of forming an explosive mixture, appropriate safety measures must be applied, e.g. working under a fume hood.
- > To avoid body injury and property damage, observe the relevant safety and accident prevention measures when processing hazardous materials.

Warning!

- > Process pathogenic material only in closed vessels under a suitable fume hood. Please contact IKA application support if you have any question.
- > When using critical or hazardous materials in your processes, IKA recommends to use additional appropriate measures to ensure safety in the experiment. For example, users can implement comprehensive monitoring equipment.

Caution!

- > Do not start up the device if:
 - It is damaged or leaking
 - Cable (not only the power cable) is damaged.

Notice!

- › Covers or parts that can be removed from the device without tools must later be refitted to ensure safe operation. This will prevent the infiltration of foreign objects, liquids and other contaminants.
- › The feet of the device must be clean and undamaged.
- › Always empty the bath if the unit is unused for long time.

- › A laboratory circulator refrigerates and circulates fluid according to specified parameters. This involves hazards due to high and low temperatures and general hazards due to the device of electrical energy. The user safety can not be ensured simply with design requirements on the part of the device. Further hazard sources may arise due to the type of tempering fluid, e.g. by exceeding or undercutting certain temperature thresholds or by the breakage of the container and reaction with the carrier fluid. It is not possible to consider all eventualities. They remain largely subject to the judgment and responsibility of the operator. For this reason, it become necessary for user to take other precautionary safety measures.
- › The device may only be used as prescribed and as described in these operating instructions. This includes operation by instructed specialist personnel.
- › When device is used for external circulation, extra precaution must be taken for cold fluid leakage due to damaged hose.
 - Use suitable hoses for connection.
 - Secure hoses and tubes against slippage and avoid kinks.
 - Check hoses, tubes and bath at regular intervals for possible material fatigue (cracks/leaks).
- › Do not transport or empty the bath while it is still cold.
- › Always empty the bath before moving the unit.

/// Accessories

- › Protect the device and accessories from bumping and impacting.
- › Check the device and accessories beforehand for damage each time when you use them. Do not use damaged components.
- › Safe operation is only guaranteed with the accessories described in the "Accessories" section.

/// Power supply / switching off the device

- › The voltage stated on the type plate must correspond to the power voltage.
- › The device can only be disconnected from the power supply by pulling out the power plug or the connector plug.
- › The device must only be operated with the original power cord set.
- › The socket for the power cord must be easily accessible.
- › Socket must be earthed (protective ground contact).
- › After a power failure during operation, the device may start automatically (depending on operating mode).
- › Disconnect the power plug before attaching or changing any accessories.
- › Disconnect the power plug before cleaning, maintenance and transportation of the unit.

/// Maintenance

- › The device must only be opened by trained specialists, even during repair. The device must be unplugged from the power supply before opening. Live parts inside the device may still be live for some time after unplugging from the power supply.

/// Disposal instructions

- › The device, accessories and packaging must be disposed of in accordance with local and national regulations.

/// Fluid

Warning!

- › Only use the fluids, which fulfill the requirements for safety, health and device compatibility. Be aware of the chemical hazards that may be associated with the bath fluid used. Observe all safety warning for the fluids.
- › Depending on the bath fluid used and the type of operation or toxic can arise. Ensure suitable extraction.
- › Do not use any fluid which may cause dangerous reactions during processing.

Notice!

- › Only use recommended bath fluid. Only use non-acid and non corroding fluid.
- › Never operate the device without sufficient fluid! Check the fluid level detection at a regular basis.

- › Continuous monitoring of the filling level of the bath fluid is required.
- › To ensure a sufficient fluid circulation, the viscosity of the bath fluid must not exceed of 50 mm²/s at the lowest operating temperature.
- › Untreated tap water is not recommended. It is recommended to use distilled water or high purity water (ion exchangers) and add 0.1 g soda (sodium carbonate Na₂CO₃) /liter, to reduce corrosive properties.
- › Don't use following fluids:
 - Untreated tap water
 - Acids or bases
 - Solutions with halides: chlorides, fluorides, bromides, iodides or sulfur
 - Bleach (Sodium Hypochlorite)
 - Solution with chromates or chromium salts
 - Glycerine
 - Ferrous water.

/// Refrigerant

Warning!

- > The device is not allowed to be used in an ATEX (Atmosphere Explosible) zone.
- > Flammable refrigerant used. Do not use mechanical devices or other means to accelerate the defrosting process, other than those recommended by the manufacturer.
- > Do not damage the refrigerant circuit.
- > Insufficient ventilation may result in the formation of explosive mixtures. Only use the device in well ventilated areas and keep all air vents open.
- > Flammable refrigerant used. Only the manufacturer may open and repair the closed refrigeration system.
- > The minimum room size is 9.5 m³, so that in the event of an error (exit of the complete refrigerant filling quantity) no flammable refrigerant-air mixture can arise (R 290 max. 8 g/m³ ambient air).
- > The device has been constructed to the requirements of EU (European Union) and EFTA (European Free Trade Association) countries.

Intended use

/// Use

Use **RC 2 lite** for tempering and circulating fluids.

Intended use: Tabletop device

/// Range of use

Indoor environments similar to that a laboratory of research, teaching, trade or industry.


The safety of the user cannot be guaranteed:

- > if the device is operated with accessories that are not supplied or recommended by the manufacturer.
- > if the device is operated improperly or contrary to the manufacture's specifications.
- > if the device or the printed circuit board are modified by third parties.

Unpacking


- > Please unpack the device carefully.
- > Any damage should be notified immediately to the shipping agent (post office, railway network or logistics company).

1	RC 2 lite	5	Hose olives DN 8 (2 x)
2	Power cord set	6	User guide
3	USB cable	7	Warranty card
4	Hexagon socket offset screw key		

Operator panel and display


A	On/Off button	Switch on/off the device.
B	"Set" LED	The LED lights up at the same time as the set value is displayed.
C	Cooling LED	The light LED indicated the cooling function is activated.
D	LED display	Display the settings and actual temperature values.
E	"Menu" button	Press it once: menu option is display. Press it once again: back to working screen.
F	Minus (-) button	Decrease the temperature setting. Navigation, selecting the settings in the menu. Change the menu value setting.
G	Plus (+) button	Increase the temperature setting. Navigation, selecting the settings in the menu. Change the menu value setting.
H	"OK/Temp" button	Start/Stop the tempering function. Confirm the menu options.



Installation



- › Place the unit on an even, stable, clean, non-slip, dry and fireproof surface.
- › Keep at least 20 cm of open space on the front and rear side.
- › The place for installation should be large enough and provide sufficient air ventilation to ensure the room does not warm up excessively because of the heat from device radiates to the environment.
- › Do not set up the device in the immediate vicinity of heat sources and do not expose to sun light.
- › Cooling machine, pump motor and electronics produce intrinsic heat that is dissipated via the venting grids! Never cover these venting grids!

/// Connecting the tubings and external device



/// Filling the device

/// Connecting to power


Operation


⚠ Notice!

Before commissioning, make sure that the device has not been moved in one hour.
The socket used must be earthed (fitted with earth contact).
Observe the ambient conditions (temperature, humidity etc.) listed under "Technical Data".


/// Switching on


/// Setting temperature


/// Starting tempering function

/// Stopping tempering function

/// Switching off

Menu navigation and structure

/// Menu navigation

- › Enter the menu by pressing "Menu" button (E).
- › Press the plus (+) button (G) or minus (-) button (F) to change the menu options.
- › Press "OK/Temp" button (H) to enter the select menu option.
- › Press the plus (+) button (G) or minus (-) button (F) to change the menu value settings.
- › Confirm the menu settings by pressing the "OK/Temp" button (H).

/// Menu structure

		Default settings
Menu	MODE (M)	
	\bar{n} A	A ----- Activated
		B ----- -
		C ----- -
	Maximum temperature (H)	40 °C
	H 40	
	Minimum temperature (L)	-10 °C
	L -10	
	Fluid type (F)	8
	F 8	
	The type of temperature PID control arithmetic (A)	
	A 1	AUTO 1 ----- Activated
		AUTO 0 -----
		Kp 1 (P) ----- 40
		Ti 1 (i) ----- 10.0
	Td 1 (d) ----- 3.0	
	d 3.0	
Alarm and key tone (b)		
b 1	0 ----- -	
	1 ----- Activated	
Calibration (C)		
C 0	0 ----- Activated	
	2 ----- -	
	3 ----- -	
Reset to factory settings (r)	0	
r 0		

/// Menu details

Operating mode (M):

A **A**

Mode A: After power on / power failure no automatic restart of functions.

Mode B: After power on / power failure automatic restart of functions, depending on previous settings.

Mode C: Set values (set in mode A or B) cannot be changed. After power on/power failure automatic restart of functions, depending on previous settings.

Maximum temperature (H):

H **40**

The maximum adjustable value: 40 °C.

Minimum temperature (L):

L **-10**

The minimum adjustable value: -10 °C.

Fluid type (F):

F **8**


The chosen fluid (No.) limits the maximum and minimum temperature. Customized fluid (No.: 8) enables the whole working temperature range. See Section "Fluids (Standard information for IKA fluid)".

Temperature control type (A 1):

A **1**

AUTO 1: AUTO 1 is the default setting. The default PID settings are used automatically.


Cooling down curve in AUTO 1 control (AUTO 1):

AUTO 0: For special requirements, AUTO 0 can be used with manual adjustment of the control parameters.

If the PID control (AUTO 0) option is selected, the following items including **Kp 1**, **Ti 1** and **Td 1** could be set. Otherwise they are not displayed in the menu list.

Note: Inappropriate settings may produce the following cooling down curves:

Proportional coefficient of PID (Kp 1):

P **40**

The proportional coefficient **Kp** is the controller amplification and determines how strongly the control deviation (the difference between the target temperature and actual temperature) directly affects the control variable (cooling down time). **Kp** values that are too large can lead to the controller overshooting.

Integral time of PID (Ti 1):

I **10.0**


The integral time **Ti (s)** is the correction time and determines how strongly the duration of the control deviation affects the control variable. **Ti** compensates an existing control deviation. A high **Ti** means a smaller and slower effect on the control variable. **Ti** values that are too small can lead to instability of the controller.

Differential time of PID (Td 1):

d **3.0**


The differential time **Td (s)** is the derivative time and determines how strongly the rate of change of the control deviation affects the control variable. **Td** compensates for rapid control deviations. A high **Td** means a smaller and slower effect on the control variable. **Td** values that are too small can lead to instability of the controller.

Alarm and button tone (b):

- b 0:** no alarm and button tone
- b 1:** alarm and button tone is activated

Calibration and adjustment (C):

- C 0:** Reset calibration
- C 2:** 2-point calibration
- C 3:** 3-point calibration.

Example: 2-point calibration:

Dip the temperature sensor of the reference measuring instrument into the bath fluid. Select 2-point calibration with plus (+) button / minus (-) button and "OK/Temp" button in the menu.

Press the "OK / Temp" button to start the 2-point calibration. The last calibration temperature value flashes.

Set the first point temperature (step 1, e.g. 10 °C) with plus (+) button / minus (-) button.

Confirm the setting by pressing "OK/Temp" button. The unit now starts and controls to the set value.

When the temperature has reached the set value and has become constant, the value in the screen changes from the set temperature value to "0.0" and flashes.

Input the calibration value from the reference measuring instrument (e.g. 10.3 °C) with plus (+) button / minus (-) button.

Note: the difference between the input calibration value from reference measuring instrument and set value is limited to less than +/-3 K.

Confirm the value by pressing the "OK/Temp" button. The first point calibration is finished now. The screen displays the set value for the first point calibration and flashes.

Calibration of the other points is performed in the same way.


Reset to factory settings (r):

- r 0:** Without resetting to factory settings.
- r 1:** Resetting to factory settings.

Transportation and storage

- › Remove the power cord before transporting the device.
- › The device must be protected against mechanical impact, vibrations, dust deposits and corrosive ambient air during transportation and storage.
In addition, take care to ensure that the relative humidity does not exceed 80 %.
- › The device must be completely emptied before transportation and storage.

- › Beware of the high dead weight of the device when transporting.
- › Ensure that your fingers do not get crushed when setting down the device.
- › The device can be transported:
 - by hand / unaided (2 to 4 persons).
 - with a pallet truck.

⚠ Caution!

Beware of tipping over because of the high center of gravity when lift and move the device!


Fluids (Standard information for IKA fluid)

No.	IKA Designation	Operating temperature range for open bath application (°C)	Operating temperature range for closed bath applications (°C)	Safety temperature (°C)	Flash point (°C)
0	CF.EG28.N10.80.8	-10 ... 80	-10 ... 80	90	115
1	CF.EG39.N20.80.16	-20 ... 80	-20 ... 80	90	115
2	CF.EG44.N25.80.19	-25 ... 80	-25 ... 80	90	115
3	CF.EG48.N30.80.22	-30 ... 80	-30 ... 80	90	115
4	UF.Si.N30.150.10LV	-30 ... 130	-30 ... 150	145 ①	>170
5	HF.Si.20.200.50	20 ... 200	20 ... 200	255	>280
6	HF.Si.20.250.50A	20 ... 200	20 ... 250	255	>280
7	Water ②	5 ... 95	5 ... 95	-	-
8	Customized ③				

Check the suitability of the liquid according to your application.

Nomenclature for IKA fluids:

CF.EG28.N10.80.8 --

- (1) Classification:
 HF: Heating Fluid
 CF: Cooling Fluid
 UF: Universal Fluid
- (2) Chemical composition:
 Si: Silicone oil
 EG: Ethylene Glycol
- (3) Temperature range: (Minimum temperature. Maximum temperature)
 N: Negative Temperature
- (4) Viscosity:
 Viscosity at 25 °C for Heating Fluid (HF)
 Viscosity at -20 °C for Cooling Fluid (CF)
 Viscosity at 25 °C for Universal Fluid (UF)
 Dynamic viscosity [mPa•s] is a product of kinematic viscosity [mm²/s] and density [kg/m³] of the fluid divided by 1000.
- (5) Additional information:
 A: Oil Additives
 LV: Low Viscosity


- ① **Note:** for open bath application!
- ② **Note:** Tap water may be unsuitable for operation because the calcium carbonate content may cause calcification.
 High purity water (from ion exchangers) and distilled or bi-distilled water are unsuitable for operation due to corrosive properties of these media. High purity water and distillates are suitable as a medium after adding 0.1 g soda (Na₂CO₃, sodium carbonate) per liter of water.
- ③ **Note:** The temperature limit values are adjustable in accordance with the fluid used.

Interfaces and outputs

The device can be connected to a PC and operated with the laboratory software labworldsoft® through USB or RS 232 interface.

The device software can also be updated with a PC via the RS 232 or USB port.

Note: Please comply with the system requirements together with the operating instructions and help section included with the software.

/// USB interface

The Universal Serial Bus (USB) is a serial bus for connecting the device to the PC. Equipped with USB devices can be connected to a PC during operation (hot plugging). Connected devices and their properties are automatically recognized.

/// USB device drivers

First, download the latest driver for IKA devices with USB interface from:

www.ika.com/ika/lws/download/usb-driver.zip.

Install the driver by running the setup file. Connect the IKA device through the USB data cable to the PC. The data communication is via a virtual COM port. Configuration, command syntax and commands of the virtual COM ports are as described in RS 232 interface.

Note: the USB driver is unnecessary for Windows 10 system. Don't install it for Windows 10 system!

/// RS 232 interface

Configuration:

- › The functions of the interface connections between the device and the automation system are chosen from the signals specified in EIA standard RS 232 in accordance with DIN 66 020 Part 1.
- › For the electrical characteristics of the interface and the allocation of signal status, standard RS 232 applies in accordance with DIN 66 259 Part 1.
- › Transmission procedure: asynchronous character transmission in start-stop mode.
- › Type of transmission: full duplex.
- › Character format: character representation in accordance with data format in DIN 66 022 for start-stop mode. 1 start bit; 7 character bits; 1 parity bit (even); 1 stop bit.
- › Transmission speed: 9600 bit/s.
- › Data flow control: none
- › Access procedure: data transfer from the device to the computer takes place only at the computer's request.

/// Command syntax and format

The following applies to the command set:

- › Commands are generally sent from the computer (Master) to the lab device (Slave).
- › The lab device only sends at the computer's request. Even fault indications cannot be sent spontaneously from the lab device to the computer (automation system).
- › Commands are transmitted in capital letters.
- › Commands and parameters including successive parameters are separated by at least one space (Code: hex 0x20).
- › Each individual command (incl. parameters and data) and each response are terminated with CR LF (Code: hex 0x0d hex 0x0A) and have a maximum length of 80 characters.
- › The decimal separator in a number is a dot (Code: hex 0x2E).

The above details correspond as far as possible to the recommendations of the NAMUR working party (NAMUR recommendations for the design of electrical plug connections for analogue and digital signal transmission on individual items of laboratory control equipment, rev. 1.1).


The NAMUR commands and the additional specific IKA commands serve only as low level commands for communication between the device and the PC. With a suitable terminal or communications program these commands can be transmitted directly to the device. The IKA software *labworldsoft®*, provides a convenient tool for controlling the device and collecting data under MS Windows.

Commands	Function
IN_PV_2	Read the internal actual temperature
IN_SP_1	Read the internal setting temperature
OUT_SP_1 xxx	Set the internal setting temperature XXX
OUT_SP_12@n	Set the WD safety temperature with echo of the set (defined) value.
OUT_WD1@n	Start the watchdog mode 1 and set the time for the watchdog to n (20...1500) seconds. Echo of the Watchdog time. During a WD1-event, the tempering and pump functions are switched off. This command needs to be sent within the watchdog time.
OUT_WD2@n	Start the watchdog mode 2 and set the watchdog time to n (20...1500) seconds. Echo of the watchdog time. During a WD2-event, the set temperature is changed to the WD safety temperature and the pump set speed is set to the WD safety speed. This command needs to be sent within the watchdog time.
RESET	Reset the PC control and stop the device functions.
START_1	Start the tempering function.
STOP_1	Stop the tempering function.

/// Connections between device and external devices


PC 1.1 Cable:

This cable is required to connect RS 232 interface to a PC.

USB 2.0 Cable (A - B):

This cable is required to connect USB interface to a PC.

/// Device software update

Keep your device up-to-date with the IKA Firmware update tool.

The firmware update can be done with a computer connected through USB-Interface.

For this, you need register on our website MyIKA first.

After registering your device IKA will inform you about available updates for your devices.


Please download the software "FWUToolSetup.zip" from our IKA service website www.ika.com.

Maintenance and cleaning

To avoid contamination, it is essential to change the bath fluid regularly. If water is used for bath fluid, we recommend to add Water bath protective media. The protective media stops the growth of algae, bacteria and other microorganisms.

To keep the full cooling performance, the dust filter of the cooler must be checked regularly and cleaned if necessary.

- › Switch off the device and disconnect power cable.
- › Take out the dust filter.

- › Clean condenser dust protection sheet with a vacuum cleaner or wash the sheet in the water and dry the sheet before assembly.

Note: Don't touch the condenser surface with hard parts.

/// Cleaning

For cleaning disconnect the mains plug!

Use only cleaning agents which have been approved by IKA to clean the devices:
Water containing surfactant / isopropyl alcohol.

- › Wear protective gloves during cleaning the devices.
- › Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- › Do not allow moisture to get into the device when cleaning.
- › Before using another than the recommended method for cleaning or decontamination, the user must ascertain with that this method does not destroy the device.

/// Ordering spare parts

When ordering spare parts, please indicate:

- › Device type.
- › Serial number, see type plate.
- › Item and designation of the spare part, see: www.ika.com, spare parts diagram and spare parts list.
- › Software version (briefly visible in the display when the device is switched on).

/// Repairs

Please send in device for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For repair, please request the "Decontamination Certificate" from **IKA**, or download printout of it from the **IKA** website www.ika.com.

If you require servicing, return the device in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Error codes

Any malfunctions during operation will be identified by an error message on the display.

Proceed as follows in such cases:

- › Turn off the device by using the power switch.
- › Carry out corrective measures.
- › Restart the device.

E 15

Causes	› Safety shutdown of the compressor.
Effect	› Pump off › Tempering function off
Solutions	› Check operating conditions according to the "Technical data". › Wait until over-temperatures / -pressures are equilibrated. › Restart the device.

E 16

Causes	› Compressor suction temperature fault.
Effect	› Pump off › Tempering function off
Solutions	› Check operating conditions according to the "Technical data". › Wait until over-temperatures / -pressures are equilibrated. › Restart the device.

E 17

Causes	› Compressor discharge temperature fault.
Effect	› Pump off › Tempering function off
Solutions	› Check operating conditions according to the "Technical data". › Wait until over-temperatures / -pressures are equilibrated. › Clean the filter.

E 19

Causes	› Fluid temperature too high.
Effect	› Pump off › Tempering function off
Solutions	› Check operating conditions according to the "Technical data". › Wait until over-temperatures / -pressures are equilibrated. › Switch the device off and let the fluid cool down.

E 22

Causes	› Motor is not connected
Effect	› Pump off › Tempering function off
Solutions	› Send device for repair.

If the actions described fails to resolve the fault or another error code is displayed then take one of the following steps:

- › contact the service department.
- › send the device for repair, including a short description of the fault.

Accessories

/// Tubings, hoses and adapters

LT 5.20	Metal hose (isolated M16 x 1)
LT 5.21	PTFE hose (isolated M16 x 1)
H.PVC.8	PVC tube (nominal width 8 mm)
H.PVC.12	PVC tube (nominal width 12 mm)
H.SI.8	Silicone tube (nominal width 8 mm)
H.SI.12	Silicone tube (nominal width 12 mm)
H.PUR.8	PUR tube (nominal width 8 mm)
H.PUR.12	PUR tube (nominal width 12 mm)
H.FKM.8	FKM tube (nominal width 8 mm)
H.FKM.12	FKM tube (nominal width 12 mm)

Hose olive DN 6 (2 pieces)
Hose olive DN 8 (2 pieces)
Hose olive DN 10 (2 pieces)
Hose olive DN 12 (2 pieces)

/// Tubing insulations

ISO. 8	Insulation (8 mm)
ISO.12	Insulation (12 mm)

/// Additional accessories

PC 1.1	Cable (RS 232)
USB 2.0 cable (A - B) Labworldsoft®	

See more accessories on www.ika.com.

Technical data

Operating voltage	VAC	230 ± 10 % 100 ... 115 ± 10 %
Frequency	Hz	50/60
Max. input power	W	250
Fuse		T 10 A 250 V (2 x)
Working temperature range	°C	- 10 ... RT
Operating temperature range (with external heating)	°C	- 10 ... + 70
Temperature stability – internal temperature control, ethylene glycol (according to DIN 12876)	K	± 0.5
Temperature control		PID
Temperature measurement, absolute accuracy	K	± 0.5
Temperature setting		Buttons
Temperature setting resolution	K	0.1
Temperature display		LED
Temperature display resolution		0.1
Cooling capacity according to DIN 12876: + 20 °C + 10 °C 0 °C - 10 °C	W	400 350 250 140
Refrigerant		R 290 ●
Refrigerant quantity	g	76
Max. refrigerating system pressure	bar	21
Max. pump pressure / suction	bar	0.35/0.15
Max. flow rate (at 0 bar)	l/min	18
Bath volume	l	1 ... 3.5
Fluid maximum viscosity	mm ² /s	50
Low fluid level protection		Yes
Interface		USB, RS 232
Permitted on-time	%	100
IP code according to EN 60 529		IP 21
Protection class		I
Overvoltage category		II
Contamination level		2
Permissible ambient temperature	°C	+5 ... +32
Permissible relative humidity	%	80
Dimensions (W x D x H)	mm	225 x 430 x 385
Weight	kg	24.5
Operation at a terrestrial altitude	m	max. 2000

● **Note:** Refrigerant must be disposed of in accordance with local and national regulations.

Subject to technical changes!


Warranty

In accordance with IKA warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs. The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Pump performance curve

Pump performance curve measured with water:

(Measurements done according DIN 12876-2 with water at 20 °C; pump in a closed-loop circuit).

designed for scientists

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10, 79219 Staufen, Germany
Phone: +49 7633 831-0, Fax: +49 7633 831-98
eMail: sales@ika.de

USA

IKA Works, Inc.
Phone: +1 910 452-7059
sales@ika.net

KOREA

IKA Korea Ltd.
Phone: +82 2 2136 6800
sales-lab@ika.kr

BRASIL

IKA Brasil
Phone: +55 19 3772-9600
sales@ika.net.br

MALAYSIA

IKA Works (Asia) Sdn Bhd
Phone: +60 3 6099-5666
sales.lab@ika.my

CHINA

IKA Works Guangzhou
Phone: +86 20 8222 6771
info@ika.cn

POLAND

IKA Poland Sp. z o.o.
Phone: +48 22 201 99 79
sales.poland@ika.com

JAPAN

IKA Japan K.K.
Phone: +81 6 6730 6781
info_japan@ika.ne.jp

INDIA

IKA India Private Limited
Phone: +91 80 26253 900
info@ika.in

UNITED KINGDOM

IKA England LTD.
Phone: +44 1865 986 162
sales.england@ika.com

VIETNAM

IKA Vietnam Company Limited
Phone: +84 28 38202142
sales.lab-vietnam@ika.com

THAILAND

IKA Works (Thailand) Co. Ltd.
Phone: +66 2059 4690
sales.lab-thailand@ika.com

TURKEY

IKA Turkey A.S
Phone: +90 216 394 43 43
sales.turkey@ika.com

KENYA

IKA Works Kenya Ltd.
Phone: +254 112 323 745
sales.kenya@ika.com

UGANDA

IKA Works Kampala Limited
Phone: +254 112 323 745
sales.uganda@ika.com

SPAIN

IKA Works Spain, S. L.
Barcelona
sales.spain@ika.com

Discover and order the fascinating products of IKA online:
www.ika.com

IKAworlwide

IKAworlwide /// #lookattheblue


@IKAworlwide

IKA

designed for scientists

RV 3 eco

ENGLISH

	EU Declaration of conformity	6
	Explication of warning symbols	6
	Safety instructions	7
	Intended use	11
	Unpacking	12
	Useful information	13
	Assembly	14
	Operator panel	21
	Operation	22
	Maintenance and cleaning	26
	Accessories	26
	Troubleshooting	27
	Technical Data	28
	Warranty	28
	Solvent table (excerpt)	29

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Explication of warning symbols

/// Symbols used in this instruction manual

“Safety instructions” chapter

“Unpacking” chapter

“Accessories” chapter

“Assembly” chapter

“Operation” chapter

“Troubleshooting” chapter

“Maintenance and cleaning” chapter

“Technical data” chapter

Item numbers
Indicates relevant device components for actions.

Correct / result
Indicates correct implementation or the result of an action step.

Incorrect
Indicates incorrect implementation of an action step.

Please note
Indicates action steps for which particular attention must be paid to a particular detail.

Beep
Indicates action steps, for which beep sounds are to be heard.

Power outage
Indicates device restart after a power outage.

Internet
Indicates action steps for which further information is included in the online instruction manual.

Indicates lists.

Indicates action steps.

Safety instructions

/// Warnings used in this instruction manual

Risk of burning at hot surfaces

- ▷ Do not touch hot surfaces without protective gloves.
- ▷ Only fill or empty the device when it is cold.
- ▷ Let the components cool off before continuing work with the device.

Falling glass vessels

- ▷ Check that the glass vessel is securely fixed with the clamp before operating the device.

Caution - glass may break

For the evaporating equipment (evaporating flask plus contents) the maximum permissible weight is 3.0 kg. Loads exceeding this risk breakage of glass at the steam pipe.

- ▷ When working with large loads, always use low speeds. Unbalanced loads can result in breakage of the steam pipe.

Risk of crushing when setting the height

- ▷ Check that no objects or devices are located in the area below the hand lever.

- ▷ Before commissioning, read and observe the safety instructions.
- ▷ Keep the operating instruction manual accessible for all users.

/// Other risks

Risk of injury from breakage of glass

Incorrect assembly, mechanical influences from outside or local temperature peaks may produce critical stress levels in the glass.

- ▷ Make sure that the glassware assembly is free of stress.
- ▷ Do not subject the device to overpressure during operation.
- ▷ The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at condenser) in order to prevent a pressure build-up.
- ▷ When using oil as the tempering medium, make especially sure that personal protective equipment is fitted properly.

Risk of injury from breakage of glass with vacuum operation

The glassware is designed to be used at a vacuum to 2 mbar. Failure to observe this may cause breakage of glass.

- ▷ Condense out uncondensed vapours or discharge without risk.
- ▷ If there is a risk that the distillation residue will decompose in the presence of oxygen, only ventilate with inert gases.

Risk of injury through operation with negative pressure

A power failure may produce negative pressure in the glassware and result in injury following breakage of glass.

- ▷ Ventilate glassware manually.

Risk of injury from breakage of glass if weight is exceeded

For the evaporating equipment (evaporating flask plus contents) the maximum permissible weight is 3.0 kg. Loads exceeding this risk breakage of glass at the steam pipe.

- ▷ When working with large loads, always use low speeds. Unbalanced loads can result in breakage of the steam pipe.

Risk of injury from breakage of glass with collision

Depending on the size of the flask, the setting angle of the rotation drive and the position of the heating bath and lift, the evaporator flask may be in contact with the bottom of the heating bath.

- ▷ Set variable end stop to limit the lower lift position.

Risk of injury through toxic / corrosive gases or suspended particles

The open outlet at the cooler may allow gases, vapours or suspended particles to escape in hazardous concentrations.

- ▷ Make sure that any risk is ruled out, e.g. with downstream cold traps, gas wash bottles or an efficient extraction system.

Risk of burns from glass vessels

The evaporator flask and glassware assembly can become hot during operation and remain so for a long time.

- ▷ Do not touch hot glass vessels without personal protective equipment.
- ▷ Allow components to cool off before continuing work with the device.

Risk of burns from splashing with hot tempering media

Do not operate the device when the evaporator flask is rotating and the lift is raised.

- ▷ Adjust the speed of the drive to prevent tempering medium being ejected from the heating bath by the rotating evaporator flask.

Risk of crushing from rotating parts

- ▷ Do not reach into rotating parts.

Risk of injury from delay in boiling

Heating the evaporator flask in the heating bath without connecting the rotary drive may delay boiling.

- ▷ Connect rotary drive to avoid delay in boiling.
- ▷ Sudden foaming or exhaust gases indicate that flask content is beginning to decompose. Switch off heating immediately.
- ▷ Lift the evaporator flask out of the heating bath using the lifting mechanism and ventilate the system to normal pressure as appropriate.
- ▷ Clear the hazardous area and warn any persons in the vicinity.

Risk of injury caused by electric shock

- ▷ Only put together accessories with the power supply cable disconnected.

Risk of injury caused by explosion

- ▷ Do not use the device in explosive atmospheres, as it is not EX-protected.
- ▷ When using tempering media that could create explosive mixtures, take appropriate safety measures (e.g. work under an extractor hood).

Risk of injury through explosive decomposition of peroxides

Organic peroxides can accumulate in distillation and exhaust residues and explode while decomposing. Numerous organic compounds are liable to form peroxides, e.g. dekalin, diethyl ether, dioxane, tetrahydrofuran (also unsaturated hydrocarbons such as tetralin, diene, cumene and aldehydes, ketones and solutions of these substances).

- ▷ Avoid peroxide formation.
- ▷ Keep liquids that are liable to form organic peroxides away from light (in particular UV rays).
- ▷ Check liquids for peroxides before distillation and exhaust.
- ▷ Eliminate any peroxides present.

Risk through occurrence of electrostatic discharge

Electrostatic discharge may occur between the medium and drive unit, possibly leading to a risk from static electricity (e.g. risk of explosion).

- ▷ Rule out any interaction between the medium and drive unit prior to operation.

Risk of damage through imbalance

Imbalance may result in uncontrolled resonance behaviour of the device or assembly. Glass apparatus may be damaged or destroyed.

- ▷ Immediately switch off the device with any imbalance or unaccustomed noise or reduce speed.

Risk of damage to glass vessels

If the glass vessel is not properly secured, it may fall off.

- ▷ Check that the glass vessel is securely fixed before operating the device.

Risk of damage through incorrect use of distillate

- ▷ Adapt the quantity and type of the distillation material to the size of the distillation apparatus.
- ▷ Check that the cooler is sufficiently effective.
- ▷ Monitor the flow of coolant at the cooler outlet.

/// Requirements for operators / users

- ▷ Ensure that only qualified persons use the device.
- ▷ Ensure that the applicable safety and accident prevention measures are observed when processing hazardous materials.
- ▷ Wear personal safety equipment in accordance with the hazard class of the medium that is used.
- ▷ Inform the user about possible danger due to contact with or inhalation of media (e.g. poisonous liquids, gases, mist, vapours, dusts or biological or microbiological materials).
- ▷ Only work with media for which the energy input during processing is harmless. This also needs to take into account other sources of energy such as, for example, light irradiation.
- ▷ Ensure that the specified settings on the rating plate coincide with the actual power supply.
- ▷ Ensure that the power socket used is earthed (grounded contact).
- ▷ The outlet for the mains plug must be easily accessible.
- ▷ Place device on an even, stable, clean, non-slip, dry, fire-resistant surface.
- ▷ Ensure that there is sufficient space above the device as the glass assembly may exceed the height of the device.
- ▷ Set up the device so that ventilation slots are not covered.
- ▷ Position the positive pressure outlet of the vacuum pump under a fume hood.
- ▷ Make sure that the device is not operated while unattended.
- ▷ Only operate the device under a hood closed on all sides or a comparable protective device.
- ▷ Refer to the operating instructions for the accessories.
- ▷ Always check the device and accessories for damage prior to each use. Never use damaged parts.

Intended use

/// Use

- ▷ The rotary evaporator is a tabletop device. It is suitable for the following applications:
 - ▷ Quick and gentle distillation of liquids,
 - ▷ Evaporation of solutions and suspensions,
 - ▷ Crystallization, synthesis or cleaning of fine chemicals,
 - ▷ Drying of powder and granulate material,
 - ▷ Recycling of solvents.

We recommend operation with IKA heating baths. Safe operation is guaranteed only with the use of original IKA accessories. Ensure that the device is only opened by a qualified technician.

/// Area of use

- ▷ Indoor environments similar to that a laboratory of research, teaching, trade or industry area.
- ▷ The safety of the user cannot be guaranteed:
 - ▷ if the device is operated with accessories that are not supplied or recommended by the manufacturer,
 - ▷ if the device is operated improperly or contrary to the manufacturer's specifications,
 - ▷ if the device or the printed circuit board are modified by third parties.

/// Further remarks

- ▷ Legal bases and country-specific characteristics, especially for medicine and food:
The operator is responsible for compliance with the laws and regulations applicable to him.


Unpacking


/// Unpacking

▷ Unpack the device carefully. Any damage should immediately be reported to the carrier (mail, rail or freight forwarding company).

/// Scope of delivery

A	lift
B	power cord
C	seal
D	vapour tube
E	(cylindrical) screw
F	(serrated) washer
G	base
H	annular spring / union nut
I	user guide

J	connection	N	receiving flask
K	clamp NS29 plastic	O	evaporator piston
L	condenser	P	ball joint clamp
M	(PTFE) tube	Q	stopcock

Useful information

Distillation is a thermal separating process for liquid compounds based on substance-specific, pressure-dependent boiling points through evaporation and subsequent condensation.

The boiling point temperature decreases with decreasing external pressure which means that work is usually done under reduced pressure. In this way the heating bath can be maintained at a constant temperature (e.g. 60 °C).

Using the vacuum, the boiling point is set with a steam temperature of approx. 40 °C. The cooling water for the condensation condenser should not be warmer than 20 °C (60-40-20 rule).

A chemical resistant membrane pump with a vacuum controller should be used to create the vacuum. The pump is protected from solvent residue by the addition of a condensate flask and/ or a vacuum separator.

Working with a jet pump to create a vacuum can only be recommended to a limited extent as the solvents may contaminate the environment when using these systems.

Speed, temperature, flask size and system pressure all affect the evaporator air capacity.

The optimum capacity of the flow-through condenser is approx. 60%.

This corresponds to condensation on approx. 2/3 of the cooling coil. With larger capacities there is the risk that the uncondensed solvent vapor will be extracted.

The device is designed for operation with a cooling water supply system (e.g. laboratory thermostat), but can also be run off a water supply line. Please refer to the "Technical Data" for information on cooling water pressure, temperature constancy, and flow rate.


Assembly

/// RV 3 eco drive

▷ Mount the base to the lift

Apply the cylindrical M5x25 (E) (4x) and the serrated washer M5 (F) (4x) to mount the base (G) with the lift unit (A) after you unpack the package.


▷ Fix the switching adapter

Attach the switching adapter on the rear side and secure it with the brackets.


- › Insert the switching adapter into the upper bracket.
- › Bend down the lower bracket and push the switching adapter into it.
- › Connect the low-voltage connector to the connection socket and fix it with the clamp.
- › Plug the power cord in the switching adapter and connect it to the power supply of the valid voltage.

After connecting to the power supply, the power indicator will light up.




▷ Fix the holding bracket (not included in delivery)

▷ Insert the condensate flask (not included in delivery) and attach the hose connectors to it

▷ Adjust the angle of the drive

1. Remove the clamping device for the angle setting of the rotary drive by rotating the angle adjustable knob anticlockwise according to the indication mark.
2. Set the desired angle for the rotary drive.
3. Then, secure the rotary drive against rotation by tightening the angle adjustable knob clockwise.

/// Glassware






Notice: Please read the operating instructions of the glassware for the safe handling of laboratory glassware!

▷ Locking device

- › There is a locking knob for locking or unlocking the locking device.
- › If a red mark is shown, the locking device is unlocked, otherwise it is locked.
- › To lock or unlock the locking device, push the locking knob to the end position.

▷ Feed the vapour tube

- › Unlock the locking device by turning it 60° anticlockwise according to the indication mark. The red mark should be visible.
- › Feed the vapour tube in until it stops.
- › Then, lock the locking device by turning it clockwise by 60°.
- › Push the lock knob to the end position. The red mark should be covered and invisible.
- › The vapour tube is not allowed to be pulled out!
- › Check the correct axial locking device on the vapour tube.
- › Keep the red mark invisible.

▷ **Mount & remove the evaporation flask**

Turn clockwise for locking.
Turn anticlockwise for unlocking.



Mount the evaporation flask:

- › Place the evaporation flask on the vapour tube.

- › Turn the plastic nut with the clip anticlockwise until the clip can be put on the flange.
- › Now turn the plastic nut clockwise until the clamp is hard up against the flange. When doing so, hold the locking device on the vapour tube.

Remove the evaporation flask:

- › Hold the locking device and loosen the plastic nut by turning it anticlockwise. This releases a tightly clamped evaporation flask.
- › Hold the evaporation flask at the flange and open the clip.
- › Remove the evaporation flask.
- › Check that the locking of vapour tube is still closed!

▷ **Mount the condenser seal**

- › Place the RV 10.8001 condenser seal in the condenser receptacle and fit the glassware to the device according to the assembly instructions.






▷ **Mount & remove the condenser**

Mount the condenser:

- › Place the cap nut into the condenser, and then put the spring ring to the flange of the condenser.
- › Place the condenser on the rotary drive and tighten the cap nut by hand.
- › Fix the receiving flask and hose connectors as shown.
See also operating instructions of glassware.


Remove the condenser:

- › Use the ring spanner provided to loosen union nuts that are tightly fitted.
- › Loosen the union nut by turning anticlockwise.

▷ **Mount the washer**

Notice: Pay attention to the correct position of the washer!


▷ **Description of special condensers**

(See products/accessories at www.ika.com)

- RV 10.3 Vertical-intensive condenser with manifold.
This vertical-intensive condenser features a double jacket design for particularly efficient condensation.
Also available with coating (RV 10.30)
- RV 10.4 Dry ice condenser
Dry ice condenser is for distilling low-boiling solvents.
Cooling by dry ice, no cooling water is required.
Maximum condensation thanks to low temperatures.
Also available with coating (RV 10.40)
- RV 10.5 Vertical-condenser with manifold and cut-off valve for reflux distillation.
Also available with coating (RV 10.50)
- RV 10.6 Vertical-intensive condenser with manifold and cut-off valve for reflux distillation
This vertical-intensive action condenser features a double jacket design for particularly efficient reflux distillation.
Also available with coating (RV 10.60)

/// Hose system


- ▷ Connect the water hoses (**H₂O**) to the condenser according to the counter-flow principle.
- ▷ Install the vacuum connections to the condenser, condensate flask, vacuum controller with valve and vacuum pump.
- ▷ Always connect the vacuum hose (**vac**) to the condenser at the highest point to minimize solvent losses during suctioning.
- ▷ Use standard laboratory vacuum hoses with an inside diameter of 8 mm and a wall thickness of 5 mm.

/// Heating bath

- ▷ Refer to the chapter entitled "Operation" in the operating instructions of heating bath.

Operator panel

A	Mains switch	E	Base
B	Angle adjusting knob	F	Height adjusting handle
C	Drive	G	Glassware
D	Speed adjusting knob	H	Heating bath


Operation

/// Start-up

- ▷ The unit is ready for service when the mains plug has been plugged in.
- ▷ Use the power switch to switch on or off the device.
Adjust the speed by turning the speed adjusting knob. The scale from 1 to 10 corresponds to a speed range of 20 to 300 rpm under nominal load.

Notice: Make sure that the speed adjusting knob is in the "0" position when switching on the power switch, otherwise a warning sound will be given.

/// Fill the evaporation flask

Manual filling: You can fill the evaporation flask manually prior to creating the vacuum. The evaporation flask should not be filled more than half its volume.

Automatic filling: Prior to filling the evaporation flask, a vacuum controller is used to regulate the glass apparatus to the target pressure.

- ▷ Now fill the evaporation flask using the back feed line.
- ▷ Due to the vacuum present, the solvent is suctioned into the evaporation flask. This enables you to keep solvent loss due to suctioning to a minimum.

Notice: The maximum allowed load (evaporation flask and contents) is 3 kg.


/// Set up the heating bath

Notice: Refer to the operating instructions of the heating bath!

- ▷ Move the drive to the bottom position and check the position of the heating bath in relation to the evaporation flask. When using larger evaporation flasks or depending on the angle of the rotary drive, you can move the heating bath to the right.
- ▷ Fill the heating bath with the tempering medium until the evaporation flask is surrounded by tempering medium to 2/3 of its volume.
- ▷ Switch on the rotary drive and slowly increase the speed.
Notice: Avoid creating waves.
- ▷ Switch on the heating bath using the power switch.
Notice: Avoid stress on the glass due to different evaporation flask and heating bath temperatures when lowering the evaporation flask into the heating bath!

/// Adjust the drive height

- ▷ To adjust the drive height, grasp the height adjusting handle and hold the trigger. The lift is unlocked and you can move the drive up and down. Release the trigger, the lift is locked in the adjusted height.


/// Adjust the lower end stop

Notice: Depending on the size of the evaporation flask, the setting angle of the rotary drive and the position of the heating bath and lift, the evaporation flask can be in contact with the bottom of the heating bath.

Notice: Glass may break!

- ▷ Limit the lower lift position with the variable end stop.
- ▷ Turn the knob anticlockwise to unlock the stopper, and then you may adjust the lower end stop by moving it up and down.
- ▷ Turn the knob clockwise, the stopper is locked in the adjusted lower end stop.

/// Assemble the stand pole RV 3.1 (Accessories)

Assemble the condenser locking device according to the diagram.

- ▷ Mount the plate (A) using the thumb screw (B).
- ▷ Put the support rod (D) on the plate (A) and attach it with the screw nut (C).
- ▷ Attach the rubber protector (E).
- ▷ Fasten the Velcro® band (F) to the support rod (D).
- ▷ Secure the vertical glassware with the Velcro® band (F).


No stand is required if the condenser has been properly fitted and secured by tightening the condenser union nut on the rotary drive.

The purpose of the stand is solely to prevent the condenser from twisting.

Notice:

If the stand is fitted incorrectly, the strong forces at the Velcro strap can cause stresses within the glass, which in turn can cause damage to the glass condenser.

After the condenser has been fitted, the stand is attached to the rotary drive. Make sure that the condenser is installed parallel to the lift body.

Maintenance and cleaning

- ▷ The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.
- ▷ The seal on the glass condenser should be checked at regular intervals and replaced if necessary.

/// Cleaning:

- ▷ For cleaning disconnect the mains plug!
- ▷ Use only cleaning agents which have been approved by IKA to clean IKA devices. These are water (with tenside) and isopropanol.
- ▷ Wear protective gloves during cleaning the devices.
- ▷ Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- ▷ Do not allow moisture to get into the device when cleaning.
- ▷ Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

/// Spare parts ordering:

- ▷ When ordering spare parts, please give:
 - › device type.
 - › serial number, see type plate.
 - › position number and description of spare part, see www.ika.com.
 - › software version.

/// Repair:

- ▷ Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.
- ▷ For repair, please request the “**Safety Declaration (Decontamination Certificate)**” from IKA or use the download printout of it from IKA website at www.ika.com.
- ▷ If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Note: For glassware unit, please contact your local dealer and do not send the glassware back to our factory.

Accessories

- ▷ For accessories see www.ika.com.

Troubleshooting

- ▷ Any malfunctions during operation will be identified by the error sound. Proceed as follows in such cases:
 - › Turn off the device by using the main switch.
 - › Carry out corrective measures.
 - › Restart the device.

Error sound | Causes | Effect | Solutions

Continuously beeping

Causes	› The speed has not been adjusted to 0 rpm when starting the device › PCBA failure
Effect	› Motor can't run
Solutions	› Turn the speed adjusting knob and set the speed to 0 rpm, after which you may set the required speed value.

Alternatively beep 1 sec and stop 1 sec

Causes	› Overload › Motor blocked › Defective cable › PCBA failure
Effect	› Motor stops running
Solutions	› Reduce the volume in the evaporator › Contact the service department

Alternatively beep 0.5 sec and stop 1.5 sec

Causes	› Temperature inside the device is too high › PCBA failure
Effect	› Motor stops running
Solutions	› Switch off the device › Let the device cool down › Switch on the device › Contact the service department

Alternatively beep 1.5 sec and stop 0.5 sec

Causes	› Overvoltage › Undervoltage › PCBA failure › Desktop power supply failure
Effect	› Motor stops running
Solutions	› Use the original power supply to supply the device › Contact the service department

- ▷ If the actions described fail to resolve the fault or the other kind of error occurred, then:
 - › contact the service department.
 - › send the device for repair, including a short description of the fault.

Technical data

Operating voltage range	Vac	(100...240) ± 10%
Rated voltage	Vac	100...240
Frequency	Hz	50/60
Rated power without heating bath	W	40
Rated power (max.) of switching adaptor	W	90
Speed	rpm	20...300
Speed tolerance	Set speed: < 100 rpm Set speed: ≥ 100 rpm	rpm ± 1 % ± 1
Speed display		Plastic part + scale
Smooth start		yes
Lift		manual
Stroke	mm	150
Head angle adjustable		0°...60°
Cooling surface	cm²	1500
Cooling water flow rate min.	l/h	30
Cooling water flow rate max.	l/h	100
Cooling water pressure max.	bar	1
Perm. On-time	%	100
Perm. ambient temperature	°C	5...40
Perm. relative humidity	%	80
Protection acc. to DIN EN 60529		IP 20
Protection class		I
Contamination level		2
Weight (no glassware; no heating bath)	kg	8
Dimensions (W x D x H)	mm	440 x 330 x 530
Operation at a terrestrial altitude	m	max. 2000

Subject to technical changes!

Warranty

- ▶ In accordance with IKA Terms and Conditions of Sale, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the device direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.
- ▶ The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating instructions.

Solvent table (excerpt)

Solvent	Formula	Pressure for boiling point 40 °C in mbar (For HB approx. 60 °C)
Acetic acid	C ₂ H ₄ O ₂	44
Acetone	C ₃ H ₆ O	556
Acetonitrile	C ₂ H ₃ N	226
N-Amylalcohol	C ₅ H ₁₂ O	11
n-Pentanol	C ₅ H ₁₀ O	11
n-Butanol	C ₄ H ₁₀ O	25
tert. Butanol	C ₄ H ₁₀ O	130
2-Methyl-2-Propanol	C ₄ H ₁₀ O	130
Butylacetate	C ₆ H ₁₂ O ₂	39
Chlorobenzene	C ₆ H ₅ Cl	36
Chloroform	CHCl ₃	474
Cyclohexane	C ₆ H ₁₂	235
Dichloromethane	CH ₂ Cl ₂	atm. press.
Methylenechloride	CH ₂ Cl ₂	atm. Press.
Diethylether	C ₄ H ₁₀ O	atm. press.
1,2,-Dichloroethylene (trans)	C ₂ H ₂ Cl ₂	751
Diisopropylether	C ₆ H ₁₄ O	375
Dioxane	C ₄ H ₈ O ₂	107
Dimethylformamide (DMF)	C ₃ H ₇ NO	11
Ethanol	C ₂ H ₆ O	175
Ethylacetate	C ₄ H ₈ O ₂	240
Ethylmethylketone	C ₄ H ₈ O	243
Heptane	C ₇ H ₁₆	120
Hexane	C ₆ H ₁₄	335
Isopropylalcohol	C ₃ H ₈ O	137
Isoamylalcohol	C ₅ H ₁₂ O	14
3-Methyl-1-Butanol	C ₅ H ₁₂ O	14
Methanol	CH ₄ O	337
Pentane	C ₅ H ₁₂	atm. press.
n-Propylalcohol	C ₃ H ₈ O	67
Pentachloroethane	C ₂ HCl ₅	13
1, 1, 2, 2, -Tetrachloroethane	C ₂ H ₂ Cl ₄	35
1, 1, 1, -Trichloroethane	C ₂ H ₃ Cl ₃	300
Tetrachloroethylene	C ₂ Cl ₄	53
Tetrachloromethane	CCl ₄	271
Tetrahydrofuran (THF)	C ₄ H ₈ O	357
Toluene	C ₇ H ₈	77
Trichloroethylene	C ₂ HCl ₃	183
Water	H ₂ O	72
Xylene	C ₈ H ₁₀	25

designed for scientists

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10,
79219 Staufen, Germany
Phone: +49 7633 831-0
eMail: sales@ika.de

USA

IKA Works, Inc.
Phone: +1 910 452-7059
eMail: sales@ika.net

KOREA

IKA Korea Ltd.
Phone: +82 2 2136 6800
eMail: sales-lab@ika.kr

BRAZIL

IKA Brasil
Phone: +55 19 3772 9600
eMail: sales@ika.net.br

MALAYSIA

IKA Works (Asia) Sdn Bhd
Phone: +60 3 6099-5666
eMail: sales.lab@ika.my

CHINA

IKA Works Guangzhou
Phone: +86 20 8222 6771
eMail: info@ika.cn

POLAND

IKA Poland Sp. z o.o.
Phone: +48 22 201 99 79
eMail: sales.poland@ika.com

JAPAN

IKA Japan K.K.
Phone: +81 6 6730 6781
eMail: info_japan@ika.ne.jp

INDIA

IKA India Private Limited
Phone: +91 80 26253 900
eMail: info@ika.in

UNITED KINGDOM

IKA England LTD.
Phone: +44 1865 986 162
eMail: sales.England@ika.com

VIETNAM

IKA Vietnam Company Limited
Phone: +84 28 38202142
eMail: sales.lab-vietnam@ika.com

THAILAND

IKA Works (Thailand) Co. Ltd.
Phone: +66 2059 4693
eMail: sales-lab.thailand@ika.com

TURKEY

IKA Turkey A.Ş.
Phone: +90 216 394 43 43
eMail: sales.turkey@ika.com

KENYA

IKA Works Kenya Ltd.
Phone: +254 112 323 745
eMail: sales.kenya@ika.com

UGANDA

IKA Works Kampala Limited
Phone: +254 112 323 745
eMail: sales.uganda@ika.com

SPAIN

IKA Works Spain, S. L.
Barcelona
eMail: sales.spain@ika.com

Register now: www.ika.com/register

www.ika.com

IKAworlwide // #lookattheblue

Technical specifications may be changed without prior notice.

IKA

designed for scientists

HB eco

ENGLISH


Fig. 1

	EU Declaration of conformity	6
	Explication of warning symbols	6
	Safety instructions	7
	Intended use	9
	Unpacking	10
	Operator panel and display	11
	Operation	12
	Interfaces and outputs	18
	Maintenance and cleaning	19
	Error codes	20
	Technical Data	22
	Warranty	23

EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-010; EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete Declaration of Conformity or further declarations of conformity can be requested at sales@ika.com.

Explication of warning symbols

/// Symbols used in this instruction manual

"Safety instructions" chapter

"Unpacking" chapter

"Accessories" chapter

"Assembly" chapter

"Operation" chapter

"Error codes" chapter

"Maintenance and cleaning" chapter

"Technical data" chapter

A — Item numbers
Indicates relevant device components for actions.

Correct / result
Indicates correct implementation or the result of an action step.

Incorrect
Indicates incorrect implementation of an action step.

Please note
Indicates action steps for which particular attention must be paid to a particular detail.

Beep
Indicates action steps, for which beep sounds are to be heard.

Power outage
Indicates device restart after a power outage.

Internet
Indicates action steps for which further information is included in the online instruction manual.

Indicates lists.

Indicates action steps.

Safety instructions

/// Warnings used in this instruction manual

Risk of burning at hot surfaces

- ▷ Do not touch hot surfaces without protective gloves.
- ▷ Only fill or empty the device when it is cold.
- ▷ Let the components cool off before continuing work with the device.

Burns hazard due to hot tempering media

- ▷ Do not touch hot tempering media.
- ▷ Note the danger of burns in case of tempering media that have a flashpoint of less than 260 °C.
- ▷ Note the danger resulting from flammable materials.

- ▷ Before commissioning, read and observe the safety instructions.
- ▷ Keep the operating instruction manual accessible for all users.

/// Other risks

Risk of burning at hot surfaces

The top edge of the heating bath can heat up to temperatures above 65 °C during operation.

- ▷ Do not touch hot device components without protective gloves.
- ▷ Only fill or empty the device when it is cold.

Risk of injury caused by explosion

- ▷ Do not use the device in explosive atmospheres, as it is not EX-protected.
- ▷ When using tempering media that could create explosive mixtures, take appropriate safety measures (e.g. work under an extractor hood).

Risk of injury caused by self-igniting tempering media

- ▷ The safety temperature must be set in accordance with EN 61010-2-010 Chapter "Requirements for devices containing or using flammable liquids".
The surface temperature of the flammable medium that is exposed to air may not exceed its flash point.
 - A danger usually arises if a medium is heated in open vessels.
The surface temperature of the heating device may not exceed the value of $(t - 25)$ °C (= set value of the safety circuit) on the surface of the flammable medium and in contact with air, whereby t is the fire point of the liquid.
 - A danger usually arises if a medium is heated in glass vessels (glass breakage).
If a setting made by the user (medium temperature or safety temperature) could bring a flammable medium into a state in which the conditions mentioned above could be exceeded, additional measures must be introduced that will protect the user from this danger.
- ▷ Do not touch hot tempering media.

Risk of burns due to self-igniting tempering media (also in connection with rotary evaporator or other devices)

- ▷ For tempering media with a flashpoint ≤ 260 °C: Note the danger resulting from ignition of the medium.
- ▷ Note the danger resulting from flammable materials.

Risk of injury caused by electric shock

- ▷ Only put together accessories with the power supply cable disconnected.

Damage to the device due to improper use

- ▷ Empty the heating bath before transportation.
- ▷ Never operate the heating bath without tempering medium.
- ▷ Ensure that the specified settings on the rating plate coincide with the actual power supply.
- ▷ Ensure that the power socket used is earthed (grounded contact).
- ▷ Avoid knocks or impacts to the device and accessories.

Hazards resulting from uncontrolled restart

- ▷ After an interruption to the power supply, the device starts up again with previously set operating mode C.
- ▷ To cut the power supply, press the device main switch or unplug the mains or device plug.
- ▷ After interrupting the power supply, ensure that the device is not unintentionally restarted.

/// Requirements for operators / users

- ▷ Ensure that only qualified persons use the device.
- ▷ Ensure that the device is only opened by a qualified technician.
- ▷ Ensure that the applicable safety and accident prevention measures are observed when processing hazardous materials.
- ▷ Wear personal safety equipment in accordance with the hazard class of the medium that is used.
- ▷ Inform the user about possible danger due to contact with or inhalation of media (e.g. poisonous liquids, gases, mist, vapours, dusts or biological or microbiological materials).
- ▷ Only work with media for which the energy input during processing is harmless. This also needs to take into account other sources of energy such as, for example, light irradiation.
- ▷ Place device on an even, stable, clean, non-slip, dry, fire-resistant surface.
- ▷ Refer to the operating instructions for the accessories.
- ▷ Always check the device and accessories for damage prior to each use. Never use damaged parts.
- ▷ During operation with rotary evaporator: Observe the maximum filling quantity also for the displaced piston volume.

Intended use

/// Use

- ▷ The heating bath is used for tempering media. It is the user's responsibility to determine which medium is tempered. The heating bath is not used for tempering foods. The heating bath is suitable for use with rotary evaporators from IKA.

Recommended tempering media

- ▷ Water (up to 80 °C)
- ▷ Low viscosity silicone oils (50 mPas) with a flashpoint ≥ 260 °C
- ▷ Untreated tap water is not recommended. It is recommended to use distilled water or high purity water (ion exchangers) and add 0.1 g soda (sodium carbonate Na_2CO_3)/liter, to reduce corrosive properties.

/// Area of use

- ▷ Indoor environments similar to that a laboratory of research, teaching, trade or industry area.
- ▷ The safety of the user cannot be guaranteed:
 - ▷ if the device is operated with accessories that are not supplied or recommended by the manufacturer,
 - ▷ if the device is operated improperly or contrary to the manufacture's specifications,
 - ▷ if the device or the printed circuit board are modified by third parties.


Unpacking

/// Unpacking

Unpack the device carefully. Any damage should immediately be reported to the carrier (mail, rail or freight forwarding company).


/// Scope of delivery

A	heating bath
B	power cord
C	user guide



Operator panel and display

/// Operator panel

A	Power socket	F	Bath insert
B	Mains switch	G	LED heating
C	Rotating / pressing knob	H	LED (set = setpoint)
D	USB interface	J	Display
E	Handle		

/// Display

Operation

/// Start-up

- ▷ Observe the ambient conditions (temperature, humidity, etc.) listed under “Technical data”.
- ▷ The unit is ready for service when the mains plug has been plugged in.
- ▷ Prior to commissioning the device, determine the filling quantity of the tempering medium due to the used evaporator piston size (approx. 2.5 l tempering medium by using a standard 1 l evaporator piston).


- ▷ Switch on the device using the mains switch.
- ▷ At every start, the display will show all display segments (Fig. 2), the software version (Fig. 3), and the operating mode (Fig. 4).
- ▷ The heating bath temperature is held constant by the control circuit and is additionally monitored by the safety circuit.

In the event of an error occurred in the control circuit, the heating bath is switched off permanently by the safety circuit.

An error in the control or safety circuits is shown on the display.

The heating function can no longer be started.



- ▷ When switching off the heating function and the medium temperature is higher than 50 °C, the display will alternately show the set temperature (Fig. 5) and “HOT” (Fig. 9). Otherwise, the display will alternately show the set temperature (Fig. 5) and OFF (Fig. 8). When switching on the heating function, the display will alternately show set temperature (Fig. 5) and the actual temperature (Fig. 6 or Fig. 7).

/// Setting the target temperature

Activate the SET function by turning and pressing the rotating / pressing knob.

- ▷ Set the desired heating bath temperature by turning the rotating / pressing knob.
- ▷ The LED setpoint temperature will be always on whenever the set temperature is shown. (Fig. 5)
- ▷ Confirm the value by pressing the rotating / pressing knob.
- ▷ Start the heating process by pressing the rotating / pressing knob.

/// Setting the operating modes

Operating mode A

- ▷ Switch on the device at the mains switch.
- ▷ The heating function is switched off.
- ▷ The target value is set to 20 °C.
- ▷ The target value can be adjusted.
- ▷ The heating function must be restarted following a power failure.
- ▷ Mode A appears on the display when switching on the device.

Factory setting: Mode A

Operating mode b


- ▷ Switch on the device at the mains switch.
- ▷ The heating function is switched off.
- ▷ The target value is set to 20 °C or the most recent set temperature.
- ▷ The target value can be adjusted.
- ▷ The heating function must be restarted following a power failure.
- ▷ Mode b appears on the display when switching on the device.

Operating mode C

- ▷ Switch on the device at the mains switch.
- ▷ The heating function is switched on or off depending on the last selected setting.
- ▷ The last target value set in mode b applies.
- ▷ The target value cannot be adjusted.
- ▷ The heating function will be automatically restarted following a power failure.
- ▷ Mode C appears on the display when switching on the device.

Switching the operating modes

- ▷ Switch the device off at the mains switch.
- ▷ Press and hold rotating / pressing knob and switch the device on at the mains switch. Release the rotating / pressing knob after 2 seconds.
- ▷ Sequence A, b, C, A etc.
- ▷ The operating mode appears on the display when switching on the device.

i

	A	b	C
	20 ... 99	20 ... 99	—

/// Safety temperature

The safety circuit prevents the temperature of the heating bath from getting too high as a result of controller error. Once the safety temperature has been reached, the device switches off permanently.

Furthermore, dry running of the heating bath is detected. The device switches off permanently. The dry running function detects both accidental heating of the heating bath without liquid and dry running caused by water evaporation starting at a set temperature of 60 °C. The error message E26 (Fig. 11) appears and the heating bath is switched off permanently. To remedy this error, see "Error codes".

/// Controlling the medium temperature

The medium temperature is limited by the set safety temperature controlled by means of a PID controller and measured by the PT 1000 temperature sensor, and the medium is heated up as quickly as possible without overshooting to the set temperature.

The PID controller adapts to the different tempering media and ensures optimum temperature control with minimum temperature drift and fluctuation.

Optimum control is only achieved when the medium is mixed using a rotating evaporator piston.

- ▷ Set the desired medium temperature between room temperature and 99 °C by turning the rotating / pressing knob.
- ▷ Start the heating function by pressing the rotating / pressing knob.
- ▷ The LED heating appears.
- ▷ The heating bath is heated up to the set temperature.
- ▷ The set and actual temperature referred to the medium are shown on the display.

Interfaces and outputs

- ▷ Please comply with the system requirements together with the operating instructions and help section included with the software.

/// USB interface:

The connected devices and their properties are detected automatically. The USB interface is used in conjunction with software control for “remote” operation and can also be used for software update of the device.

/// USB device drivers:

First, download the latest driver for IKA devices with USB interface from:

<http://www.ika.com/ika/lws/download/usb-driver.zip>

Install the driver by running the setup file. Then connect the IKA device to the PC via the USB data cable and follow the instructions.

The data communication is via a virtual COM port.

/// USB cable A - B:

This cable is required to connect USB port to a PC.

Maintenance and cleaning

- ▷ The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

/// Cleaning:

- ▷ For cleaning disconnect the mains plug!
- ▷ Use only cleaning agents which have been approved by IKA to clean the devices:
Water containing surfactant / isopropyl alcohol.
- ▷ Wear protective gloves during cleaning the devices.
- ▷ Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
- ▷ Do not allow moisture to get into the device when cleaning.
- ▷ Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

/// Spare parts order:

- ▷ When ordering spare parts, please give:
 - › device type.
 - › serial number, see type plate.
 - › position number and description of spare part, see **www.ika.com**.
 - › software version.

/// Repair:

- ▷ Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.
- ▷ For repair, please request the “**Safety Declaration (Decontamination Certificate)**” from IKA or use the download printout of it from IKA website at **www.ika.com**.
- ▷ If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Error codes

- ▷ The fault is shown by an error code on the display as following if the error occurs.
Proceed as follows in such cases:
- › Turn off the device by using the main switch.
 - › Carry out corrective measures.
 - › Restart the device.

Error code | Causes | Effect | Solutions

E 3

Causes	› Temperature inside the device is too high
Effect	› Heating switched off
Solutions	› Switch off the device › Let the device cool down › Switch on the device

E 9

Causes	› Fault in storing the safety circuit set temperatures › Memory chip (EPROM) defective
Effect	› Heating switched off
Solutions	› Switch off the device › Let the device cool down › Switch on the device

E21

Causes	› Safety relay does not open
Effect	› Heating switched off
Solutions	› Switch off the device and switch on again

E24

Causes	› Set safety temperature exceeded
Effect	› Heating switched off
Solutions	› Switch off the device › Let the device cool down › Switch on the device

E25

Causes	› The heater control circuit switch (TRIAC) is faulty. › The heater or the supply line is disconnected.
Effect	› Heating switched off
Solutions	› Switch off the device and switch on again

E26

Causes	› Dry running
Effect	› Heating switched off
Solutions	› Switch off the device › Let the device cool down › Fill up the tempering medium › Switch on the device

E27

Causes	› Error in calibration
Effect	› Heating switched off
Solutions	› Switch off the device and switch on again

E28

Causes	› Breakage of the controller sensor
Effect	› Heating switched off
Solutions	› Switch off the device and switch on again

E29

Causes	› Short-circuit in the safety sensor
Effect	› Heating switched off
Solutions	› Switch off the device and switch on again

E30

Causes	› Short-circuit in the controller sensor
Effect	› Heating switched off
Solutions	› Switch off the device and switch on again

E31

Causes	› Breakage of the safety sensor
Effect	› Heating switched off
Solutions	› Switch off the device and switch on again

E32

Causes	› Temperature deviation is too big
Effect	› Heating switched off
Solutions	› Switch off the device and switch on again

- ▷ If the actions described fails to resolve the fault or another error code is displayed then take one of the following steps:
- › contact the service department.
 - › send the device for repair, including a short description of the fault.

Technical data

Power supply		
Voltage	VAC	200...240 ±10 % 100...120 ±10 %
Frequency	Hz	50 / 60
Power input	W	1450
Power input standby	W	3
Heating function		
Heat control	-	LED
Set temperature resolution	K	± 1
Heating temperature min. ... max.	°C	Room temperature ... 99
Heat output	W	1400
Controller oscillation (3l water / 90 °C)	K	± 1
Absolute deviation / mean value (3l water / 90 °C)	K	± 2
Fixed safety temperature	°C	109
Filling point min.	mm	60
Filling volume max.	l	4
General information		
USB interface	-	yes
Protection class according to DIN EN 60529	-	IP 21
Overvoltage category	-	II
Permissible relative humidity	%	80
Permissible ambient temperature min. ... max.	°C	5 ... 40
Material in contact with medium	-	Stainless steel 1.4404
Dimensions (W x D x H)	mm	330 x 325 x 190
Outer diameter	mm	257
Weight	kg	3.9
Inner diameter	mm	225
Inner height	mm	130
Outer height	mm	190
Class designation acc. DIN 12876	-	I
Operation at a terrestrial altitude	m	max. 2000 above sea level

Subject to technical changes!

Warranty

- ▷ In accordance with IKA Terms and Conditions of Sale, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the device direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.
- ▷ The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

designed for scientists

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10,
79219 Staufen, Germany
Phone: +49 7633 831-0
eMail: sales@ika.de

USA

IKA Works, Inc.
Phone: +1 910 452-7059
eMail: sales@ika.net

KOREA

IKA Korea Ltd.
Phone: +82 2 2136 6800
eMail: sales-lab@ika.kr

BRAZIL

IKA Brasil
Phone: +55 19 3772 9600
eMail: sales@ika.net.br

MALAYSIA

IKA Works (Asia) Sdn Bhd
Phone: +60 3 6099-5666
eMail: sales.lab@ika.my

CHINA

IKA Works Guangzhou
Phone: +86 20 8222 6771
eMail: info@ika.cn

POLAND

IKA Poland Sp. z o.o.
Phone: +48 22 201 99 79
eMail: sales.poland@ika.com

JAPAN

IKA Japan K.K.
Phone: +81 6 6730 6781
eMail: info_japan@ika.ne.jp

INDIA

IKA India Private Limited
Phone: +91 80 26253 900
eMail: info@ika.in

UNITED KINGDOM

IKA England LTD.
Phone: +44 1865 986 162
eMail: sales.english@ika.com

VIETNAM

IKA Vietnam Company Limited
Phone: +84 28 38202142
eMail: sales.lab-vietnam@ika.com

THAILAND

IKA Works (Thailand) Co. Ltd.
Phone: +66 2059 4690
eMail: sales.lab-thailand@ika.com

TURKEY

IKA Turkey A.Ş.
Phone: +90 216 394 43 43
eMail: sales.turkey@ika.com

Discover and order the fascinating products of IKA online:
www.ika.com