







# Ammonium 400mg/L

# **DATA SHEET**

# Instructions for use

- 1. Fill the test vessel with 5mL of sample.
- 2. Add 10 drops of Reagent  $NH_4^+$  A and mix thoroughly.
- 3. Immerse an indicator strip in the test solution for 3 seconds.
- 4. Shake off excess solution.
- 5. Wait 10 seconds.
- 6. Compare the colour produced on the strip against the colour chart.

### **Measuring range**

0 - 10 - 25 - 50 - 100 - 200 - 400mg/L

# **Application**

This Ammonium test kit detects aqueous ammonium ions in wastewater as well as aquariums, where ammonium can be easily converted to ammonia which is highly toxic to fish.

1

#### **Shelf Life**

2 years from manufacturing date.



# **Contents**

100 strips 1x test vessel 2x 30mL reagent NH<sub>4</sub>+ - A

#### **Dimensions**

Strip dimension: 83mm x 5mm

# Notes

Keep any unused strips in the original container. Do not store where they are liable to be exposed to acid or alkaline gases or vapours. Store in a cool dry place. Reagent  $NH_4^+$ - A is corrosive.