



2019
EDITION



SURFACES

INSTRUMENTS

ASPIRATING UNITS &
DENTAL IMPRESSIONS

HAND HYGIENE

LAUNDRY

Swiss Made Disinfectants Since 1969



OCC Switzerland is a Swiss, family-owned manufacturer of disinfectants.


OCC Switzerland's success story started in 1969 when Heinz Suter developed his first cleaner for dental suction units. ORO-CLEAN, a phenol-free cleaner with antimicrobial properties was an instant hot seller which gave Heinz Suter the confidence not only to establish Dentalchemie Heinz Suter but also to invest a large portion of his earnings into the development of disinfectants. In 1977, together with his wife, Heinz Suter incorporated Oro Clean Chemie AG, also known as OCC Switzerland. Since then, the company has grown into a respected infection control player with subsidiaries in Germany, the Czech Republic, and Malaysia and customers in more than 60 countries around the world. OCC Switzerland's product portfolio includes disinfectants and cleaners for surfaces, instruments, hands, laundry, and industry specific areas.

Killing Germs for a safer world



1969

Chemist Heinz Suter develops ORO-CLEAN in his father's garage, a phenol-free cleaner for dental suction units.



1977

The legal form was changed to a corporation and the company renamed to Oro Clean Chemie AG.



1980

Introduction of DENTIRO®, a fast-acting surface disinfectant.



1996

OCC Switzerland became the first Swiss dental company to be certified according to EN 46001.



“ Since our early days, our goal has been to manufacture fast-acting disinfectants that are safe to use and compatible with a wide range of materials. We constantly improve our disinfectants to ensure that they are scientifically up-to-date and meet the requirements of today's infection control professionals. We, at OCC Switzerland strive to manufacture disinfectants in the most environmentally friendly manner and align ourselves with entities that support environmental sustainability. Our skilled workforce is the driver of our business and allows us to support our partners by increasing the quality and efficiency of their infection control and prevention endeavours.

Ursula Suter, CEO



1998

Relocation to newly erected production and research premises in Fehraltorf.



2006

Expansion of production and research premises in Fehraltorf.



2007

Expansion of analytical and microbiological laboratory with state-of-the-art instruments and equipment.



2014

Launch of DENTIRO® Wipes and DENTIRO® Light, two surface disinfectants that kill 101 microorganisms including poliovirus in 60 seconds.



2018

Launch of OROLIN® Multisept Plus and ISORAPID® OP-Forte AF together with the newly improved ORO CLEAN® Plus and ORO CLEAN® Liquid



Our vision

Our vision is to be one of the world's leading manufacturers of fast-acting, reliable, and safe disinfectants that are effective against a broad spectrum of viruses, bacteria, and fungi.



Our mission

At OCC Switzerland, 'Killing germs' isn't just our slogan. It's our mission.



Our values

As an organization, our core values are driven by a desire to improve human life, achieve scientific excellence, operate within the highest standards of integrity, and employ a diverse, multicultural workforce that values collaboration and acts ethically at all times.



Instruments



Careful and thorough disinfection and cleaning of medical and surgical instruments are essential to protect patients and staff from infections.

Page 14

Hands



Proper hand hygiene is the single most important measure to prevent the spread of Healthcare Acquired Infections (HAIs).

Page 21

Surfaces



Thorough cleaning and disinfection of surfaces are critical to prevent the spread of disease within healthcare settings.

Page 7

Aspirating units and dental impressions



Suction systems are heavily contaminated and pose a real danger of causing infection.

Dental impressions and casts are highly contaminated with patient saliva or blood upon removal from the patient's mouth.

Page 19

Laundry



Contaminated laundry in healthcare settings must be processed correctly to eliminate the risk of cross infection.

Page 27

DENTIRO® Light Kills 101 microorganisms

active in
60
seconds

Candida albicans

Hepatitis C
virus

Mycobacterium bovis

*Yersinia
enterocolitica*

Klebsiella pneumoniae

Burkholderia cepacia

*Enterococcus
hirae*

***Mycobacterium
terrae***

Proteus mirabilis

Bordetella pertussis

Influenza A virus

*Mycobacterium
avium*

*Neisseria
meningitidis*

Streptococcus pyogenes

Adenovirus

*Mycobacterium
smegmatis*

*Proteus
vulgaris*

Moraxella catarrhalis

*Escherichia
coli*

Listeria monocytogenes

Poliovirus

*Saccharomyces
cerevisiae*

Staphylococcus epidermidis

*Salmonella
typhimurium*

Bovine Viral Diarrhea Virus

*Vibrio
parahaemolyticus*

MRSA

*Propionibacterium
acnes*

Vancomycin-resistant
Enterococcus

Acinetobacter baumannii

Vaccinia virus

Mycobacterium chelonae

HIV

*Staphylococcus
aureus*

Herpes simplex virus

Norovirus

Hepatitis B virus

Trichoderma viride

*Enterobacter
aerogenes*



CONTACT US FOR A COMPLETE LIST
TEL. +41 (0)44 226 44 44
EMAIL: INFO@OROCLEAN.COM

READY-TO-USE SURFACE DISINFECTANTS



DENTIRO® Sensitive

Fully virucidal, ready-to-use surface disinfectant for the quick disinfection and cleaning of non-invasive medical devices (medical equipment, installations, and instruments). Lightly scented.

- Kills 101 microorganisms
- Active in 60 seconds
- Compatible with 250 materials including artificial leather
- Low alcohol content of less than 50%

Composition in 100 g: 20 g ethanol, 28 g 1-propanol, 0.056 g quaternary ammonium compounds

Delivery forms: 1 litre bottle with flip top cap, 5 litre canister

Complies with EN 13624, EN 13697, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563, EN 16615 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	60 sec
Antibiotic-resistant bacteria	60 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	60 sec
Antibiotic-resistant yeasts	60 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	30 sec
Poliovirus	60 sec



DENTIRO® Light

Fully virucidal, ready-to-use surface disinfectant for the quick and residue-free disinfection of non-invasive medical devices (medical equipment, installations, and instruments). Perfume-free.

- Kills 101 microorganisms
- Active in 60 seconds
- Compatible with more than 275 materials including artificial leather
- Leaves no residues or smearing

Composition in 100 g: 31.3 g ethanol, 16.7 g isopropanol

Delivery forms: 1 litre bottle with flip top cap, 2 litre bottle, 5 litre canister

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	60 sec
Poliovirus	60 sec
Rotavirus	30 sec



ISORAPID® Liquid

Ready-to-use, aldehyde-free solution for the quick disinfection and cleaning of surfaces of non-invasive medical devices (medical equipment, installations and instruments). Active against norovirus in 30 seconds.

- Kills 101 microorganisms
- Broad effective spectrum
- Low alcohol content of 50%
- Excellent material compatibility
- Tested under high organic load

Composition in 100 g: 50 g 1-propanol, 0.057 g aminoalkylglycine

Delivery forms: 5 litre canister

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Fungicidal	2 min
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	30 sec
Poliovirus	60 sec
Polyomavirus	30 sec



ISORAPID® Spray

Ready-to-use, aldehyde-free spray for the quick disinfection and cleaning of surfaces of non-invasive medical devices (medical equipment, installations, and instruments).

- Kills 101 microorganisms
- Active in 60 seconds
- Compatible with 250 materials including artificial leather
- Contains less than 50% alcohol

Composition in 100 g: 20 g ethanol, 28 g 1-propanol, 0.056 g quaternary ammonium compounds

Delivery forms: 1 litre bottle with flip top cap, 2 litre bottle, 5 litre canister

Complies with EN 13624, EN 13697, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563, EN 16615 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	60 sec
Antibiotic-resistant bacteria	60 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	60 sec
Antibiotic-resistant yeasts	60 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	30 sec
Poliovirus	60 sec



DENTIRO® Zero

Alcohol-free, ready-to-use surface disinfectant for the quick disinfection and cleaning of non-invasive medical devices (medical equipment, installations, and instruments).

- Kills 101 microorganisms
- Short exposure time
- Alcohol-free
- Compatible with more than 250 materials including artificial leather

Composition in 100 g: 0.4 g benzalkonium chloride, 0.4 g aminoalkylglycine

Delivery forms: 1 litre bottle with flip top cap, 5 litre canister

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14563 and EN 16615.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 min
Mycobactericidal	60 min
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 min
Adenovirus	30 min
Norovirus	60 min
Poliovirus	60 min
Rotavirus	30 sec



DENTIRO® Foam

Alcohol-free, ready-to-use disinfectant foam for the disinfection, cleaning, and care of non-invasive medical devices with sensitive surfaces such as gurneys and medical chairs.

- Kills 101 microorganisms
- Short exposure time
- Compatible with more than 250 materials including artificial leather
- Tested under high organic load

Composition in 100 g: 0.4 g benzalkonium chloride, 0.4 g aminoalkylglycine

Delivery forms: 200 ml foamer bottle, 2 litre bottle

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14563 and EN 16615.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 min
Mycobactericidal	60 min
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 min
Adenovirus	30 min
Norovirus	60 min
Poliovirus	60 min
Rotavirus	30 sec

DENTIRO® Wipes

Kill 101 microorganisms

active in
60
seconds

Candida albicans

Hepatitis C
virus

Mycobacterium bovis

Yersinia
enterocolitica

Klebsiella pneumoniae

Burkholderia cepacia

Enterococcus
hirae

Mycobacterium
terrae

Proteus mirabilis

Bordetella pertussis

Influenza A virus

Mycobacterium
avium

Neisseria
meningitidis

Streptococcus pyogenes

Adenovirus

Mycobacterium
smegmatis

Proteus
vulgaris

Moraxella catarrhalis

Escherichia
coli

Listeria monocytogenes

Poliovirus

Saccharomyces
cerevisiae

Staphylococcus epidermidis

Salmonella
typhimurium

Bovine Viral Diarrhea Virus

Vibrio
parahaemolyticus

MRSA

Propionibacterium
acnes

Vancomycin-resistant
Enterococcus

Acinetobacter baumannii

Rotavirus

Clostridium difficile

Campylobacter
jejuni

Coronavirus

Vaccinia virus

Mycobacterium chelonae

HIV

Staphylococcus
aureus

Herpes simplex virus

Norovirus

Hepatitis B virus

Trichoderma viride

Enterobacter
aerogenes



CONTACT US FOR A COMPLETE LIST
TEL. +41 (0)44 226 44 44
EMAIL: INFO@OROCLEAN.COM

SURFACE DISINFECTANT WIPES



DENTIRO® Wipes Extra Comfort

Highly saturated disinfectant wipes with a unique texture for a residue-free disinfection and cleaning of non-invasive medical devices (medical equipment, installations, and instruments). Bactericidal, tuberculocidal, mycobactericidal, yeasticidal, and fully virucidal in 60 seconds.

- Kills 101 microorganisms
- Active in 60 seconds
- Compatible with more than 275 materials including artificial leather
- Residue-free drying

Composition of the moisturizing solution in 100 g:

31.3 g ethanol, 16.7 g isopropanol

Delivery forms: Refill pouch with 165 disinfectant wipes (14.5 x 18 cm), lemon scent

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines. Dermatologically tested.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	60 sec
Poliovirus	60 sec
Rotavirus	30 sec



DENTIRO® Wipes Smart Saver

Cheap disinfectant wipes for the disinfection and cleaning of non-invasive medical devices (medical equipment, installations, and instruments). Bactericidal, tuberculocidal, mycobactericidal, yeasticidal, and fully virucidal in 60 seconds.

- Kills 101 microorganisms
- Active in 60 seconds
- Compatible with more than 275 materials including artificial leather
- Highly saturated wipes

Composition of the moisturizing solution in 100 g:

31.3 g ethanol, 16.7 g isopropanol

Delivery forms: Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), lemon scent, Refill pouch with 250 disinfectant wipes (14.5 x 20 cm), lemon scent

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines. Dermatologically tested.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	60 sec
Poliovirus	60 sec
Rotavirus	30 sec



DENTIRO® Wipes Ultra Soft

Ultra soft disinfectant wipes with maximized strength for the simultaneous disinfection, cleaning, and deodorizing of non-invasive medical devices (medical equipment, installations, and instruments). Bactericidal, tuberculocidal, mycobactericidal, yeasticidal, and fully virucidal in 60 seconds.

- Kills 101 microorganisms
- Active in 60 seconds
- Compatible with more than 275 materials including artificial leather
- Soft yet highly tear-resistant wipes

Composition of the moisturizing solution in 100 g:

31.3 g ethanol, 16.7 g isopropanol

Delivery form: Refill pouch with 120 disinfectant wipes (19.3 x 21 cm), lemon scent

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines. Dermatologically tested.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	60 sec
Poliovirus	60 sec
Rotavirus	30 sec



DENTIRO® Wipes

Fully virucidal disinfectant wipes for the disinfection of non-invasive medical devices (medical equipment, installations, and instruments).

- Kills 101 microorganisms
- Active in 60 seconds
- Compatible with more than 275 materials including artificial leather
- Residue-free drying

Composition of the moisturizing solution in 100 g:

31.3 g ethanol, 16.7 g isopropanol

Delivery form: Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), classic scent; Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), lemon scent; Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), floral scent

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines. Dermatologically tested.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	60 sec
Poliovirus	60 sec
Rotavirus	30 sec



DENTIRO® Jumbo Wipes

Fully virucidal disinfectant wipes for the disinfection of non-invasive medical devices (medical equipment, installations, and instruments).

- Kills 101 microorganisms
- Active in 60 seconds
- Compatible with more than 275 materials including artificial leather
- Residue-free drying

Composition of the moisturizing solution in 100 g:

31.3 g ethanol, 16.7 g isopropanol


Delivery form: Refill pouch with 100 disinfectant wipes (14.5 x 40 cm), classic scent; Refill pouch with 100 disinfectant wipes (14.5 x 40 cm), lemon scent

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines. Dermatologically tested.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	60 sec
Poliovirus	60 sec
Rotavirus	30 sec

11 reasons why you should use DENTIRO® Wipes




- 1 Kills 101 microorganisms
- 2 Active in 60 seconds
- 3 Doesn't foam, smear or leave any residue
- 4 Compatible with more than 275 materials including artificial leather
- 5 Dermatologically tested
- 6 Highly saturated wipes
- 7 Pleasant lemon scent
- 8 Available in different sizes
- 9 Available in refill pouches
- 10 Convenient and durable tub with flip lid
- 11 Silicone cross-slit opening prevents wipes from drying out

DENTIRO® Wipes

Compatible with more than 275 MATERIALS including artificial leather

We have tested the stability of artificial leather according to the EN ISO 21530 norm. This test takes two weeks and simulates 20,000 disinfections with one minute contact time. Cotton balls are soaked with disinfectant and placed on the leather while being changed everyday for two weeks. The test is passed when no changes on the colour of the surface are visible.



Material: KaVo®
Time: 2 weeks (20,000 reprocessing cycles)
Result: No change.
Comparable to water.



Material: Sirona®
Time: 2 weeks (20,000 reprocessing cycles)
Result: No change.
Comparable to water.



Material: A-dec®
Time: 2 weeks (20,000 reprocessing cycles)
Result: No change.
Comparable to water.




Material: Planmeca®
Time: 2 weeks (20,000 reprocessing cycles)
Result: No change.
Comparable to water.



Material: Castellini®
Time: 2 weeks (20,000 reprocessing cycles)
Result: No change.
Comparable to water.



Material: Anthos®
Time: 2 weeks (20,000 reprocessing cycles)
Result: No change.
Comparable to water.



Material: Victor®
Time: 2 weeks (20,000 reprocessing cycles)
Result: No change.
Comparable to water.



Material: Stern Weber®
Time: 2 weeks (20,000 reprocessing cycles)
Result: No change.
Comparable to water.

KAVO® is a registered trademark of Kaltenbach & Voigt GmbH.
Sirona® is a registered trademark of Sirona Dental Systems GmbH.
A-dec® is a registered trademark of A-dec Inc.
Planmeca® is a registered trademark of Planmeca Group.
Castellini®, Anthos®, Victor® and Stern Weber® are registered trademarks of Cefla Group.

CONTACT US FOR A COMPLETE LIST
TEL. +41 (0)44 226 44 44
EMAIL: INFO@OROCLEAN.COM

SURFACE DISINFECTANT CONCENTRATES



ISORAPID® Floor

Aldehyde- and phenol-free concentrate for the disinfection and cleaning of large water-washable surfaces of non-invasive medical devices (medical appliances and installations). Also suitable for immersion disinfection.

- Kills 101 microorganisms
- Powerful cleanser
- High-yielding concentrate
- Long-lasting antimicrobial action on treated surfaces
- Good material compatibility
- Simple application

Composition in 100 g: 13.9 g N,N-bis(3-aminopropyl)dodecylamine, anionic surfactants, non-ionic surfactants, cleaning booster, fragrances, dyes

Delivery forms: 1 litre bottle with measuring cap, 5 litre canister

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	0.5% in 10 min
Antibiotic-resistant bacteria	0.5% in 10 min
Tuberculocidal	0.5% in 30 min
	2% in 15 min
Yeasticidal	0.5% in 5 min
Antibiotic-resistant yeasts	0.5% in 5 min
Enveloped viruses	1% in 5 min
BVDV	0.25% in 5 min
Coronavirus	1% in 5 min
Hepatitis B virus	1% in 5 min
Hepatitis C virus	0.25% in 5 min
Herpes simplex virus	1% in 5 min
HIV	1% in 5 min
Influenza A virus	1% in 1 min
Vaccinia virus	1% in 5 min
Adenovirus	1% in 60 min

NEW



ISORAPID® OP-Forte AF

Mild disinfectant concentrate for sensitive water-washable surfaces and non-invasive medical devices such as bedpans and urine sample containers.

- Kills 101 microorganisms
- Effective against bacteria, yeasts and viruses
- Better value for money
- Excellent material compatibility
- Simple application
- Biodegradable

Composition in 100 g: 10 g quaternary ammonium compounds, perfumes, dyes

Delivery forms: 2 litre bottle with measuring cap, 5 litre canister

Complies with EN 13624, EN 13727, EN 14476 and EN 16615.

EFFECTIVE SPECTRUM

Bactericidal	1% in 1 min
Antibiotic-resistant bacteria	1% in 1 min
Yeasticidal	1% in 1 min
Antibiotic-resistant yeasts	1% in 1 min
Enveloped viruses	1% in 5 min
BVDV	1% in 5 min
Coronavirus	1% in 5 min
Hepatitis B virus	1% in 5 min
Hepatitis C virus	1% in 5 min
Herpes simplex virus	1% in 5 min
HIV	1% in 5 min
Influenza A virus	1% in 5 min
Vaccinia virus	1% in 5 min
Adenovirus	2% in 1 min
Norovirus	2% in 5 min

The stronger and more powerful ISORAPID® OP-Forte AF

OROLIN® OP-Forte AF

1% in 60 min

1% in 60 min

5% in 1 min

N/A

5% in 1 min

1% in 15 min

N/A

N/A

N/A

OLD

ISORAPID® OP-Forte AF

1% in 1 min

1% in 1 min

1% in 5 min

1% in 5 min

1% in 5 min

1% in 5 min

1% in 5 min

1% in 5 min

1% in 5 min

NEW





NEW



OROLIN® Multisept Plus

Powerful concentrate for the rapid disinfection and cleaning of medical and surgical instruments.

- Kills 101 microorganisms
- Fully virucidal within 15 minutes
- 2 litre of concentrate yields 100 litres of working solution
- Powerful cleanser
- Excellent material compatibility
- Suitable for use in ultrasonic cleaning devices

Composition in 100 g: 3.6 g didecylmethylammonium chloride, 9.1 g alkylamine, cleaning booster, auxiliaries

Delivery forms: 1 litre bottle with measuring cap, 2 litre bottle with measuring cap, and 5 litre canister

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	2% in 5 min
Antibiotic-resistant bacteria	2% in 5 min
Tuberculocidal	2% in 5 min
Mycobactericidal	2% in 5 min
Yeasticidal	2% in 5 min
Antibiotic-resistant yeasts	2% in 5 min
Enveloped viruses	2% in 5 min
BVDV	2% in 5 min
Coronavirus	2% in 5 min
Hepatitis B virus	2% in 5 min
Hepatitis C virus	2% in 5 min
Herpes simplex virus	2% in 5 min
HIV	2% in 5 min
Influenza A virus	2% in 1 min
Vaccinia virus	2% in 5 min
Non-enveloped viruses	2% in 15 min
Adenovirus	2% in 5 min
Norovirus	2% in 10 min
Poliovirus	2% in 15 min



ULTRACLEAN-3

Alkaline-based concentrate for the cleaning of surgical and dental instruments. Particularly suitable for use in ultrasonic cleaning devices.

- Powerful cleanser
- Suitable for use in ultrasonic cleaning devices
- Simple application

Composition in 100 g: 16 g potassium hydroxide, complexing agent, special surfactants

Delivery form: 1 litre bottle with measuring cap

NEW



OROCID-MULTISEPT® Plus has been replaced with **OROLIN® Multisept Plus**, a more powerful formulation for the same price

OROCID-MULTISEPT® Plus (Old formulation)

Active in 60 minutes at 2%
Limited virucidal
Tuberculocidal
With corrosion inhibitors

OROLIN® Multisept Plus (Improved formulation)

Active in only 15 minutes at 2%
Fully virucidal
Tuberculocidal & mycobactericidal
Even less corrosive than water



OROLIN® Burbath

Ready-to-use, aldehyde-free solution for the disinfection and cleaning of medical and dental rotating instruments.

- Kills 101 microorganisms
- Fully virucidal within 30 seconds
- Highly effective even in the presence of blood and proteins
- Contains corrosion inhibitors
- Time-saving application

Composition in 100 g: 32 g isopropanol, 0.7 g phenoxyethanol, 0.2 g quaternary ammonium compounds, alkali, corrosion inhibitors, surfactants

Delivery forms: 1 litre bottle, 2 litre bottle, 5 litre canister

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	30 sec
Antibiotic-resistant bacteria	30 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	30 sec
Antibiotic-resistant yeasts	30 sec
Fungicidal	60 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	30 sec
Adenovirus	30 sec
Norovirus	30 sec
Poliovirus	30 sec
Polyomavirus	30 sec



OROLIN® Burbath Container

Soaking tray for the immersion disinfection of rotating instruments.

- Effective spectrum and exposure time are indicated on the lid
- White plastic with a fine sieve and lid with insertion hole and plug
- Ideal for the disinfection of small objects such as rotating instruments



OROLIN® Perasept

Concentrated powder for the disinfection, cold sterilization, and cleaning of medical and dental instruments.

- Kills 101 microorganisms
- Cold sterilization
- Tested under high organic load
- Aldehyde-free
- No discoloration or fixation of blood
- Free of carcinogenic substances

Composition in 100 g: 45 g sodium percarbonate, 25 g tetraacetylenediamine, 24 g citric acid, surfactants, corrosion inhibitors.

Delivery form: 1 kg package with measuring spoon, 5 kg package with measuring spoon

Complies with EN 13624, EN 13704, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562 and EN 14563.

EFFECTIVE SPECTRUM

Bactericidal	1% in 5 min 2% in 30 sec
Antibiotic-resistant bacteria	1% in 10 min 2% in 1 min
Tuberculocidal, Mycobactericidal	1% in 15 min 2% in 5 min
Yeasticidal	1% in 10 min 2% in 1 min
Antibiotic-resistant yeasts	1% in 10 min 2% in 1 min
Fungicidal	1% in 30 min 2% in 10 min
Enveloped viruses	1% in 5 min 2% in 1 min
BVDV, Coronavirus	1% in 5 min 2% in 1 min
HBV, HCV, HSV, HIV	1% in 5 min 2% in 1 min
Influenza A virus	1% in 1 min 2% in 30 sec
Vaccinia virus	1% in 5 min 2% in 1 min
Non-enveloped viruses	1% in 5 min 2% in 1 min
Adenovirus, Norovirus	1% in 30 sec 2% in 1 min
Bovine enterovirus, Poliovirus	1% in 5 min 2% in 1 min
Sporicidal	1% in 5 min 2% in 1 min



1 litre

OROCLEAN® Soaking Tray

OROCLEAN® Soaking Tray for the immersion disinfection of medical instruments.

- Requires minimal space
- Simple application
- Recessed metal plate in lid perfect for fitting a timer
- Removable strainer for easy rinsing



3 litre



OROCLEAN® Dosing Pump 20 ml

Dosing pump with 20 ml output for 1 litre and 2 litre bottles.



Scan
QR code to
view video!

OROCLEAN® Hygiene Club

Want infection control tips?
Join **OROCLEAN®** hygiene club!

OROCLEAN® Hygiene Club is not just about providing information. It's about creating a dialogue between you and OCC Switzerland.

Become a member now!

Not yet a Member? [Register for FREE](#) ↓

Salutation <input type="radio"/> Mr <input type="radio"/> Ms <input type="radio"/> Mrs	Title <input type="text"/>
First name*	Last name*
Email*	Confirm Email*
Clinic /Surgery /Company	Street *
Street 2	City*
Country*	Industry*

[Register for FREE](#)



Scan
QR code to
view video!



OROLIN® Intra

Fully virucidal, ready-to-use suspension for the disinfection, cleaning, and lubrication of dental handpieces, contra-angles, and turbines.

- 3-in-1
- Maximum cost-effectiveness
- Fully virucidal within 30 seconds
- Tested under high organic load

Composition in 100 g: 53.7 g ethanol, 19.0 g 1-propanol, auxiliaries, medical oil

Delivery form: 575 ml bottle with trigger pump

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	60 sec
Antibiotic-resistant bacteria	60 sec
Tuberculocidal	5 min
Mycobactericidal	5 min
Yeasticidal	60 sec
Antibiotic-resistant yeasts	60 sec
Fungicidal	2 min
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	30 sec
Adenovirus	30 sec
Norovirus	30 sec
Poliovirus	30 sec
Polyomavirus	30 sec

NUMEROUS ADAPTERS GUARANTEE FULL COMPATIBILITY



OROLIN® Adapter
(NSK Mach)



OROLIN® Adapter
(KAVO® Multiflex®)



OROLIN® Adapter
(Midwest® ISO-C)
4 or 6 holes



OROLIN® Adapter
(NSK QD-J®)



OROLIN® Adapter
(Siemens® 4-holes)



OROLIN® Adapter
(BienAir Unifix®)



OROLIN® Adapter
(Borden® ISO-A)
2-3 holes



OROLIN® Adapter
(Sino®)



OROLIN® Adapter
(Sirona® A,B,F,R)



OROLIN® Adapter
(W&H® Roto Quick®)



OROLIN® Adapter
(Plastic ISO 3964)



DISINFECTION, CLEANING, AND LUBRICATION WITH OROLIN® INTRA PROCESS YOUR INSTRUMENTS IN 5 MINUTES



① Connect the OROLIN® Adapter to the OROLIN® Intra bottle.



② Fully insert the OROLIN® Adapter into the connection end of your handpiece or turbine.



③ Hold the OROLIN® Intra bottle upside down. Activate the trigger until the liquid is coming out from the instrument head.



④ Wait for 5 minutes to allow the disinfectant to act.







⑤ Hold the OROLIN® Intra bottle upright and cover the instrument with a single-use cloth. Activate the trigger several times to flush out the solution and to blow-dry the instrument.



⑥ Use the single-use cloth for the outer disinfection of your handpiece or turbine.

AEROSOL-FREE APPLICATION. COMPLETE SURFACE COVERAGE.



Complete Coverage – OROLIN® Intra is the only product on the market that completely fills your handpiece, contra-angle or turbine with the solution, thereby, ensuring that all surfaces are covered with the disinfectant for the entire exposure time. There will never be any gaps.

No Further Processing – Your handpiece, contra-angle or turbine is ready for use right after cleaning, disinfection, and lubrication. You do not need to sterilize your instrument if used for a non-invasive procedure.

Scan QR code to know more.





NEW



ORO CLEAN® Plus

Broad spectrum and non-foaming 3-in-1 concentrate for the disinfection, cleaning and deodorisation of dental suction units.

- Kills 101 microorganisms
- Amazingly fast-acting
- 3-in-1 solution
- Non-foaming formula
- 2 litres yield 50 applications
- Fully virucidal

Composition in 100 g: 2.475 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine, 7.5 g quaternary ammonium compound, anti-foaming agents, special surfactants

Delivery forms: 2 litre bottle with measuring cap, 5 litre canister

Complies with EN 13624, EN 13727, EN 14476, EN 14561, EN 14562 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	2% in 5 min
Antibiotic-resistant bacteria	2% in 5 min
Yeasticidal	2% in 5 min
Antibiotic-resistant yeasts	2% in 5 min
Enveloped viruses	2% in 5 min
BVDV	2% in 1 min
Coronavirus	2% in 5 min
Hepatitis B virus	2% in 5 min
Hepatitis C virus	2% in 1 min
Herpes simplex virus	2% in 5 min
HIV	2% in 5 min
Influenza A virus	2% in 1 min
Vaccinia virus	2% in 5 min
Non-enveloped viruses	2% in 60 min
Adenovirus	2% in 30 min
Norovirus	2% in 60 min
Poliovirus	2% in 30 min

NEW



ORO CLEAN® Liquid

Powerful non-foaming concentrate for the cleaning and deodorisation of dental suction units.

- Excellent cleaning power
- Non-foaming formula
- Excellent material compatibility
- 2 litres yield 50 applications
- Great value for money
- Simple application

Composition in 100 g: alkaline cleaning compounds, anti-foaming agents, special surfactants

Delivery forms: 2 litre bottle with measuring cap, 5 litre canister



OROCLEAN® Cup

Aspirating cup for the quick and easy disinfection as well as cleaning of suction units.

- Simple application
- Closed system for your safety
- 2 litre capacity



ASEPTOPRINT® Liquid

Aldehyde-free concentrate for the disinfection and cleaning of dental impressions.

- Kills 101 microorganism
- Broad effective spectrum
- Short exposure time
- Excellent material compatibility with alginate, polyether, polysulfide, A- and C-silicone
- Aldehyde-, phenol- and chlorine-free
- Simple application

Composition in 100 g: 30.0 g benzalkonium chloride, surfactants, cleaning booster, auxiliaries

Delivery form: 1 litre bottle with measuring cap

Complies with EN 13624, EN 13727, EN 14476, EN 14561, EN 14562 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	1% in 2 min
Antibiotic-resistant bacteria	3% in 1 min
Yeasticidal	1% in 2 min
Antibiotic-resistant yeasts	3% in 1 min
Enveloped viruses	1% in 1 min
BVDV	1% in 1 min
Coronavirus	1% in 1 min
Hepatitis B virus	1% in 1 min
Hepatitis C virus	1% in 1 min
Herpes simplex virus	1% in 1 min
HIV	1% in 1 min
Influenza A virus	1% in 1 min
Vaccinia virus	1% in 1 min
Adenovirus	1% in 2 min
Norovirus	1% in 2 min
Rotavirus	3% in 1 min



1 litre

ASEPTOPRINT® Liquid Soaking Tray

ASEPTOPRINT® Liquid soaking tray for the immersion disinfection of dental impressions.

- Requires minimal space
- Simple application
- Recessed metal plate in lid perfect for fitting a timer
- Removable strainer for easy rinsing



OROCLEAN® Dosing Pump 20 ml

Dosing pump with 20 ml output for 1 litre, 2 litre and 2.5 litre bottles.



ASEPTOPRINT® Spray

Ready-to-use, aldehyde-free disinfectant spray for the quick disinfection of dental impressions made of polyether, polysulfide, alginate, as well as A- and C-silicone.

- Kills 101 microorganisms
- Short exposure time
- Excellent material compatibility
- Biodegradable

Composition in 100 g: 20 g ethanol, 28 g 1-propanol, 0.056 g quaternary ammonium compounds

Delivery forms: 1 litre bottle with sprayer head

Complies with EN 13624, EN 13697, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563 and RKI/DVV guidelines.

EFFECTIVE SPECTRUM

Bactericidal	60 sec
Antibiotic-resistant bacteria	60 sec
Tuberculocidal	60 sec
Mycobactericidal	60 sec
Yeasticidal	60 sec
Antibiotic-resistant yeasts	60 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus	60 sec
Norovirus	30 sec
Poliovirus	60 sec



OROMED® Gel

Clear hydroalcoholic gel with low viscosity for surgical and hygienic hand disinfection.

- Kills 101 microorganisms
- Broad effective spectrum
- Free of fragrance and colorants
- Excellent skin compatibility
- Contains anti-aging ingredients
- Dermatologically tested

Composition in 100 g: 52.8 g ethanol, 19.0 g 1-propanol, 2.8 g isopropanol, moisturizers, emollients, vitamin E acetate, panthenol, bisabolol, gelling agents

Delivery forms: 100 ml bottle with flip top cap, 500 ml bottle with pump, 500 ml bottle (EURO), 1 litre bottle (EURO), 1 litre refill cartridge for myArtisto dispenser, 5 litre canister

Complies with EN 1500, EN 12791, EN 13624, EN 13727, EN 14348, EN 14476 and RKI/DVV guidelines. Dermatologically tested.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	15 sec
Mycobactericidal	30 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Fungicidal	60 sec
Enveloped viruses	15 sec
BVDV, Coronavirus	15 sec
HBV, HCV, HSV, HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	5 min
Adenovirus, Norovirus	60 sec
Poliovirus	5 min
Rotavirus	15 sec
Biocidal product authorisation number: CHZN1698	

APPLICATION AREA

Hygienic hand disinfection EN 1500	30 sec
Surgical hand disinfection EN 12791	90 sec

NEW



OROMED® Liquid

Fast-acting alcohol-based hand rub for hygienic and surgical disinfection. Bactericidal, tuberculocidal, mycobactericidal, yeasticidal, and fully virucidal.

- Kills 101 microorganisms
- Free of fragrance and colorants
- Excellent skin compatibility
- Contains anti-aging ingredients
- Dermatologically tested

Composition in 100 g: 52.8 g ethanol, 19.0 g 1-propanol, 2.8 g isopropanol, moisturizers, emollients, vitamin E acetate, panthenol, bisabolol

Delivery forms: 100 ml bottle with flip top cap, 500 ml bottle with pump, 500 ml bottle (EURO), 1 litre bottle (EURO), 1 litre refill cartridge for myArtisto dispenser, 5 litre canister

Complies with EN 1500, EN 12791, EN 13624, EN 13727, EN 14348, EN 14476.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	15 sec
Mycobactericidal	15 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV, Coronavirus	15 sec
HBV, HCV, HSV, HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Non-enveloped viruses	60 sec
Adenovirus, Norovirus	30 sec
Poliovirus	60 sec
Biocidal product authorisation number: CHZN4499	

APPLICATION AREA

Hygienic hand disinfection EN 1500	30 sec
Surgical hand disinfection EN 12791	90 sec

WANT TO KNOW THE COLD, HARD TRUTH ABOUT GERMS?

Like our Facebook page
@Oroclean and follow us

Scan here



When it comes to infections and diseases, our posts bare it all so you don't become a statistic.





OROSEPT® Solution

Alcohol-based hand rub for hygienic and surgical hand and forearm disinfection.

- Kills 101 microorganisms
- Short exposure time
- Contains skincare, moisturizing, and nourishing ingredients
- Excellent skin compatibility
- Pleasant lemon fragrance

Composition in 100 g: 44.9 g isopropanol, 30.0 g 1-propanol, 0.15 g dialkyldimethylammonium chloride, water, natural moisturizing factors, perfume

Delivery forms: 100 ml bottle with flip top cap, 500 ml bottle with pump, 500 ml bottle (EURO), 1 litre bottle (EURO), 1 litre refill cartridge for myArtisto dispenser, 2 litre bottle, 5 litre canister

Complies with EN 1500, EN 12791, EN 13624, EN 13727, EN 14348 and EN 14476. Dermatologically tested.

EFFECTIVE SPECTRUM

Bactericidal	15 sec
Antibiotic-resistant bacteria	15 sec
Tuberculocidal	30 sec
Mycobactericidal	30 sec
Yeasticidal	15 sec
Antibiotic-resistant yeasts	15 sec
Enveloped viruses	15 sec
BVDV	15 sec
Coronavirus	15 sec
Hepatitis B virus	15 sec
Hepatitis C virus	15 sec
Herpes simplex virus	15 sec
HIV	15 sec
Influenza A virus	15 sec
Vaccinia virus	15 sec
Norovirus	60 sec
Biocidal product authorisation number: CHZB0364	

APPLICATION AREA

Hygienic hand disinfection EN 1500	30 sec
Surgical hand disinfection EN 12791	90 sec



OROMED® Hand Wash

Alkali-free hand washing liquid with high-quality moisturizing agents.

- Mild cleansing lotion
- Contains skincare substances
- Maximum cost-effectiveness

Ingredients (INCI): AQUA, SODIUM LAURETH SULFATE, COCAMIDOPROPYL BETAINE, DECYL GLUCOSIDE, GLYCERIN, SODIUM CHLORIDE, COCO-GLUCOSIDE, GLYCERYL OLEATE, SODIUM BENZOATE, CITRIC ACID, LEONTOPODIUM ALPINUM EXTRACT, PARFUM, C.I. 14720

Delivery forms: 500 ml bottle with pump, 500 ml bottle (EURO), 1 litre bottle (EURO), 1 litre refill cartridge for myArtisto dispenser, 5 litre canister

Dermatologically tested.



OROCREME

Hand cream with protective and nourishing agents for skin care. Easily absorbed and non-greasy.

- Contains skin-moisturizing ingredients
- Excellent skin compatibility
- Non-greasy

Ingredients (INCI): AQUA, DECYL OLEATE, CAPRYLIC/CAPRIC TRIGLYCERIDE, CETEARYL ALCOHOL, GLYCERIN, GLYCERYL STEARATE SE, PHENOXYETHANOL, PARFUM, CITRIC ACID

Delivery form: 230 ml dispenser bottle with pump

HOW TO DISINFECT YOUR HANDS

Apply the disinfectant onto dry hands and rub in for 30 seconds according to the following six steps. It is important that your hands remain moist during the entire application time.



Rub palm to palm.



Run right palm over the back of the left hand - and vice versa.



Rub palm to palm with interlaced spread fingers.



Rub the backs of the interlocked fingers on the opposing palms.



Rub the right thumb in the closed palm of the left hand using a circular motion - and vice versa.



Rub the closed fingertips of the right hand on the left palm using a circular motion - and vice versa.

CLEAN HANDS SAVE LIVES



|O|C|C|

Switzerland

Oro Clean Chemie AG
Allmendstrasse 21
8320 Fehraltorf
Switzerland

Tel. +41 (0)44 226 44 44
Fax +41 (0)44 226 44 00
info@oroclean.com

www.oroclean.com



OROCLEAN® myArtisto MS-3 Dispenser

The OROCLEAN® myArtisto MS-3 Sanitizer Dispenser is a push-button soap and hand disinfectant dispenser ideal for use in health care settings like hospitals and clinics as well as kitchens and restaurants, day-care centres, fitness gyms, or public buildings. The dispenser features three extra-large windows that allow users to customize the appearance simply by replacing the graphics. The sanitizer dispenser is compatible with alcohol-based hand disinfectants, hand sanitizers, liquid hand soaps, as well as hand creams. The OROCLEAN® myArtisto MS-3 Sanitizer Dispenser can be effortlessly wall-mounted by using the included double-sided adhesive tape or wall brackets with the corresponding screws and dowels. The unit has a vandal resistant security lock and key. An optional arm/elbow lever is available for use in areas where cross contamination may be of concern.



OROCLEAN® myArtisto AS-3 Dispenser

The OROCLEAN® myArtisto AS-3 Sanitizer Dispenser is a touch-free soap and hand disinfectant dispenser ideal for use in health care settings like hospitals and clinics where the risk of cross-infections is particularly high. The dispenser features three extra-large windows that allow users to customize the appearance simply by replacing the graphics. The device automatically dispenses 1.5 ml or 3 ml of liquid depending on the setting. The sanitizer dispenser is compatible with hand disinfectants, hand sanitizers, liquid hand soaps, as well as hand creams. The dispenser is powered by four AA batteries that last for more than 10,000 cycles (batteries are not included in the delivery). A power adapter is optionally available. The OROCLEAN® myArtisto AS-3 Sanitizer Dispenser can be effortlessly wall-mounted by using the included double-sided adhesive tape or wall brackets. The unit has a vandal resistant security lock and key.



OROCLEAN® myArtisto Drip Tray

The myArtisto Drip Tray is designed to protect your floor and walls by keeping them free of excess sanitizer. The clip-on drip tray easily attaches to the bottom of any myArtisto Sanitizer Dispenser without the need for tools.



OROCLEAN® myArtisto Elbow Lever

The arm/elbow lever for myArtisto MS-3 Sanitizer Dispenser converts the push-button dispenser into a lever-operated dispenser. The lever is made of stainless steel and can easily be attached without the need for tools.

Medical uniform disinfection made as easy as 1-2-3

NEW




Step 1

Place laundry in
the machine



Step 2

Add 100 g of **OROCLEAN® Splendid**
and laundry detergent



Step 3

Set temperature to 40°C
and wash as usual

That's it!

Your medical clothes are now free of
bacteria, fungi, viruses and spores.



Available in
5 kg and
1 kg



Make a change and simplify your life.
Visit www.orocean.com for more details.

Use disinfectants safely. Always read the label and product information before use

Fully virucidal in 60 seconds Nourishes skin for hours

OROMED® Liquid

- ✓ Contains anti-aging ingredients
- ✓ Free of fragrance and colorants
- ✓ Excellent skin compatibility



Call your dealer today to order or
visit www.orocean.com for more information.

Use disinfectants safely. Always read the label and product information before use

Fully
Virucidal





OROCLEAN® Latex Gloves

Very thin, tear-resistant, and powder-free latex examination gloves for unprecedented comfort and precise handling. Fit either hand. Non-sterile.

- Low-allergenic ingredients
- Optimal wearing comfort
- Excellent tactile sensitivity
- Powder-free

Delivery form: Box of 100 gloves
Glove sizes: XS, S, M, L



OROCLEAN® Nitrile Gloves

Ergonomically engineered, powder- and latex-free nitrile examination gloves for maximum protection against germs, punctures, and chemicals. Non-sterile.

- Advanced barrier protection
- Optimum grip
- Powder-free
- 100% latex-free

Delivery form: Box of 100 gloves
Glove sizes: S, M, L, XL **NEW**



OROCLEAN® Face Masks

Comfortably fitting surgical earloop face masks for single-use. Free of latex, lint, fibreglass, and gluten. Non-irritating.

- Excellent wearing comfort
- Padded nosepiece
- Low breathing resistance
- Enriched with Aloe Vera

Delivery form: Box of 50 masks



OROCLEAN® Splendid

Laundry detergent additive for the disinfection of all kinds of textiles (except for delicate and surgical textiles) in medical or public institutions such as medical and dental practices or beauty salons as well as in private sector.

- Kills 101 microorganisms
- Broad effective spectrum
- Biodegradable
- Complies with European standards

Composition in 100 g: 30 – 50 % oxygen-based bleaching agent, < 5 % non-ionic surfactants, phosphonates. Other components: TAED, citric acid

Delivery form: 1 kg package with measuring spoon, **NEW**
5 kg package with measuring spoon


Complies with EN 13624, EN 13704, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562 and EN 14563.

EFFECTIVE SPECTRUM


Bactericidal	30 min
Tuberculocidal	30 min
Mycobactericidal	30 min
Yeasticidal	30 min
Fungicidal	30 min
Virucidal	30 min
Sporicidal	30 min




MICROBIAL CONTAMINATION OF PRACTICE CLOTHING



Dental clinic in Bern, Switzerland
2015




Dental clinic in Zurich, Switzerland
2014



Dental clinic in Wintherthur, Switzerland
2014



Dental clinic in Basel, Switzerland
2015



Dental clinic in Aarau, Switzerland
2014



Dental clinic in St. Gallen, Switzerland
2014



DENTIRO® Wipes Kill 101 microorganisms

active in
60
seconds



Our best disinfectant wipes ever!

- Compatible with more than 275 materials including artificial leather




DENTIRO® Wipes
14.5 cm x 20 cm, 120 pcs
14.5 cm x 40 cm, 100 pcs



**DENTIRO® Wipes
Extra Comfort**
14.5 cm x 18 cm, 165 pcs



**DENTIRO® Wipes
Ultra Soft**
19.3 cm x 21 cm, 120 pcs



**DENTIRO® Wipes
Smart Saver**
14.5 cm x 20 cm, 120 pcs
14.5 cm x 20 cm, 250 pcs