


LM-LiftOut luxating instruments

Angulations and sizes



LM-TwistOut elevators

Angulations and sizes



LM-Instruments Oy
PL 88 (Norrbyn rantatie 8)
FI-21601 Parainen
Tel. +358 2 4546 400
Fax +358 2 4546 444
info@lminstruments.com
www.lminstruments.com

Contact your local dealer

LM LiftOut & LM TwistOut



LM

feel the
difference

Extract with confidence

LM LiftOut

Luxating instruments for atraumatic extraction

LM-LiftOut instruments allow tooth extraction to be done as atraumatically as possible. This is important in order to enable rapid healing and future implant placement. The thin tip of the instrument is introduced into the periodontal space and slowly advanced toward the apex of the root while twisting the instrument gently back and forth. LM-LiftOut is not suitable for use as an elevator.



- Ergonomic design – comfortable to hold and easy to rotate
- Secure LM-ErgoTouch surface – light, non-slip grip
- Well-balanced – controlled grip
- Lightweight – less strain on the hands
- Clear color-coding – easy identification
- LM-DuraGradeMAX supersteel – strong and sharp blades
- Tight seam at blade juncture – hygienic

LM TwistOut

Elevators for forceful extraction



LM-TwistOut instruments are suitable for tooth extraction in situations where strong force or torque must be applied.



Easy to maintain

Cassette LM 810001

LM-LiftOut and LM-TwistOut instruments are supplied in a convenient cassette that protects both the instrument and the handler during the maintenance cycle. The cassette keeps the instrument from puncturing the sterilization pouch, so sterility is more likely to be assured. The cassette can be color-coded.

Sharpening stone LM 818002

The blades, hand-finished and made from LM-DuraGradeMAX supersteel, will stay sharp and last a long time when properly used and maintained. Sharpen the extracting instruments occasionally using this round sharpening stone as a simple and quick way to keep them in top condition and always ready to use.

Familiar LM ergonomics and tactile sensitivity, now in extracting instruments