VOCO**VOCO Ionofil Molar**

Instructions for use

0482

in accordance with EN 29917/ISO 9917/1994

Product description:

VOCO Ionofil Molar is a radiopaque glass polyalkenoate (ionomer) cement for fillings with chemical adhesion to enamel and dentine and continuously high release of fluoride ions.

VOCO Ionofil Molar is available in the shades A1, A2 and A3.

Fields of application:

- treatment of V-shaped defects and enamel or root erosions in the dental neck area incl. class V fillings
- restoration of deciduous teeth (especially class I)
- cavity fillings of class II (not occlusion-bearing) by means of tunnel preparation
- build-up fillings and linings, core build-up
- temporary long-term treatment of class I and II cavities

Application:

Preparation:

The adhesion of **VOCO Ionofil Molar** to enamel and dentine prevents the formation of marginal leakage and creates retentions. Undercuts give additional stability to the filling. The layer thickness of the fillings should not be less than 1 mm. Avoid feather-shaped margins. In case of doubt a conventional preparation should be chosen.

Clean the cavity (e.g. 3% hydrogen peroxide), thoroughly rinse with water and dry (warm air). Do not overdry dentine. In case of a remaining dentine thickness of more than 1.5 mm no additional lining is necessary.

To avoid pulp irritations which may occur in case of insufficient dentine thickness, a covering with calcium hydroxide of areas in proximity of the pulp is recommended.

Apply **VOCO Ionofil Molar** at room temperature (19°C - 23°C); in case of high temperatures the use of a cooled glass plate is advantageous. Low temperatures prolong the working time.

Furthermore, **VOCO Ionofil Molar** is suitable for use for the purposes of ART-technique (Atraumatic restorative treatment).

Mixing powder/liquid:


Shake powder well before use!

Mix one level measuring scoop of powder with one drop of liquid. This corresponds to a mixing ratio of 4.0 parts powder to 1 part liquid by weight. Hold dropping bottle in vertical position and place 1 drop of liquid onto a dry glass or plastic plate. Hold bottle high enough so that a single drop can fall down freely. Only use drops without enclosed air bubbles!

Mix the powder into the liquid in two equal portions using a clean metal spatula, or preferably with a clean **plastic spatula**. Only mix in the second portion when the first portion has been mixed in homogeneously. The mixed cement has a firm, packable consistency.

Avoid any contact with humidity!

Never add new liquid during mixing!

The mixing time takes approx. 50 - 60 s and should not exceed 1 min.

Application:

After application into the cavity **VOCO Ionofil Molar** will set within approx. 2 - 3 min. During placement of the filling the cavity has to be kept dry. Use of a rubber dam is recommended. Linings or build-up fillings can be carefully finished after setting (approx. 4 - 6 min after application).

During the setting phase, especially immediately after application, glass ionomer cements are sensitive to humidity (breath, saliva) and drying-out respectively. Therefore, if possible, cover immediately with a strip or matrix or apply a protective varnish.

Finishing and polishing (flexible polishing discs) of the fillings should follow 7 min after application at the earliest. Afterwards, cover again with varnish. The best finishing results will be achieved approx. 24 h after the filling has been inserted.

Additional information/precautions:

Close powder tightly after use in order to keep it dry and to maintain its high physical strength. Remove any remaining liquid from the outside of the liquid dispenser (by tissue) to keep opening clean for a correct mixing ratio.

VOCO Ionofil Molar should only be mixed immediately before use to ensure that the entire working time is available for placing and modelling. The development of fissures due to delayed use of the mixed cement can thus be prevented.

VOCO Ionofil Molar contains polyacrylic acid, fluorosilicate and parabenes. It should not be applied in case of known hypersensitivities or allergies.

Storage:

Store **VOCO Ionofil Molar** powder and liquid containers tightly closed at room temperature (15°C - 23°C) and protected from light. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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