



VOCO Profluorid® Varnish

DESENSITISING VARNISH WITH FLUORIDES

VOCO Profluorid® Varnish

HIGHLY EFFECTIVE

Hypersensitivity


Sensitive teeth are a widespread problem for a large number of people. Usually, the condition causes only transient pain which patients do not always report to their dentists. Hypersensitivity is mainly caused by exposed dentinal tubules. Therefore, the aim of treatment is to seal these tubules.

Indications

Profluorid Varnish is a fluoride-containing varnish for treating hypersensitivity. Dentinal tubules exposed on the necks of teeth or after cavity preparation are effectively sealed.


Components and effect of VOCO Profluorid® Varnish

Profluorid Varnish is a colophony-based varnish containing 5 % sodium fluoride (22,600 ppm fluoride). The fluoride ion, together with the calcium ions accumulated in the tubules, causes a precipitation of calcium fluoride, which effectively seals the tubules. According to a user survey, over 85 % of patients were completely free of pain after only one treatment with Profluorid Varnish.



As well as rapidly sealing the dentinal tubules, Profluorid Varnish also causes calcium fluoride deposits to form on the tooth surface. These deposits protect the tooth from acid attacks, promote remineralisation and contribute to the long-term formation of fluoroapatite.

Xylitol, which has a proven cariostatic effect in addition to its taste-enhancing properties, is also added to the varnish.



VOCO Profluorid® Varnish





SIMPLE AND AESTHETIC

Easy to apply

Thanks to its colophony matrix, Profluorid Varnish adheres well even to wet surfaces. It is not necessary to dry the tooth surfaces. This improves cooperation of patients with low levels of compliance (e.g. children).

Aesthetics and taste

Profluorid Varnish is a tooth-shaded varnish which does not impair aesthetics. Profluorid Varnish is available in four pleasant flavours: melon and cherry – aimed especially at children and younger patients, mint – a classic, particularly for adults, and caramel – to provide a slightly more exotic option. These four refined flavours will enable you to provide an even better service by tailoring it to the individual patient.



Survey of testing dentists

Opinion of dentists


Dentists using the product all confirmed its good handling and aesthetic properties. A vast majority also reported positive differences between Profluorid Varnish and competing products. Profluorid Varnish was equally well accepted by the patients. After the treatment, nearly all of the patients reported that they were either symptom-free or had a significantly reduced level of sensitivity.

Profluorid Varnish: much more than just good taste!

The *SingleDose* advantages

Profluorid Varnish *SingleDose* comes in a blister pack containing a brush and the varnish. This enables the product to be applied under perfectly hygienic conditions.

Simple and aesthetic desensitising with VOCO Profluorid® Varnish



VOCO Profluorid® Varnish

YOUR CHOICE

Indications

Treatment of hypersensitive teeth

Sealing of the dentinal tubules on sensitive root surfaces or following cavity preparation

Advantages

- Easy to apply on wet surfaces
- Aesthetic, white-transparent tooth colour
- Fast desensitisation and fluoride release (5 % NaF \pm 22,600 ppm fluoride)
- Available in four flavours: melon, mint, cherry and caramel

Advantages of the *SingleDose*

- Product and applicator in hygienic disposable packing
- No additional devices required
- No spilling in any position



Presentation

Flavor	SingleDose 48 x 0.40 ml	SingleDose 50 x 0.25 ml	SingleDose 50 x 0.40 ml	SingleDose 200 x 0.40 ml	Tube 10 ml
Melon	–	REF 1268	REF 1269	REF 1271	REF 1271
Caramel	–	–	REF 2226	REF 2229	REF 2229
Cherry	–	–	REF 2227	REF 2230	REF 2230
Mint	–	–	REF 2228	REF 2231	REF 2231
Mixed	REF 2225	–	–	–	–

VOCO GmbH
Anton-Flettner-Straße 1-3
27472 Cuxhaven
Germany

Tel.: +49 (0) 4721-719-0
Fax: +49 (0) 4721-719-140

info@voco.com
www.voco.com

Available from: