

# Certification Pack.

# N95 BYD Face Mask. This pack contains TGA Certificate **EN Test Report ASTM 2100 Test Report** Date Released November 2020 Certification Pack 73700 | Veridia Australia



#### **Australian Government**

#### **Department of Health**

#### Therapeutic Goods Administration

#### **Public Summary**

Summary for ARTG Entry: 337767 BYD AUSTRALIA PTY LTD - Surgical/medical respirator, single-use

ARTG entry for Medical Device Included Class 1

Sponsor BYD AUSTRALIA PTY LTD

Postal Address Unit 11 9-11 Chaplin Drive, Lane Cove West, NSW, 2066

Australia

ARTG Start Date 12/06/2020

Product Category Medical Device Class 1

Status Active

Approval Area Medical Devices

#### Conditions

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

#### Manufacturers

Name Address

BYD Precision Manufacture Co Ltd

No 1 Baoping Road

Long gang District, Shenzhen,

China

#### **Products**

#### 1 . Surgical/medical respirator, single-use

Product Type Single Device Product Effective Date 12/06/2020

**GMDN** 57794 Surgical/medical respirator, single-use

Intended Purpose The N95 Healthcare particulate respirator is designed to cover the nose and mouth of a healthcare worker to permit normal

breathing while protecting the wearer and patient from particles during medical, surgical, dental, and isolation procedures. It is a

single-use device.

#### **Specific Conditions**

No Specific Conditions included on Record

© Commonwealth of Australia. This work is copyright. You are not permitted to re-transmit, distribute or commercialise the material without obtaining prior written approval from the Commonwealth. Further details can be found at http://www.tga.gov.au/about/website-copyright.htm.





中国认可 国际互认 检测 TESTING CNAS L1499

# National Quality Supervision and Testing Center for Personal Protective Equipment (Beijing) (Testing Laboratory for Labour Protection Products of Politics

# (Testing Laboratory for Labour Protection Products of Beijing Municipal Institute for Labour Protection)

No.55 Taoranting Street, Xicheng District, Beijing, China. Phone: +86 10 63519250 +86 10 63520770 +86 10 83530311

Fax: +86 10 63519250 +86 10 63520770

The Testing Center is accredited for compliance with ISO/IEC 17025.

The results of tests, calibrations and/or measurements included in this document are traceable to Chinese/national standards.

CNAS is a signatory to the ILAC mutual recognition arrangement for the mutual recognition of the equivalence of testing, calibration and inspection reports.

#### **TEST REPORT**

#### Particulate respirator-half facepiece

EN 149: 2001 +A1: 2009 Respiratory protective devices — Filtering half masks to protect against particles — Requirements, testing, marking

**Product:** Particle filtering half mask

**Report No:** 2020 (D) - 0615

Client: CCQS Certification Services Limited

**Model (s): DE2323** 

Date(s) of tests: 2020.05.03-2020.05.16

#### DESCRIPTION OF SAMPLES

**General Information**Classification
FFP2 NR

Main Components
White folding mask

Manufacturer BYD Precision Manufacture Co., Ltd.

Manufacturer Address

No.3001 Baohe Road, Baolong Industrial City, Longgang Street, Longgang District,

Shenzhen, P. R China

Signed:

陈强为

Issued: 2020.5.16

陈倬为 Chen Zhuowei Authorized Signatory, Lab Director

Page 1 of 10

Report No: 2020 (D) - 0615 Page 2 of 10

## **Conditions:**

The test results presented in this report relate to the samples tested only.

This report may be reproduced and distributed to your clients, provided that it is reproduced and distributed in full.

The authenticity of this test report and its contents can be verified by contacting the laboratory.

Report No: 2020 (D) - 0615

#### **Test Results**

7.3 Visual inspection Not tested<sup>1</sup>

The visual inspection shall include the marking and information supplied by the manufacturer.

Note1: As requested by the client, marking and information supplied by the manufacturer was not inspected.

7.4 Package Pass<sup>2</sup>

Particle filtering half masks shall be offered for sale packaged in such a way that they are protected against mechanical damage and contamination before use.

Note2: In accordance with the requirement.

7.5 Material Pass<sup>3</sup>

Materials used shall be suitable to withstand handling and wear over the period for which the particle filtering half mask is designed to be used.

Any material from the filter media released by the air flow through the filter shall not constitute a hazard or nuisance for the wearer.

After undergoing the conditioning described in 8.3.1 none of the particle filtering half masks shall have suffered mechanical failure of the facepiece or straps.

When conditioned in accordance with 8.3.1 and 8.3.2 the particle filtering half mask shall not collapse.

Note3: No mechanical failure after undergoing the conditioning described in 8.3.1. No collapse when conditioned in accordance with 8.3.1 and 8.3.2.

#### 7.6 Cleaning and disinfecting

 $N/A^4$ 

If the particle filtering half mask is designed to be re-usable, the materials used shall withstand the cleaning and disinfecting agents and procedures to be specified by the manufacturer.

Note4: Single shift use only.

#### 7.7 Practical performance

Pass<sup>5</sup>

The particle filtering half mask shall undergo practical performance tests under realistic conditions.

Note5: No imperfections.

7.8 Finish of parts Pass<sup>6</sup>

Parts of the device likely to come into contact with the wearer shall have no sharp edges or burrs. Note6: No sharp edges or burrs.

# 7.9.1 Total inward leakage

Pass<sup>7</sup>

For particle filtering half masks fitted in accordance with the manufacturer's information, at least 46 out of the 50 individual exercise results (i.e. 10 subjects x 5 exercises) for total inward leakage shall be not greater than: 25% for FFP1, 11% for FFP2, 5% for FFP3

and, in addition, at least 8 out of the 10 individual wearer arithmetic means for the total inward leakage shall be not greater than

22% for FFP1, 8% for FFP2, 2% for FFP3

Note7: FFP2 respirator. Test results are shown in Annex A Table 7.9.1-A&B.

#### 7.9.2 Penetration of filter material

Pass8

The penetration of the filter of the particle filtering half mask shall meet the requirements of Table 1.

Sodium chloride test 95 l/min

Paraffin oil test 95 l/min

FFP1 ≤20% ≤20%

Report No: 2020 (D) - 0615

FFP2  $\leq 6\%$   $\leq 6\%$  FFP3  $\leq 1\%$   $\leq 1\%$ 

Note8: FFP2 respirator. Test results are shown in Annex A Table 7.9.2.

#### 7.10 Compatibility with skin

Pass9

Materials that may come into contact with the wearer's skin shall not be known to be likely to cause irritation or any other adverse effect to health.

Note9: No irritation or any other adverse effect to health.

7.11 Flammability Pass<sup>10</sup>

When tested, the particle filtering half mask shall not burn or not to continue to burn for more than 5 s after removal from the flame.

Note10: Test results are shown in Annex A Table 7.11.

#### 7.12 Carbon dioxide content of the inhalation air

Pass<sup>11</sup>

The carbon dioxide content of the inhalation air (dead space) shall not exceed an average of 1,0 % (by volume) Note11: Test results are shown in Annex A Table 7.12.

7.13 Head harness Pass<sup>12</sup>

The head harness shall be designed so that the particle filtering half mask can be donned and removed easily. The head harness shall be adjustable or self-adjusting and shall be sufficiently robust to hold the particle filtering half mask firmly in position and be capable of maintaining total inward leakage requirements for the device.

Note12: Head harness can be donned and removed easily, adjustable or self-adjusting and have sufficiently robust to hold the particle filtering half mask firmly.

7.14 Field of vision Pass<sup>13</sup>

The field of vision is acceptable if determined so in practical performance tests.

Note13: Pass the practical performance tests.

7.15 Exhalation valve N/A<sup>14</sup>

A particle filtering half mask may have one or more exhalation valve(s), which shall function correctly in all orientations.

If an exhalation valve is provided it shall be protected against or be resistant to dirt and mechanical damage and may be shrouded or may include any other device that may be necessary for the particle filtering half mask to comply with 7.9.

Exhalation valve(s), if fitted, shall continue to operate correctly after a continuous exhalation flow of 300 l/min over a period of 30 s.

When the exhalation valve housing is attached to the faceblank, it shall withstand axially a tensile force of 10 N applied for 10 s.

Note14: No exhalation valve.

#### 7.16 Breathing resistance

Pass<sup>15</sup>

Classification	Maximum permitted resistance (mbar)						
	Inhalation	Exhalation					
	30 l/min	95 l/min	160 l/min				
FFP1	0.6	2.1	3.0				
FFP2	0.7	2.4	3.0				
FFP3	1.0	3.0	3.0				

Note15: FFP2 respirator. Test results are shown in Annex A Table 7.16.

Report No: 2020 (D) - 0615 Page 5 of 10

7.17 Clogging N/A<sup>16</sup>

#### 7.17.2 Breathing resistance

Valved particle filtering half masks:

After clogging the inhalation resistances shall not exceed:

FFP1: 4 mbar, FFP2: 5 mbar, FFP3: 7 mbar at 95L/min continuous flow

The exhalation resistance shall not exceed 3 mbar at 160 L/min continuous flow

Valveless particle filtering half masks

After clogging the inhalation and exhalation resistances shall not exceed:

FFP1: 3 mbar, FFP2: 4 mbar, FFP3: 5 mbar at 95L/min continuous flow

#### 7.17.3 Penetration of filter material

So	dium chloride test 95 l/min	Paraffin oil test 95 l/min
FFP1	€20%	<b>≤</b> 20%
FFP2	≤6%	≤6%
FFP3	≤1%	≤1%
Note16: Single	e shift use only.	

#### 7.18 Demountable parts

 $N/A^{17}$ 

All demountable parts (if fitted) shall be readily connected and secured, where possible by hand Note17: No demountable parts.

9 Marking Not tested

## 9.1 Packaging

The following information shall be clearly and durably marked on the smallest commercially available packaging or legible through it if the packaging is transparent.

- 9.1.1 The name, trademark or other means of identification of the manufacturer or supplier.
- **9.1.2** Type-identifying marking.
- **9.1.3** Classification

The appropriate class (FFP1, FFP2 or FFP3) followed by a single space and then: "NR" if the particle filtering half mask is limited to single shift use only. Example: FFP3 NR, or "R" if the particle filtering half mask is re-usable. Example: FFP2 R D.

- **9.1.4** The number and year of publication of this European Standard.
- **9.1.5** At least the year of end of shelf life. The end of shelf life may be informed by a pictogram as shown in Figure 12a, where yyyy/mm indicates the year and month.
- **9.1.6** The sentence 'see information supplied by the manufacturer', at least in the official language(s) of the country of destination, or by using the pictogram as shown in Figure 12b.
- **9.1.7** The manufacturer's recommended conditions of storage (at least the temperature and humidity) or equivalent pictogram, as shown in Figures 12c and 12d.
- **9.1.8** The packaging of those particle filtering half masks passing the dolomite clogging test shall be additionally marked with the letter "D". This letter shall follow the classification marking preceded by a single space.

#### 9.2 Particle filtering half mask

Particle filtering half masks complying with this European Standard shall be clearly and durably marked with the following:

9.2.1 The name, trademark or other means of identification of the manufacturer or supplier.

Report No: 2020 (D) - 0615

- 9.2.2 Type-identifying marking.
- **9.2.3** The number and year of publication of this European Standard.
- 9.2.4 Classification

The appropriate class (FFP1, FFP2 or FFP3) followed by a single space and then: "NR" if the particle filtering half mask is limited to single shift use only. Example: FFP3 NR, or "R" if the particle filtering half mask is re-usable. Example: FFP2 R D.

- **9.2.5** If appropriate the letter D (dolomite) in accordance with clogging performance. This letter shall follow the classification marking preceded by a single space
- **9.2.6** Sub-assemblies and components with considerable bearing on safety shall be marked so that they can be identified.



Report No: 2020 (D) - 0615 Page 7 of 10

# **Annex A: Summarization of Test Data**

# Table 7.9.1-A Inward leakage test data

Test specification: EN 149-2001 Clause 8.5

Subject	Sample No.	Condition	Walk(%)	Head Side/side(%)	Head up/down(%)	Talk(%)	Walk(%)	Mean(%)
Yi	1	A.R.	9.13	9.29	9.61	9.41	9.22	9.3
Gong	2	A.R.	6.33	6.61	6.61	6.61	6.43	6.5
Yu	3	A.R.	7.10	7.15	7.12	7.40	7.36	7.2
Hu	4	A.R.	7.32	7.47	7.64	7.71	7.51	7.5
Xu	5	A.R.	7.89	8.48	7.91	7.91	8.33	8.1
Deng	6	T.C.	6.32	6.75	6.78	6.45	6.46	6.6
Zhang	7	T.C.	6.62	6.70	6.70	7.05	6.86	6.8
Liu	8	T.C.	6.11	6.58	6.13	6.14	6.13	6.2
Fang	9	T.C.	7.11	7.12	7.43	7.43	7.22	7.3
Lv	10	T.C.	7.20	7.72	7.61	7.60	7.50	7.5
	All <u>50</u> individual exercise results were not greater than <u>11</u> % 8 out of <u>10</u> individual wearer arithmetic means were not greater than $\leq$ <u>8</u> %							Pass

#### **Table 7.9.1-B Facial dimension**

Table 7.3.1-D Factal difficusion								
Subject	Face length	Face Width	Face Depth	Mouth Width				
Yi	120	130	109	59				
Gong	122	140	115	65				
Yu	119	160	139	55				
Hu	112	122	119	63				
Xu	110	130	118	60				
Deng	115	119	110	59				
Zhang	112	123	113	55				
Liu	103	130	100	50				
Zhi	118	139	130	63				
Fang	115	129	120	50				
Chen	116	150	132	56				
Lv	110	121	110	53				

Report No: 2020 (D) - 0615 Page 8 of 10

# **Table -7.9.2 Penetration of filter material**

Test specification: EN 149-2001 Clause 8.11

Aerosol	Condition	Condition Sample Penetration No. (%)				
		11	0.245			
	As received	12	0.226			
		13	0.385			
		14	0.431			
Sodium chloride test	Simulated wearing treatment	15	0.444			
		16	0.517			
		17	0.498			
	Mechanical strength+ Temperature conditioned	18	0.572			
		19	0.616			
		20	2.24	Pass		
	As received	21	2.15			
		22	2.16			
		23	2.51			
Paraffin oil test	Simulated wearing treatment	24	2.54			
		25	2.47			
		26	2.81			
	Mechanical strength+ Temperature conditioned	27	2.84			
		28	2.73			
Flow condition	Flow conditioning: Single filter: 95.0 L/min					

# **Table 7.11 Flammability**

Test specification: EN 149-2001 Clause 8.6

Condition	Sample No.	Result	Assessment
A 1	29	Burn for 1 s	
As received	30	Burn for 1 s	D
Temperature	31	Burn for 2 s	Pass
conditioned	32	Burn for 1 s	

Report No: 2020 (D) - 0615 Page 9 of 10

Table 7.12 Carbon dioxide content of the inhalation air

Test specification: EN 149-2001 Clause 8.7

Condition	Sample No.	Result	Assessment	
	33	0.41%		
As received	34	0.42%	Mean value 0.4%	Pass
	35	0.41%		

# **Table 7.16 Breathing resistance (mbar)**

Test specification: EN 149-2001 Clause 8.9


•	F1	4			36					37					38		
	Flow	rate	Α	В	С	D	Е	A	В	С	D	Е	Α	В	С	D	Е
As received	Inhalation	30 l/min	0.4	0.5	0.5	0.6	0.6	0.5	0.5	0.5	0.5	0.5	0.4	0.4	0.6	0.5	0.6
	Illialation	95 l/min	1.5	1.6	1.5	1.7	1.7	1.7	1.6	1.7	1.6	1.6	1.6	1.6	1.6	1.5	1.5
	Exhalation	160 l/min	1.6	1.7	1.8	1.7	1.8	1.7	1.7	1.8	1.7	1.8	1.6	1.6	1.7	1.8	1.6
	Flow	rato			39					40					41		
Simulated	FIOW	Tate	Α	В	С	D	Е	Α	В	С	D	Е	Α	В	С	D	Е
wearing	Inhalation	30 l/min	0.6	0.5	0.6	0.6	0.5	0.5	0.4	0.4	0.5	0.6	0.4	0.5	0.6	0.4	0.6
treatment		95 l/min	1.6	1.6	1.5	1.6	1.6	1.6	1.6	1.6	1.6	1.5	1.5	1.7	1.6	1.6	1.7
	Exhalation	160 l/min	1.7	1.8	1.8	1.7	1.8	1.6	1.7	1.8	1.8	1.7	1.8	1.7	1.7	1.7	1.7
	Flow	rato			42					43					44		
Tomas onotions	FIOW	Tate	Α	В	С	D	Е	Α	В	С	D	Е	Α	В	С	D	Е
Temperature conditioned	Inhalation	30 l/min	0.5	0.6	0.4	0.6	0.5	0.5	0.5	0.6	0.5	0.4	0.5	0.6	0.4	0.4	0.4
Conditioned		95 l/min	1.7	1.6	1.7	1.6	1.6	1.5	1.7	1.5	1.7	1.7	1.6	1.6	1.7	1.6	1.7
	Exhalation	160 l/min	1.8	1.7	1.6	1.7	1.7	1.6	1.7	1.6	1.7	1.7	1.7	1.8	1.7	1.7	1.8
Assessment	Pass																

A: facing directly ahead; B: facing vertically upwards; C: facing vertically downwards; D: lying on the left side; E: lying on the right side



Report No: 2020 (D) - 0615 Page 10 of 10

# ANNEX B PHOTOS OF SAMPLES





**End of Annex B** 



# **TEST REPORT**

Prepared For :	BYD Precision Manufacture Co., Ltd.
	No.3001 Baohe Road, Baolong Industrial City, Longgang Street,Longgang District,Shenzhen, Guangdong Province, China
Trade Mark :	N/A
Product Name :	BYD N95 Particulate Respirator Mask
Model(s):	N/A
Number:	One hundred
Prepared By:	Shenzhen CCT Testing Technology Co., Ltd. TING TECHNOR Sth Floor, Area I, Building 1, Hanhaida Science and Technology Innovation Park, Guangming New District, Shenzhen, Guangdong, China Tel: +86 755-3315 7675 Fax: +86 755-33609691 Email: cct@cct-prc.com Web: www.cct-prc.com Testing Service
Test Date:	Mar. 24, 2020 - Apr. 2, 2020
Date of Report:	Apr. 2, 2020
Report No. :	CCT20032405TRS

**Note:** This test report is limited to the above client company and the product model only. It may not be duplicated without prior written consent of Shenzhen CCT Testing Technology Co., Ltd.



# **TEST REPORT ASTM F-2100**

# **Standard Specification for**

#### Performance of Materials Used in Medical Face Masks 1

Date of issue...... Apr. 2, 2020

Contents ...... 6 pages

Testing laboratory

Name ...... Shenzhen CCT Testing Technology Co., Ltd.

8th Floor, Area I, Building 1, Hanhaida Science and Technology

Address ...... Innovation Park, Guangming New District, Shenzhen,

Guangdong, China

Testing location ...... Same as above

Client

Name ...... BYD Precision Manufacture Co., Ltd.

No.3001 Baohe Road, Baolong Industrial City, Longgang

Address ...... Street, Longgang District, Shenzhen, Guangdong Province,

China

Test specification

Procedure deviation .....: N.A.

Non-standard test method .....: N.A.



#### SHENZHEN CCT TESTING TECHNOLOGY CO., LTD. NO.: CCT20032405TRS

#### **Test case verdicts**

Test case does not apply to the test object .....: N(.A.)

Test item does meet the requirement...... P(ass)

Test item does not meet the requirement...... F(ail)

Testing

Date of receipt of test item ...... Mar. 24, 2020

Date(s) of performance of test ...... From Mar. 24, 2020 to Apr. 2, 2020

#### General remarks

This test report shall not be reproduced except in full without the written approval of the testing laboratory.

The test results presented in this report relate only to the item tested.

"(see remark #)" refers to a remark appended to the report.

"(see appended table)" refers to a table appended to the report.

Throughout this report a comma is used as the decimal separator.

#### Remark:

The product complies with requirements of the standard.



## SHENZHEN CCT TESTING TECHNOLOGY CO., LTD. NO.: CCT20032405TRS

Serial number	Test item	Standard Requirements	Standard	Single item conclusion
1	Hight Fluid Resistance	160mm Hg		Р
2	Bacterial filtration efficiency (BFE)	The bacterial filtration efficiency of mask shall be ≥98%	ASTM F-2100	P
3	Particulate filtration efficiency(PFE)	The Particulate filtration efficiency of mask shall be ≥98%	Level 3	P
4	Flame spread	Class 1		Р
5	Filtration efficiency	The Filtration efficiency≥95%	NIOSH 42 CED DADT 94	P
6	Exhalation Resistance Test	≤35mm H <sub>2</sub> O	NIOSH 42 CFR PART 84	Р
7	Inhalation Resistance Test	≤25mm H <sub>2</sub> O		Р

# **Test Results**

Bacterial Filtration Efficiency	Р		
Models Item	Sample 1	Sample 2	Sample 3
BFE (%)	99.96	99.98	99.96

Sub-Micron Particulate Filtrati	Р		
Models Item	Sample 1	Sample 2	Sample 3
Sub-Micron Particulate Filtration at 0.1 micron (%)	99.96	99.94	99.94

Filtration efficiency			Р
Models Item	Sample 1	Sample 2	Sample 3
Filtration efficiency (%)	99.72	99.76	99.68



## SHENZHEN CCT TESTING TECHNOLOGY CO., LTD. NO.: CCT20032405TRS

Exhalation Resistance Test			Р
Models Item	Sample 1	Sample 2	Sample 3
Exhalation Resistance (mmH <sub>2</sub> O)	16.5	16.2	17.0

Inhalation Resistance Test			Р
Models Item	Sample 1	Sample 2	Sample 3
Inhalation Resistance (mmH <sub>2</sub> O)	15.7	15.7	15.7



# **Attachments: real photos**



Photo 1



Photo 2