

Ambu[®] BlueSensor electrodes

Electrodes	Reference	Size (mm)	Features	Connectors	Recommended applications	Page
	SU	49 x 33	Offset connector Wet gel Polymer backing	A F C	Resting ECG	Page 7
	L VL* VLC	68 x 55 72 x 68*** 72 x 68	Offset connector Wet gel Double ring adhesive Microporous backing	S A (VLC: S only)	Holter Long-term monitoring Event Recording	Page 9
	R T SP	57 x 48 57 x 35 42 x 38	Offset connector Wet gel Foam backing	A S	Stress Test / EMS Short-term monitoring	Page 11
	Q** M P	53 x 40*** 40 x 34 40 x 34	Offset connector Wet gel Foam/polymer backing	A F S (Q: A, S)	Short -term monitoring	Page 13 Page 13 Page 17
	QR	40 x 40	Radiolucent Wet gel Superior adhesion Pre-attached leadwire	A	Short-term monitoring	Page 13
	NEOX	25 x 20	Tiny size Four-leaf clover shape Transparent material No-touch application system Long-term adhesion	K	Neonatal	Page 14
	NF	44 x 22***	Solid gel Offset connector or pre-attached leadwire Soft cloth backing	A F K S J	Neonatal	Page 15
	BR BRS	28 x 20 16 x 19	Solid gel Soft cloth backing Radiolucent Special shape for improved adhesion	A K F: BR only	Neonatal	Page 15
	N	44 x 22***	Wet gel Foam backing Offset connector or pre-attached leadwire	A F K S	Paediatric	Page 17

*VL is available with leadwires and A, K connectors.

**Q is available with leadwires and A, F, K connectors.

***Electrode size stated with snap connector. Electrode size is smaller for versions with leadwires. See details on product page.

Connectors and Fittings:



A = 4 mm Connector

F = 3 mm connector

J = 2 mm connector

K = 1.5 mm connector

A = 4 mm fitting

C = Tab

F = 3 mm fitting

S = Stud