

ONE-LINK Needle-free IV Connector

IV Access.
Connect to the power of choice.

Enhancing patient care. Simply. Neutrally.

- **Power injector compatible** (up to 325 psi)*
- **Straight fluid path**
- **Non-PVC and Non-DEHP**
- **Easy to clean surface**
- **Low priming volume of .08 mL**

The neutral displacement ONE-LINK connector is designed to help reduce the risk of thrombotic catheter occlusions compared to devices with higher reflux volumes.¹ Less chance of reflux means simplified training and usage – and more time to focus on your patients.

Plus, ONE-LINK is a needle-free connector that can be used throughout your facility. From OR to ER. Neonatal to Oncology.



One·Link
Needle-free IV Connector

Preventing Infections and improving Medication Safety is our business.
Let's raise the standards and help make healthcare...

Safer together

*Prior to power injection, ensure all parts of the IV System are compatible. Replace if a pressure over 325 psi is applied to the connector. The ONE-LINK IV Connector is intended for single patient use with a vascular access device for administration of drugs and solutions without needles and can be used to aspirate blood. For safe and proper use of this device please refer to the Instructions for Use.

¹ Data on file, Baxter Healthcare Corporation.
Baxter and One-Link are trademarks of Baxter International Inc.

Baxter

Baxter Healthcare Pty Ltd
1 Baxter Drive, Old Toongabbie
NSW 2146 Australia
Tel: + 61 1800 BAXTER (+61 1800 229 837)
Fax: +61 2 98481123

Baxter Healthcare Ltd
33 Vestey Drive, Mt Wellington
Auckland 1060, New Zealand
Tel: + 64 (0)800 BAXTER (+64 (0)800 229 837)
Fax: +64 (0)800 229 239

ANZ-MD35-250001