

bodigrip*

Elasticated Tubular Support Bandage, 1M

Specially designed to provide firm effective support for sprains, strains and weak joints

- Easy to apply and re-apply.
- No need for pins or tapes.
- Comfortable – provides even pressure and allows freedom of movement.
- Stays in place – close weave, seamless knitted bandage provides continuous support.
- Available in full range of sizes, fully washable.

Code Flesh	SIZE	Application	Bandage Width
13023161	B	Small wrist, elbow	6.5cm
13023162	C	Medium wrist, elbow /Small ankle	6.75cm
13023163	D	Large wrist, elbow /Medium ankle /Small knee	7.5cm
13023164	E	Large ankle /Medium knee /Small thigh	8.75cm
13023165	F	Large knee /Medium thigh	10cm
13023166	G	Large thigh	12cm

Comprehensive, easy to use directions and size guide on all retail packs



1. Cut to twice the required length for the limb. Make two small cuts for the thumb.



2. Pull bandage over wrist and push thumb through hole.



3. Double bandage back over wrist and push thumb through second hole. Pull top layer about 3cm past bottom layer.

