





Ascorbic Acid 2,000mg/L

DATA SHEET

Instructions for use

- 1. Immerse an indicator strip in the test solution for 1 second.
- 2. Shake off excess solution.
- 3. Wait 45 seconds.
- 4. Compare the colour produced on the strip against the colour chart.

Measuring range

0 - 50 - 100 - 200 - 300 - 500 - 700 - 1,000 - 2,000mg/L

Application

Ascorbic acid test strips detect the level of ascorbic acid (vitamin C) in fruits and fruit juices. These strips can also be used to measure the levels of ascorbic acid acting as preservatives or antioxidants in food stuffs.

Shelf Life

3 years from manufacturing date.

Contents

100 strips

Dimensions

Strip dimension: 83mm x 5mm

Notes

Keep any unused strips in the original container. Do not store where they are liable to be exposed to acid or alkaline gases or vapours. Store in a cool dry place.

E.&O.E.