



CODE: STR-162-50

Peracetic Acid 50mg/L

DATA SHEET

Instructions for use

1. Immerse an indicator strip in the test solution for 3 seconds.
2. Shake off excess solution.
3. Wait 45 seconds.
4. Compare the colour produced on the strip against the colour chart.

Measuring range

0 - 5 - 10 - 20 - 30 - 50mg/L

Application

Peracetic acid test strips detect residual peracetic acid, for example, in the food, cosmetic, beverage, pharmaceutical and medical industries, where peracetic acid may be used as a disinfectant.

Shelf Life

3 years from manufacturing date.

Contents

100 strips

Dimensions

Strip dimension: 83mm x 5mm

Notes

Keep any unused strips in the original container. Do not store where they are liable to be exposed to acid or alkaline gases or vapours. Store in a cool dry place.

