

2019 Product Catalog US & Canada





TABLE OF CONTENTS

ACTEON IMAGING

Intraoral Radiology	4
Extraoral Radiology	14
Intraoral Cameras	

ACTEON EQUIPMENT

Ultrasonics	26
Air Polisher	45
Surgery	50
Curing Lights	

ACTEON PHARMA

Gingival Retraction	69
Hemostasis & Tissue Regeneration	70
Hygiene	71

ACTEON INTEGRATED

Ultrasonics	72
	77
	79

INTRAORAL RADIOLOGY Wall-mounted units



X-MIND™ UNITY

X-MIND™ Unity is pre-wired for the integration of the SOPIX² digital sensor. Delivered with a long cone in standard version.

X-MIND™ UNITY - MOBILE STAND W0800140

X-MIND Unity Mobile Stand (XMU NOT INCLUDED)

X-MIND™ UNITY - TOP MOUNTING

XMU TOP 40 SHORT CONE

HXMU-TS

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 16"

XMU TOP 80 SHORT CONE

HXMU-TM

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 31"

XMU TOP 110 SHORT CONE

HXMU-TL

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 43"



INTRAORAL RADIOLOGY Wall-mounted units

X-MIND™ UNITY - BOTTOM MOUNTING

XMU BOTT 40 SHORT CONE

HXMU-BS

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 16"

XMU BOTT 80 SHORT CONE

HXMU-BM

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 31"

XMU BOTT 110 SHORT CONE

HXMU-BL

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 43"

Digital X-ray Systems to integrate into X-MIND™ Unity

CMOS sensor with optic fiber of 1.5 million of pixels (size 1) or 2.21 millions of pixels (size 2), theorical resolution of 25 pl/mm, actual resolution > 12 pl/mm, 3.70m cable. Includes ACE technology (Automatic Control Exposure).

Delivered with a sensor holder, 10 sensor sheaths, and SOPRO Imaging Software.

SOPIX® INSIDE SENSOR SIZE 1

S 802 6500

SOPIX (SIZE 1) DIGITAL SENSOR - DIRECT USB2

SOPIX® INSIDE SENSOR SIZE 2 S 802 6501

SOPIX (SIZE 2) DIGITAL SENSOR - DIRECT USB2





INTRAORAL RADIOLOGY Accessories

X-MIND UNITY ACCESSORIES



X-MIND™ Unity short cone W0800018 SHORT CONE: 8" Ø60 (20cm)



X-MIND™ Unity long cone W0800019 LONG CONE: 12" Ø60 (31cm)



X-MIND™ Unity short square cone W0800020 Square SHORT CONE: 8″ Ø60 (20cm)





X-MIND™ Unity remote exposure switch W0800021
REMOTE EXPOSURE SWITCH



WALLPLATE for X-MIND™ Unity installation SSWALLPLATE
Single Stud Wallplate



XMU EXTERNAL LIGHT 115V W0800034

INTRAORAL RADIOLOGY Digital radiology system by phosphor plate



PSPIX^{®2} S 702 0011


Digital imaging plate system

Includes 4 imaging plates (2 x size 1, 2 x size 2), 200 bag & cover for imaging plate (100 x size 1, 100 x size 2), an imaging plate storage box, 1 stylus for touchscreen,1 power supply 100-240 VAC, 1 Ethernet cable, SOPRO Imaging software and its license, and a Mac and its license.



STANDARD IMAGING PLATE





	PSPIX ² Set of 2	PSPIX 1 st generation Set of 2
Size 0 (22 x 35mm)	990215	700178
Size 1 (24 x 40mm)	990216	700179
Size 2 (31 x 41mm)	990217	700180
Size 3 (27 x 54mm)	990218	700181

INTRAORAL RADIOLOGY Accessories

IDOT™IMAGING PLATE

Numbered plate for better identification.



	PSPIX ²	PSPIX 1 st generation set of 6
Size 0 (22 x 35mm)		700182
Size 1 (24 x 40mm)		700183
Size 2 (31 x 41mm)		700184
Size 3 (27 x 54mm)		700185

PROTECTIVE BAG FOR IMAGING PLATE



BAG & COVER FOR IMAGING PLATE with integrated cardboard envelope









PSPIX²



Size 1

Size 2



PSPIX 1st generation

		4
Size 0	Size 1	
	PSPIX ² set of 300	

	PSPIX ² set of 300	PSPIX 1st generation set of 250
Size 0	700546	700338L
Size 1	700547	700339L
Size 2	700548	700340L
Size 3	700549	700341L

PSPIX² ACCESSORIES



Bag & cover for imaging plate start up kit

(set of 100 pieces of size 1 and 100 pieces of size 2)



Imaging plate storage box 990219

(delivered without imaging plates)



Stylus for PSPIX² touchscreen 005132



Cleaning plate for PSPIX² 990243

INTRAORAL RADIOLOGY Accessories



Autoclavable Insertion Guide 608019

Autoclavable insertion guide (Top) for imaging plates

Autoclavable Receptacle 608018

Autoclavable receptacle (Bottom) for imaging plates



Size 4 kit for PSPIX² 990250

1 plastic SUPPORT size 4, 100 protective bags size 4, 3 imaging plates size 3, 1 quick set-up, 1 quick start

Plastic Support Size 4 990251

Plastic case for 3 size 3 plates

BAGS&COVERS S4 - 100 700644

Bag and cover for imaging plates size 4 (set of 100)



PSPIX Positioning Set - Bitewing 700195

Includes five positoning holders with ring for bitewing and one centring devices.



PSPIX Positioning Set - Anterior & Posterior 700196

Includes two positoning holders with ring for anterior teeth, two positoning holders with ring for posterior teeth and two centering devices.



SOPROTABS 700406

Adhesive foam tab (500 pieces in a single pasteboard)



Imaging Plate Storage Box - Gen 1 700194

(delivered without imaging plates)

INTRAORAL RADIOLOGY Digital radiology sensors





CMOS sensor with optic fiber of 1.5 million of pixels (size 1) or 2.21 millions of pixels (size 2), theorical resolution of 25 pl/mm, actual resolution > 18 pl/mm, 3.70m cable. Includes ACE technology (Automatic Control Exposure).

Delivered with a sensor holder, 10 sensor sheaths, the ACTEON® Imaging Suite (AIS) software.

SOPIX®2 (SIZE 2) DIGITAL SENSOR

S 802 0005

SOPIX² Direct USB2 - comes w/ one sensor holder

SOPIX®2 (SIZE 1) DIGITAL SENSOR S 802 0006

SOPIX2 (SIZE 1) DIGITAL SENSOR - DIRECT USB2 - COMES W/ ONE SENSOR HOLDER

INTRAORAL RADIOLOGY Accessories

ACCESSORIES



Positioning set for SOPIX²/SOPIX series sensors 700042

SOPIX² Positioning Set (5 sets with positioning Ring for: BW, Ant, Post, Endo Ant, Endo Post)



Sopix²/ Sopix sensor holder 926039

Size 1 or Size 2



Sheaths for Sopix²/Sopix

(500 pieces in a single pasteboard) for size 1 SOPIX² series sensors for size 2 SOPIX² series sensors

403026 403027



Sensor Holder Combo Kit ACT-459-1+2

Size 1&2 Combo kit - includes 1 Ring & rod system, 4 size 1 bite blocks, 4 size 2 bite blocks (1 each size Anterior, Posterior, Vertical and Horizontal Bitewings)

Sensor Holder Size 2 Kit ACT-459-2

Size 2 kit - includes 1 Ring & rod system, 8 size 2 bite blocks (2 each: Anterior, Posterior, Vertical and Horizontal Bitewings)

Sensor Holder Size 1 Kit ACT-459-1

Size 1 Kit (includes 1 Ring & rod system, 8 size 1 bite blocks: 2 each: Anterior, Posterior, Vertical and Horizontal Bitewings)

INTRAORAL RADIOLOGY Accessories

Sensor Holder Bite Blocks Size 2 *ACT-462*

16 pack of bite blocks size 2 (4 each: Anterior, Posterior, Vertical and Horizontal Bitewings)

Sensor Holder Bite Blocks Size 1 ACT-461

16 pack of bite blocks size 1 (4 each: Anterior, Posterior, Vertical and Horizontal Bitewings)

*Sensors, ring and rod not included



3 pack of bite blocks size 1 Anterior *ACT-001*

3 pack of bite blocks size 2 Anterior *ACT-049*

*Sensors, ring and rod not included



3 pack of bite blocks size 1 Posterior *ACT-018*

3 pack of bite blocks size 2 Posterior *ACT-056*

*Sensors, ring and rod not included



3 pack of bite blocks size 1 Horizontal Bitewing *ACT-025*

3 pack of bite blocks size 2 Horizontal Bitewing ACT-063

3 pack of bite blocks size 1 Vertical Bitewing *ACT-032*

3 pack of bite blocks size 2 Vertical Bitewing ACT-070

*Sensors, ring and rod not included

3D imaging system

- Easily upgradeable
- Superior 3D image processing with SHARP and STAR
- Patented Ceph collimation
- Large selection of F.O.V: Ø40x40, 60x60, 80x80 and 110x80

X-MIND™ TRIUM PAN

Control panel on the right Control panel on the left

W0900001 W0900007



OrthoPantomograph (OPG), 115-220 V, With One CMOS Pan 2D sensor for panoramic imaging

- Adult and Child Panoramic programs
- Improved Orthogonality
- TMI section
- Maxillary Sinus exam

ACTEON® Imaging Suite Software, for visualisation and analysis of 2D acquisitions.

Installation and user e-manuals, 2 years standard warranty.

Not included: Workstation PC and monitor, light.





X-MIND® Trium is a class 3R laser product per IEC 60825-1:2007. Avoid direct eye exposure to laser radiation. Viewing the laser output with magnifying optical instruments (for example, surgical microscope and binocular glasses) may pose an eye hazard and thus the user should not direct the beam into an area where such instruments are likely to be used.

EXTRAORAL RADIOLOGY X-MIND® Trium - Upgradeable 3 in 1

X-MIND™ TRIUM PAN (CEPH READY) UPGRADES

X-MIND Trium Pan Ceph Ready Upgrade Ceph right (no sensor included)
X-MIND Trium Pan Ceph Ready Upgrade Ceph left (no sensor included)

W0900007

W0900007

X-MIND™ TRIUM PAN (3D READY)

Control panel on the right W0900225
Control panel on the left W0900226

UPGRADE 3D

X-MIND Trium Pan 3D Ready Upgrade 3D

W0900227

UPGRADES CEPH

X-MIND Trium Pan 3D Ready Upgrade Ceph right (1 sensor Ceph 2D included)
X-MIND Trium Pan 3D Ready Upgrade Ceph left (1 sensor Ceph 2D included)

W0900018 W0900019









^{*}Specific restrictions apply for 11 x 8cm FOV in Ontario, Canada. Contact ACTEON North America for more information.

EXTRAORAL RADIOLOGY X-MIND® Trium - Upgradeable 3 in 1



X-MIND™ TRIUM PAN CEPH

CEPH left (control panel on the right)
CEPH right (control panel on the left)
CEPH right (control panel on the right)
CEPH left (control panel on the left)

W0900003 W0900009 W0900230 W0900232



OrthoPantomograph (OPG), 115-220V, with one CMOS Ceph 2D sensor used for panoramic imaging and cephalometry.

- · Adult and Child Panoramic programs
- Improved Orthogonality
- TMJ section
- · Maxillary Sinus exam
- · Cephalometric exam

ACTEON® Imaging Suite (AIS) Software, for visualisation and analysis of 2D acquisitions. Installation and user e-manuals, 2 years standard warranty.

Not included: Workstation PC and monitor, light.

EXTRAORAL RADIOLOGY X-MIND® Trium - Upgradeable 3 in 1





X-MIND™ TRIUM PAN 3D

Control panel on the right Control panel on the left W0900004 W0900010

OrthoPantomograph (OPG), 115-220V, with one CMOS Ceph 2D sensor used for panoramic imaging and one Flat Panel Detector for CBCT.

- Adult and Child Panoramic programs
- Improved Orthogonality
- TMI section
- Maxillary Sinus exam
- CBCT FOV Ø40x40, Ø60x60, Ø80x80, Ø110x80

Workstation PC with standard monitor. Software preinstalled.

ACTEON® Imaging Suite (AIS) Software, for visualisation and analysis of 2D and 3D acquisitions and treatment planning. Includes five licenses for 3D app, the 3D viewer and implant planning application.

Installation and user e-manuals, 2 years standard warranty.

Not included: light.

UPGRADES

X-MIND™ Trium Pan 3D Upgrade Ceph right (1 2D Ceph sensor included)

W0900018

X-MIND[™] Trium Pan 3D Upgrade Ceph left (1 2D Ceph sensor included)

W0900019

*Specific restrictions apply for 11 x 8cm FOV in Ontario, Canada. Contact ACTEON North America for more information.

L ACTEON IMAG

EXTRAORAL RADIOLOGY X-MIND® Trium - Upgradeable 3 in 1

X-MIND™ TRIUM PAN CEPH 3D



CEPH left (control panel on the right)
CEPH right (control panel on the left)
CEPH right (control panel on the right)
CEPH left (control panel on the left)

W0900006 W0900012 W0900231 W0900233



OrthoPantomograph (OPG), 115-220V, with one CMOS Ceph 2D sensor used for panoramic imaging and cephalometry and one Flat Panel Detector for CBCT.

- Adult and Child Panoramic programs
- · Improved Orthogonality
- TMJ section
- · Maxillary Sinus exam
- Cephalometry exam
- CBCT FOV Ø40x40, Ø60x60, Ø80x80, Ø110x80

Workstation PC with standard monitor. Software preinstalled.

ACTEON® Imaging Suite (AIS) Software, for visualisation and analysis of 2D and 3D acquisitions and treatment planning. Includes five licenses for 3D app, the 3D viewer and implant planning application.

Installation and user e-manuals, 2 years standard warranty.

Not included: light.

*Specific restrictions apply for 11 x 8cm FOV in Ontario, Canada. Contact ACTEON North America for more information.

EXTRAORAL RADIOLOGY AIS & Accessories





ACTEON® IMAGING SUITE

AIS 3D APP STANDARD 1 dongle license	W1000014
AIS 3D APP BASIC 1 electronic license	W1000012
AIS 3D APP BASIC 5 electronic licenses	W1000013
AIS 3D APP STANDARD 5 dongle licenses	W1000015
AIS 3D APP PREMIUM 1 dongle license	W1000016
AIS 3D APP PREMIUM 5 dongle licenses	W1000017
AIS DICOM PACKAGE license (DICOM print license NOT included)	W0900153
DICOM print license	39800751
ACCESSORIES	
Trium POWER CONDITIONER 8A, 120V	61084-03R
Trium 2D ACQUISITION PC AND MONITOR	W0900246
Trium 2D ACQUISITION PC USA	W0911246

EXTRAORAL RADIOLOGY Accessories

Pan 2D sensor for OrthoPantomograph (OPG) (to add a second 2D sensor or for replacement)	W0900036
Ceph 2D sensor for OrthoPantomograph (OPG) and cephalometry (for replacement)	W0900037
3D sensor Plan captor for 3D OrthoPantomograph (OPG) (for replacement)	W0900038
115V Light Reports X-Rays emission during acquisition	29701225.CW
230V Light Reports X-Rays emission during acquisition	29701226.CW
PC workstation (no monitor) Windows 7, 64 bits. Nvidia Graphic card. Not included: network card, monitor, AIS 3D App	W0900040
PC workstation with network card and dongle key (no monitor) Windows 7, 64 bits. Nvidia Graphic card. 1 Gbit Ethernet Network card. 3 dongle. Not included: monitor, AIS 3D App.	W0900042 BD reconstruction
24" PC Monitor	W0900141
Remote trigger & 15m cable	W0900091
Ears protections X-MIND Trium only for arms Ceph	W0900110
XMT CHIN REST, BITE BLOCK, TMJ	W0900080
XMT SLIDING BITE BLOCKS (10pcs)	W0900080
Bite protections X-MIND Trium (500 pcs)	W0900142

X-MIND Trium Self Standing Base

W0900236

EXTRAORAL RADIOLOGY _____X-MIND® Prime



X-MIND™ PRIME 2D PAN

XMP 2D WALL 110-120V **W120002V**

X-MIND® Prime is a smart space-saving solution, providing high-technology capabilities and simplicity of use. This ingenious unit offers all essential diagnostic tools in a cost effective package.

X-MIND® Prime offers a full set of panoramic exams tailored to meet all your clinical applications. It provides all the panoramic exams required for general dentistry: dental panoramic, temporomandibular joints, sinuses.

X-MIND™ PRIME 3D CBCT/PAN

XMP 3D WALL 110-120V + WS + MONITOR + 5 licenses for AIS 3D app **W1200006V**

X-MIND® Prime is a complete imaging solution combining panoramic and 3D X-rays, from general examination to specific treatment planning.

X-MIND® Prime 3D offers a multitude of acquisition programs thanks to its multiple FOV (50x50mm, 85x50mm, and 85x93mm). Once you select the region of interest of the examination, the settings and dose are automatically adjusted.

X-MIND™ PRIME COLUMN

Column for X-MIND® Prime for when mounting to the wall is not possible **W1200101V**

*X-MIND® Prime is currently unavailable in Canada



INTRAORAL CAMERAS Diagnostic tools



Patented SOPRO® innovation using fluorescence and selective chromatic amplification technologies.

SOPROCARE® instantly highlights plaque, gingival inflammation and caries even in their early stages.

- Light (78 g)
- 4 pre-set selections (extraoral, intraoral, CARE, macro)
- Lighting by LED: four white three blue
- Freezing by sensitive touch "SOPROTouch®"



SOPROCARE®

Handpiece S 950 0002

Delivered with four SOPROTIPS, a handpiece holder, SOPRO Imaging software and its license. Needs to be used with a docking station.



Patented SOPRO® innovation using fluorescence technologies.

SOPROLIFE® Detects caries on occlusal and interproximal surfaces, helping to achieve minimally invasive treatments.

- Light (78 g)
- 4 pre-set selections (extraoral, intraoral, LIFE, macro)
- Lighting by LED: 4 whites (Daylight mode) or 4 blues (diagnostic or traitement mode)
- Freezing by sensitive touch "SOPROTouch®"



SOPROLIFE® Handpiece

S 900 0002

Delivered with four SOPROTIPS, a SOPROLIFE® handpiece holder, SOPRO Imaging software and its license. Needs to be used with a docking station.



SOPRO 717 First reveals micro-cracks, micro-lesions, infiltrations all that human vision can not see alone.

- Light (75 g)
- 3 pre-selections (extraoral, intraoral, macro)
- Lighting by LED (x8)
- Freezing by sensitive touch "SOPROTouch®"



SOPRO® 717 FIRST

Handpiece

Delivered with a SOPRO® 717 First handpiece holder, SOPRO Imaging software and its license. Needs to be used with a docking station.

SOPR **617**

- Very light (55 q)
- Freezing by sensitive touch "SOPROTouch®"

Needs to be used with a docking station

- No focus adjustment
- Lighting by LED (x8)
- Non-inverted image
- 105° angle of view



Handpiece Delivered with a SOPRO® 617 handpiece holder, SOPRO Imaging software and its license.

SOPRO® 617 S 617 0002

S 717 0101

INTRAORAL CAMERAS Docking Stations

DOCKING STATIONS



Mini dock USB 2.0 Docking station USB2.0 S 717 1600

Including one 2.5 meters linking cable and one 0.8 meter repeater cable

Dock M-USB-Light Grey

S 717 1036

BUILT-IN DOCKING STATION WITH MEMORY: Features Analog and USB 2.0 Outputs for Dental Unit - LIGHT GREY colored tubing - One built-in docking station with Image Memory, one 21 foot linking cable with LIGHT GREY color tubing, one S-Video Cable, one Composite



Linking cable 8ft

n cable 8ft

With connector on both sides SOPROLIFE®/SOPRO® 717/SOPRO® 617

Linking cable 15ft

922051

922050

With connector on both sides SOPROLIFE®/SOPRO® 717/SOPRO® 617

FOOTSWITCHES



Foot Pedal for M-USB Dock

009002



USB Foot Pedal

009009

SHEATHS



Sheaths for SOPROLIFE®, SOPROCARE® or SOPRO® 717 First 403017

500 sheaths per box

Sheaths for SOPRO® 617

403023

500 sheaths per box

Sheaths for SOPRO® 595 500 sheaths per box

SOPRO HC5

HANDPIECE HOLDERS



Handpiece holder for SOPROLIFE®, SOPROCARE®, SOPRO® 717 First, SOPRO® 617 and SOPRO® 595 Only for unsheathed linking cables.

926063

990135

ACCESSORIES



SOPROTIPS for SOPROLIFE® and SOPROCARE®

Protection tip against light. Sterilization by autoclave.

(Set of 4 pieces)

ULTRASONICS Ultrasonic generators





NEWTRON® BOOSTER US F60216 CANADA F60206

NEWTRON unit delivered with NEWTRON Handpiece, US power Supply

NEWTRON® BOOSTER WITH INTRO TIP KIT US F60211 CANADA F60201

NEWTRON unit delivered with NEWTRON Handpiece, intro kit (#1, 10Z, H3 tips, three color coded torque wrenches, autoclavable plastic box), and US power Supply

NEWTRO



NEWTRON® P5 B.LED US F61116 CANADA F61106

NEWTRON P5 B.LED unit delivered with NEWTRON Slim B.LED Handpiece with white LED ring, US power Supply

NEWTRON® P5 B.LED WITH INTRO KIT US F61117 CANADA F61107

NEWTRON P5 B.LED delivered with NEWTRON Slim B.LED Handpiece with white LED ring, intro kit (#1, 10Z, H3 tips, three color coded torque wrenches, autoclavable plastic box), and US power Supply

Class IIa medical devices (GMED) - CE 0459

ULTRASONICS Ultrasonic generators





NEWTRON® P5 XS B.LED US F62118 CANADA F62106

NEWTRON P5 XS B.LED unit delivered with NEWTRON Slim B.LED Handpiece with white LED ring, US power Supply

NEWTRON® P5 XS B.LED WITH INTRO KIT US F62117 CANADA F62107

NEWTRON P5 XS B.LED delivered with NEWTRON Slim B.LED Handpiece with white LED ring, intro kit (#1, 10Z, H3 tips, three blue color coded torque wrenches, autoclavable plastic box), and US power Supply

Class IIa medical devices (GMED) - CE 0459

ACCESSORIES



NEWTRON® handpiece

F12281

Autoclavable Handpiece: GII 2-Tone (Grey)
NEWTRON OFM



Nose for NEWTRON® handpiece

E12281



NEWTRON® SLIM B.LED handpiece F12900

NEWTRON Slim B.LED Handpiece with blue LED ring; Compatible with NEWTRON B.LED and NEWTRON XS B.LED devices for use with FLAG fluorescent liquid assistance guide

*F.L.A.G. is currently unavailable in Canada



NEWTRON® SLIM B.LED handpiece

F12905

NEWTRON Slim B.LED Handpiece with white LED ring for NEWTRON B.LED, XS B.LED



NEWTRON® SLIM handpiece without LED

F12901

For NEWTRON® Booster sold before 2016



Nose for NEWTRON® SLIM and NEWTRON® SLIM B.LED handpieces

E12904



NEWTRON® SLIM B.LED blue LED ringBlue LED Ring for NEWTRON Slim B.LED handpiece (to illuminate FLAG - not available in Canada)



NEWTRON® SLIM B.LED white LED ring

F12915

White LED Ring for NEWTRON Slim B.LED Handpiece

Light Guide

E12903

LED Light Guide for NEWTRON Slim B.LED Handpiece

ULTRASONICS Ultrasonic Accessories



300ml tank

F62004



500ml tank

F62005



F.L.A.G.™ for B.LED

F62209

FLAG 5, 1.5 ML vials. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED handpiece. Not available in Canada.

F.L.A.G.™ for B.LED

F62207

FLAG 100 ML bottle. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED handpiece. Not available in Canada.

F.L.A.G.™ Package

F62206

FLAG Package 100 ML bottle with Blue LED ring. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED.

*F.L.A.G. is currently unavailable in Canada



Large O-Ring for Handpiece *E15019*

Small O-Ring for Handpiece E15052

O-Ring Kit *E15121*

for P Max XS and P5 NEWTRON XS Irrigation Tanks

Medium O-Ring

E15002

Black Medium O-ring 9 x 1 for NEWTRON Slim and NEWTRON Slim B.LED Handpiece

O-Ring with Tool

F12304

Handpiece 3 gasket O-rings and tool

ACTEON EQUIPMENT

ULTRASONICS Ultrasonic generators

Water Filter, Complete F10386

Replacement water filter cartridge F10389





NEWTRON® LED handpiece F12609 Autoclavable LED Handpiece for P Max XS, P5 LED, P5 XS LED, SP NEWTRON LED and Piezotome 2 / **Implant Center 2**



Light Ring For NEWTRON® LED handpiece F12605

Light Guide E12674 Light guide for the NEWTRON® LED Handpiece

WRENCHES



NEWTRON® wrenchBlue color coded torque wrench

F81322

IMPORTANT

The NEWTRON® wrench is used with surgical endodontics tips: AS3D, AS6D, AS9D, ASLD, ASRD, S12-70D, P14D, P15LD and P15RD. For other tips with yellow ring, the universal wrench (Ref. F00406) is recommended.

Universal wrench	F00406
Blue Plastic Color Coded Wrench	F81302
Orange Plastic Color Coded Wrench	F81303
Black Plastic Color Coded Wrench	F81305

BOXES



Autoclavable plastic box (round) F12351 Clear Cover GII (Empty) Autoclavable Box - Holds four tips total - Requires round torque wrenches to hold the tips (not included)




Autoclavable box F12350

Clear Cover GII (Empty) Autoclavable Box - Holds four tips total - Requires Satelec plastic wrenches to hold the tips (not included)



Autoclavable box F00143 Metal Blue Autoclavable Box - Holds 6 tips (not included)



periodontics H2L H2R F00366 F00367 F00368 Anterior teeth, Premolars & Premolars & diamondmolars, I molars, coated eft-angled right-angled Н3 H4R F00369 F00114 F00115 Anterior teeth, Premolars & Premolars & molars,



molars, left-angled

right-angled

curette





IP1 F02121 **Implant** abutments and wide threads



F02122 threads, left-angled

IP2L



F02123 Medium implant Medium implant threads, right-angled



IP3L F02124 Narrow implant threads, left-angled



F02125 Narrow implant threads, right-angled

periomaintenance



F01001 Short probe



F01004 Long probe



F02162 Premolars & molars. left-angled



TK2-1R F02161 Premolars & molars. right-angled

perioSoft





PH1 (x4) F00702 Anterior teeth



PH2L (x4) F00705 Premolars & molars. left-angled



PH2R (x4) F00706 Premolars & molars. right-angled

perioPrecision



P₂L F00090 Narrow pockets debridement, left-angled



P2R F00091 Narrow pockets debridement. right-angled

endodontics





K10-21 (x4) F43710 Ø 0.1/length 21mm

K10-25 (x4) F43712 Ø 0.1/length 25mm





K15-21 (x4) F43715 Ø 0.15/length 21mm

K15-25 (x4) F43717 Ø 0.15/length 25mm





K25-21 (x4) F43725 Ø 0.25/length 21mm

K25-25 (x4) F43727 Ø 0.25/length 25mm





K30-21 (x4) F43730 Ø 0.3/length 21mm

K30-25 (x4) F43732 Ø 0.3/length 25mm

irrisafe





IRR 20-21 (x4) F43807 Ø 0.2mm length 21mm





IRR 25-21 (x4) F43805

Ø 0.25mm length 21mm

IRR 25-25 (x4) F43806 Ø 0.25mm length 25mm

endosurgery







P15RD F00108 Right-



S12-70D F00118 70° Universal

P14D P15LD F00106 F00107 Universel Left-angled

angled



OSUCCESS

AS6D

F00079

Length

6mm

Apical Surgery


AS3D F00065

Length 3mm

Universal

F00080

Lenath 3mm.




AS9D


9mm

F00081

Length 3mm,










Class IIa medical devices (GMED) - CE 0459







scaling

Scaling kit

F00900

Tips - #1, #2, 10X, 10Z with 4 torque Wrenches and Autoclavable Box



NEWTRON

NEWTRON Intro Kit

F00932

Tips - #1, 10Z, H3, three torque wrenches and Autoclavable Box



perio Precision

Perio Precision kit

F00939

Tips - TK1- 1S, P2L, P2R, 3 Wrenches and Autoclavable Box



periodontics

Perio kit

F00902

Tips - H3, H4L, H4R with 3 torque wrenches and



perio*m*aintenance

Bio-Film Disruption Kit

F00904

Tips - TK1- 1S /TK1- 1L / TK2- 1R / TK2- 1L with 4 Wrenches and Autoclavable Box



implant Protect

Implant Protect kit

F02120

Tips - IP1, IP2L, IP2R, IP3L, IP3R with universal wrench and Autoclavable stand



periodontics

Perio Soft Kit

F00906

Tips - 4 pk PH1, 4 pk PH2L, 4 pk PH2R with 3 Black Wrenches and Autoclavable Box. Only available until supplies last.



excavus

Excavus® kit

F00739

Tips - EX1, EX2, EX3, EXR, EXL with 1 Universal Wrench, and Autoclayable Stand



endosuccess

Endo Success™ Retreatment kit

F00737

ET18D, ETBD, ET20, ET25S, ET25 and ETPR with 1 Universal Wrench and Autoclavable Stand



endosuccess Canal Access Prep

Endo Success™ Canal Access Prep kit

F88180



CAP1, CAP2, CAP3 with universal wrench and Autoclavable stand



endosuccessApical Surgery

Endo Success™ Apical Surgery kit

F00069

AS3D, AS6D, AS9D, ASLD, ASRD Tips with 1 Universal Wrench and Autoclavable Stand

ULTRASONICS Tip kits





Perfect Margin® Rounded kit

F00738

Tips - PM1, PM2, PM3, PM4 with 1 Universal Wrench and Autoclavable Stand





Perfect Margin® Shoulder kit

F00736

Tips - PMS1, PMS2, PMS3, PM4 with 1 Universal Wrench and Autoclavable Stand



perfectmargin Veneers



Perfect Margin® Veneers kit

F02020

Tips - PMV1, PMV2, PMV3, PMV4, PMV5, PMV6 with universal wrench and Autoclavable stand



Retro Kit F00910

Tips - S12-70D, AS6D, ASLD, ASRD with 4 Wrenches and Autoclavable Box



Micro Retro Kit

F00912

Tips - P14D, P15LD, P15RD with 3 Wrenches and Autoclavable Box $\,$



Loosening Kit

F00914

Tips - 5AE with one Orange Wrench and Autoclavable Box



Box Reactor Kit

F00276

5AE tip with Universal Wrench and Reactor Handpiece



AIR-N-GO® easy

The new generation of dual purpose Air Polisher: SUPRA and PERIO

AIR-N-GO® easy with KaVo* adapter	F10130
AIR-N-GO® easy with Sirona* adapter	F10131
AIR-N-GO® easy with W&H* adapter	F10132
AIR-N-GO® easy with Bien-Air* adapter	F10133
AIR-N-GO® easy with Midwest* adapter	F10135

Delivered with a convertible AIR-N-GO® easy handpiece, a SUPRA nozzle, 1 intro kit of 20g bags of CLASSIC powders, 1 AIR-N-GO® easy maintenance kit (replacement o-rings, syringe of silicone grease with brushes, cannula, cleaning drill).

Class IIa medical devices (GMED) - CE 0459





 $^{{\}rm *KaVo},$ Sirona, W&H, Bien-Air and Midwest are registered trademarks.

AIR-N-GO® EASY ACCESSORIES



AIR-N-GO® easy PERIO KIT

F10120

Turns your handpiece into a Perio system to treatperiodontal disease (tooth & implants) on pockets from 3 to 8mm.

Delivered with a AIR-N-GO easy nozzle and a bottle of AIR-N-GO PERIO glycine powder.

Class IIa medical devices (GMED) - CE 0459



AIR-N-GO® SUPRA nozzle

F10125

For all prophylaxis supra-gingival treatments.



AIR-N-GO® PERIO easy nozzle

F10126

Turns your handpiece into a Perio system to treat non surgical periodontal disease (tooth & implants) in pockets from 3 to 8mm.



AIR-N-GO® PERIO nozzle

F10127

Turns your handpiece into a Perio system to treat non surgical and surgical periodontal disease (tooth & implants) in pockets from 8 to 10mm.



AIR-N-GO® PERIO Maintenance nozzle

F10128

Turns your handpiece into a system adapted for periodontal maintenance in pockets until 4mm.



Air-N-Go Easy O-ring kit

F10121

Air-N-Go Easy o-ring kit: includes silicone grease syringe, (3) brushes, o-rings

POWDER - BOTTLES

Class IIa medical devices (GMED) - CE 0459



AIR-N-GO® CLASSIC powder

4 bottles (250g)

Peppermint	F10250
Lemon	F10251
Raspberry	F10252
Cola	F10253
Neutral	F10254



4 flavors AIR-N-GO® CLASSIC powder

4 bottles (250g)

F10255



AIR-N-GO® PERIO powder 3 bottles (160g)

F10270

POWDER - STICKSClass IIa medical devices (GMED) - CE 0459



AIR-N-GO® CLASSIC powder

50 sticks (20g)

Dannarmint	F10256
Peppermint	
Lemon	F10257
Raspberry	F10258
Cola	F10259
Neutral	F10260

AIR-N-GO® ADAPTERS

:0	KaVo* Adapter	F10240
10	Sirona* Adapter	F10241
	W&H* Adapter	F10242
	Bien-Air* Adapter	F10243
(Midwest* Adapter	F10245

*KaVo, Sirona, W&H, Bien-Air and Midwest are registered trademarks



AIR-N-GO® & OTHER ACCESSORIES

	AIR-N-GO® SUPRA tank	F10213
	AIR-N-GO® PERIO tank	F10216
	Ø0.5 mm AIR-N-GO® SUPRA probe	F10207
─	Ø1 mm AIR-N-GO [®] SUPRA probe	F10208
	AIR-N-GO® PERIO probe	F10206
•0	AIR-N-GO® o-rings kit	F10205
	AIR-N-GO® Easy Front Body Autoclavable front body for the AIR-N-GO E	E10207 asy
	Blue ProphyPen Disconnectable polisher handpiece (bayonet) delivered with a 135° titanium no an autoclavable universal wrench.	F10554 szzle and
Homas I	White ProphyPen Disconnectable polisher handpiece (bayonet) delivered with a 135° titanium no an autoclavable universal wrench.	F10557 zzle and
	Titanium nozzle 135° Ø0.7 mm	F10555



Ultrasonic Bone Surgery provides Less Trauma and Less Drama for :

- Atraumatic Extractions
- Bone grafting
- · Sinus elevations
- Ridge Splitting
- Accelerated Orthodontics



PIEZOTOME® CUBE

US F50108 CANADA F50117

Piezotome CUBE with power cord, Piezotome CUBE handpiece, torque wrench, foot pedal, irrigation bracket, and five (5) sterile disposable irrigation lines


Class IIa medical devices (GMED) - CE 0459

NOTE

Some accessories adapted to Piezotome® 2, Piezotome® Solo, Piezotome® CUBE and ImplantCenter™ 2 are not compatible with the first generation Piezotome® and ImplantCenter™, and vice-versa

- References compatible with 1st generation devices are indicated by a blue square :
- References compatible with 2^{nd} and 3^{rd} generation devices are indicated by a grey square \Box

PIEZOTOME® CUBE ACCCESSORIES



IMPLANT CENTER 2



FUNCTIONS

- I-Surge™: implantology
- Piezotome[®]: bone surgery
- NEWTRON®: conventional treatments (perio, endo, ...)

IMPLANT CENTER 2 LED: STARTER KIT

F27168

Delivered with i-Surge LED motor, i-Surge LED motor cord, power cord, a progressive footswitch, 2 brackets, 2 handpiece holders and 2 sterilizable irrigation lines, 30 perforators, an Extraction II pack *Contra Angles not supplied

This reference does not include the PIEZOTOME 2 handpiece or tips.

IMPLANT CENTER 2 LED STARTER KIT WITH LED SURGICAL F27168-CA CONTRA ANGLE (20:1)

Delivered with i-Surge LED motor, i-Surge LED motor cord, power cord, a progressive footswitch, 2 brackets, 2 handpiece holders and 2 sterilizable irrigation lines, 30 perforators, an Extraction II pack and LED Surgical Contra Angle (20:1) - Bien Air

This reference does not include the PIEZOTOME 2 handpiece or tips.

- References compatible with 1st generation devices are indicated by a blue square :
- References compatible with 2nd and 3rd generation devices are indicated by a grey square

ACCESSORIES FOR IMPLANT CENTER™



Implant Center Pedal

Foot Pedal for use with the Implant Center





Implant Center 2 Foot Control

Implant Center 2 Progressive Foot Pedal





ImplantCenter™ Bracket

E27113

Irrigation for use with Implant Center 2, Piezotome Solo, Piezotome 2



Power Cord

E13125

110 Volt Power Cord for Implant Center, I-Surge and Piezotome



ACCESSORIES FOR IMPLANT CENTER™ 2 WITH LED MOTOR



I-Surge™ LED motor

F27210

I-SURGE LED MOTOR for use with the Implant Center 2 LED only *Contra Angles not supplied





i-Surge™ LED Motor Cord

F27211

Cord for use with the i-Surge LED Motor on Implant Center 2 LED only *Contra Angles not supplied





I-Surge™ micro-motor

I-Surge Micro-Motor (Motor Only) *Contra-Angles are not supplied For Implant Center 1 only





I-Surge Cord

F27124

I-Surge Micro-Cord (Cord Only)

*Contra-Angles are not supplied For Implant Center 1 only

SURGERY Accessories

IMPORTANT

The handpieces are equipped with keyway systems to prevent a Piezotome® insert from being screwed onto a NEWTRON® handpiece or a classic insert (scaling, periodontics, endodontics, etc.) from being screwed onto a Piezotome® handpiece.



Piezotome[®] handpiece + cord *F57321*

Non-disconnectable assembly. For use with Piezotome 1 and Implant Center 1 only.



NEWTRON® Handpiece

F12281

Autoclavable Handpiece: GII 2-Tone (Grey) NEWTRON OEM. Compatible with Piezotome 1 and Implant Center 1. Need F57305 to connect.



NEWTRON Cord to Piezotome

F57305

Sterilizable and disconnectable cord for Piezotome® and ImplantCenter™





Piezotome® LED handpiece + cord **F12815** Compatible with Piezotome Solo, Piezotome 2 and Implant Center 2 device



LED Ring F1

Light Ring for the Piezotome 2 LED Handpiece



NEWTRON® LED handpiece

F12609

Autoclavable LED Handpiece for Piezotome 2 / Implant Center 2. Need F12830 to connect.



NEWTRON® LED cord **F12830**Sterilizable and disconnectable cord to connect the NEWTRON LED handpiece to the Implant Center 2

• References compatible with 1st generation devices are indicated by a blue square :

or Piezotome 2

• References compatible with 2^{nd} and 3^{rd} generation devices are indicated by a grey square \blacksquare



Easytorque Wrench *F57350* Easytorque Wrench For Piezotome Tips only



Piezotome® Pack (without tips)Sterilization Box, Piezotome Handpiece, Piezotome Cord and Easytorque Dynamometric Wrench



Piezotome® 2 Pack (without tips) F87523Sterilization Box, Piezotome LED Handpiece,
Piezotome 2 Cord and Easytorque Wrench



Sterilization Box F27156

Piezotome (PZI, PZII) Sterilization box: can hold one handpiece with cord, 2 Tip Kits and 1 Wrench (items are not included)



Irrigation Kit (Sterile Disposable) F5737

Kit includes 10 sterile disposable lines and cartridge.



Irrigation Kit (Autoclavable) F57374

Kit includes 1 autoclavable irrigation line and 10 sterile disposable perforators. (For use with all I-Surge, Piezotome / Piezotome 2, Implant Center / Implant Center 2)



Perforators (Disposable)

10 Disposable perforators for use with the F57374 Irrigation lines



F57372

E57353



Irrigation Clip

Single Irrigation Clip for Surgical Ultrasonic and Micro Motor Devices





Piezotome 2 Foot Control

F57404

Piezotome 2 Progressive Foot Pedal





Piezotome SOLO: Foot Pedal

F57504



E57341

Foot Pedal for use with the Piezotome





Bracket

E27113

Irrigation for use with Implant Center 2, Piezotome Solo, Piezotome 2





Piezotome LED Handpiece

F12815

Piezotome 2, Piezotome Solo LED and Implant Center 2 LED Handpiece and Cord: Nondisconnectable assembly

Piezotome 2 Nose Cone

E12813

Piezotome 2 stainless steel nose cone

LED Light Guide

E12812

12012

Light Guide fo the Piezotome 2 LED Handpiece



Piezotome Bracket

E57325

~ ~

Bottle Holder Bracket for the Piezotome

Hook

E57348

Irrigation Bottle Bracket Hook for Piezotome 1 only



Bone Surgery Kit

Tips - BS1, BS2L, BS2R, BS4, BS5, BS6, Autoclavable Stand for bone grafting procedures

Bone Surgery Kit F87309

Bone Surgery-2 Kit

Tips - BS1-2, BS2L-2, BS2R-2, BS4-2, BS5-2, BS6-2, Autoclavable Stand for bone grafting procedures

Bone Surgery-2 Kit F87509



BS1 - Straight BS2L - Left-angled

BS2R - Right-angled

F87301

F87302 F87303



BS1S-2 Slim - Depth of cut 9mm

BS1L-2 Long - Depth of cut 15mm with universal wrench

F87525

F87527



BS1-RD - Rounded saw, 280° of active surface

F87557



BS₅ BS₆

BS4 BS5

BS₆

F87304 F87305 F87306

F87504 F87505 F87506

- References compatible with 1st generation devices are indicated by a blue square
- References compatible with 2nd and 3rd generation devices are indicated by a grey square



Sinus Lift Kit

Tips - SL1, SL2, SL3, SL4, SL5, Autoclavable Stand (These tips will only work with the Piezotome HP) for lateral approach sinus elevation. Sinus Lift Kit F87319

Sinus Lift-2 Kit

F7519



Sinus Lift tips

SL1	F87311	F87511
SL2	F87312	F87512
SL3	F87313	F87513
SL4	F87314	F87514
SL5	F87315	F87515

- References compatible with 1st generation devices are indicated by a blue square
- References compatible with 2^{nd} and 3^{rd} generation devices are indicated by a grey square



Intralift[™] Kit

Tips - TKW1, TKW2, TKW3, TKW4, TKW5, Autoclavable Stand to perform internal sinus elevation.

Intralift[™] Kit *F87336*

Intralift[™] 2 Kit **F87536**



Intralift[™] tips

TKW1 - Conical tip ø1.35mm (pilot drilling)

TKW2 - Cylindrical tip ø2.1mm (preliminary drilling)

TKW3 - Cylindrical tip ø2.35mm (preliminary drilling)

TKW4 - Cylindrical tip ø2.8mm (secondary drilling)

F87331 F87531

F87332 F87532

F87333 F87533

F87334 F87534



TKW5 - Non cutting "trumpet" tip (sinus membrane elevation)

F87335 F87535



EXTRACTION 2 PACK

Sterilization Box, Piezotome 2 Handpiece, Piezotome 2 Cord, Extraction 2 Kit and Easytorque Wrench for atraumatic extractions

Extraction 2 Pack F87540



Extraction Kit (1st generation) Tips - LC1, LC190, LC2, LC2R, LC2L, Ninja, Autoclavable Stand and Universal Wrench to perform atraumatic extractions

Extraction kit **F87346**

Extraction-2 Kit

Tips - LC1-2, LC1 90-2, LC2-2, LC2R-2, LC2L-2, Ninja-2, Autoclavable Stand and Universal Wrench for atraumatic extractions.

Extraction-2 Kit **F87546**



Extraction tips

LC1 LC2 LC2L LC2R F87307 F87507 F87342 F87542 F87343 F87543 F87344 F87544



Ninja®

F87345

F87545



Crown Extension Kit

BS6, CE1, CE2, CE3 Tips with 1 Universal Wrench and Autoclavable Stand to perform crown lengthening procedures. CE Kit **F87354**

CE 2 Kit **F87554**



Crown Extension tips

CE1 - Diamond-coated spherical tip ø1.75mm

CE2 - Diamond-coated spherical tip ø1.20mm

CE3 - Diamond-coated straight tip on 5mm (laser mark placed at 3mm: biologic width) F87351 F87551

F87352 F87552

F87353 F87553



BS₆

F87306

F87506



Specialist Upgrade Tip Kits

Bone Surgery 2 Tip Kit, Sinus Lift 2 Tip Kit, IntraLift 2 Tip Kit, Clinical Book and DVD F87521



Crest Splitting Kit

CS1, CS2, CS3, CS4, CS5, CS6, tips an autoclavable metal support, an autoclavable universal wrench. CS 2 Kit **F87567**



Crest Splitting tips

CS1	
CS2	
CS3	
CS4	

F87561 F87562 F87563 F87564



CS5 CS6 F87565 F87566



Piezocision™ Kit

Delivered with PZ1, PZ2L, PZ3R, PZ3 tips, an autoclavable metal support and an autoclavable universal wrench.

F87576

F87571

F87572



Piezocision™ tips

PZ1 - Anterior corticotomies.

PZ2L (left oriented) Lateral corticotomies on patient's right vestibular arch

To be used with pendulum motion

PZ2R (right oriented) Lateral corticotomies on patient's left vestibular arch To be used with pendulum motion

PZ3 - Anterior corticotomies for root proximity concerns

F87574

F87573





Essential-2 Kit

BS1S II, BS4 II, LC2 II, SL1 II, SL2 II, SL3 II tips, an autoclavable metal support, an autoclavable universal wrench

F87528



NEWTRON® PERIO PACK

Sterilization Box, NEWTRON Handpiece, NEWTRON Cord, Auto Stand, Tips H1, H2L, H2R, H3, H4L, H4R and Easytorque Wrench NEWTRON® Perio Pack F87320

NEWTRON® Perio 2 LED Pack F87520

NEWTRON® RETRO PACK

Delivered with a sterilization box, an autoclavable NEWTRON® handpiece, an autoclavable NEWTRON® cord, a NEWTRON® Retro kit and a Piezotome® autoclavable wrench.

F87322

NEWTRON PACK (without tips)

Sterilization Box, NEWTRON Handpiece, NEWTRON Cord and Easytorque Wrench

F87324

NEWTRON 2 PACK (without tips)

Sterilization Box, NEWTRON LED Handpiece, NEWTRON LED Cord and Easytorque Wrench

F87524

IMPORTANT

NEWTRON® kits are compatible with the NEWTRON® and NEWTRON® LED handpieces.



NEWTRON® Perio Kit

Tips - H1, H2L, H2R, H3, H4L, H4R, Autoclavable Stand





NEWTRON® Implant Protect Kit

Tips IP1, IP2L, IP2R, IP3L and IP3R, an autoclavable metal support, an autoclavable universal wrench.





NEWTRON® Micro-Retro Kit

Tips - P14D, P15LD, P15RD with 3 Yellow Wrenches and Autoclavable Box F87325



NEWTRON® Retro Kit

Tips - S12-70D, AS6D, ASLD, ASRD, Autoclavable Stand

F87326



Autoclavable Stand - NEWTRON Tips

NEWTRON- Silver Metal Autoclavable Stand - Holds 6 tips (not included) N88449

POLYMERIZATION Curing lights



MINILED SUPERCHARGED

Mini LED™ SuperCharged

F05214

F05217

2,200 mWcm2 Output, (Modes: Full, Pulse & Ramp) Handpiece, Battery, Charger, Power Supply, Light Guide 7.5mm, Oval Light Shield (Cordless)

Mini LED™ SuperCharged handpiece



MINILED ORTHO 2

Mini LED™ Ortho 2

F05220

3,000 mWcm2 Ouput (Modes 4, 8, 16 & 32 seconds) Handpiece, Battery, Charger, Power Supply, Light Guide 5.5mm, Oval Light Shield (Cordless)



Mini LED™ Black

F05400

Mini LED Black handpiece, 7.5mm light guide and light shield. Handpiece requires 3 AA batteries (Batteries are included)

POLYMERIZATION Curing lights





Mini LED™ Standard

F02611



Wavelength: 420 - 480nm

Li-Ion battery

ruur

Handpiece: anodized aluminium

Mini LED™ Handpiece

F02530


Mini LED Handpiece for Battery operated unit

POLYMERIZATION Accessories

For Mini LED™ SuperCharged, Mini LED™ Ortho 2

7	Mini LED Battery Charger Base	F02510
	Mini LED Battery Lithium Ion	F02520
7-1 T#	Power supply with US plug for Mini LED	F02544
	Light Shield Light shield for Mini LED, SuperCharged, Or Black	F05407 tho2 &
	Oval Light Shield Oval Light Shield for Standard Mini LED	F02555

For Mini LED™, Mini LED™ SuperCharged, Mini LED™ Ortho 2



PROSTHESIS & AESTHETICS Gingival Retraction Paste

EXPASYL



EXPASYL™ INTRO KIT

294100

Expasyl Intro Kit (Applicator Gun, 6 Capsules, 12 Applicator Tips) System for temporary gingival retraction.



1 MANUAL APPLICATOR Expasyl Applicaton Gun

260900



EXPASYL™ CAPSULE - 20 PACK (STANDARD) 261030

EXPASYL[™] CAPSULE - 20 PACK (STRAWBERRY) 261001



EXPASYL™ APPLICATOR TIP - BOX OF 40 261040

EXPASYL™ APPLICATOR TIP - BOX OF 100 261005

Manufacturer: Produits Dentaires Pierre Rolland - FRANCE

PROPHYLAXIS & PERIODONTICS Hemostasis and Splinting

HEMOSTASY



Hemostatic dressing with mechanical action

INDICATION

To be used in all cases of moderate bleeding occurring in routine dental practice (scaling, etc.)

HEMOSTASYL™ INTRO KIT - 2 Syringes, 40 Tips

246200

HEMOSTASYL™ SYRINGE - 2 Syringes

246050

HEMOSTASYL™ APPLICATOR TIP - Box of 40 Tips

246100

*HEMOSTASYL™ is currently unavailable in Canada



BIOSPLINT

Splinting ribbons

INDICATIONS

- · Stabilization of mobile teeth
- · Temporary replacement of missing teeth
- · Reinforcement or repairing of temporary or permanent bridges
- · Maintain of inter-dental spaces

PROPERTIES

Contention

- · Low thickness avoids any oral embarrassment
- · May remain in place for many years
- · Treatment with supercritical fluids: improves impregnation of the fibre by composite.

BIOSPLINT SPLINTING RIBBONS - REFILL

BIOSPLINT Ribbons: 8 Splinting Ribbons

201661

*BIOSPLINT is currently unavailable in Canada

Manufacturer: Produits Dentaires Pierre Rolland - FRANCE





Disposable tips for air/water syringes



Box of 250 non sterile tips



White	20166
Blue	20166
Yellow	20166
Green	20166
Pink	20166

Box of 250 non sterile tips



White/Mint 201700



RISKONTROL ADAPTERS

Adapter ADEC II - fits Riskontrol Perfect only	200195
Adapter BELMONT	200158
Adapter Dental EZ	200164
Adapter DCI - fits Riskontrol Perfect only	200152
Adapter Siemens 4000	200178
Adapter Sirona	200154
Adapter KAVO	200171
Adapter Luzzani - fits Riskontrol Perfect only	200150

Manufacturer: Produits Dentaires Pierre Rolland - FRANCE All information for a perfect use of these medical devices are indicated on the website: www.acteongroup.com. Products reserved for professional dental use.



MODULES



SP NEWTRON® 4K5V MODULE with connectors (10 & 4 pin)

F55297



SP NEWTRON® 4K3, 7V MODULE

F55287

Module only. Does not contain harness adapter (available seperately if needed).

PIEZO ULTRASONIC MODULE F55292

SP NEWTRON universal w/ connector, AHP, handpiece tubing, Potentiometer. Note: Does not include Air-Electric switch or power supply. Items are available separately.

PIEZO ULTRASONIC MODULE WITH HYGIENE TIPS F55292S

Reference F55292 with Autoclavable Box, 3 color coded torque wrenches, Tips: #1,10Z, H3. Note: Does not include Air-Electric switch or power supply. Items are available separately.





MODULES



SP NEWTRON® LED MODULE with connector

X55508

PIEZO ULTRASONIC MODULE LED F55521

SP NEWTRON LED (Universal), w/ LED handpiece, HP tubing Potentiometer. Note: Does not include Air-Electric switch or power supply. Items are available separately.

PIEZO ULTRASONIC MODULE LED W/ HYGIENE TIPS F55521S

Reference F55521 with Autoclavable Box, three color coded torque wrenches, Tips: # 1, 10Z, H3. Note: Does not include Air-Electric switch or power supply. Items are available separately.

HANDPIECES

ACTEON INTEGRATED



NEWTRON® handpiece autoclavable, light grey/dark grey

F12281



NEWTRON[®] **LED handpiece** autoclavable


F12609



NEWTRON® SLIM B.LED handpiece with white LED ring

F12905

HANDPIECES ACCESSORIES



ULTRASONICS Accessories

ACCESSORIES





Potentiometer F00524
Potentiometer with wire harness F55213
Potentiometer Complete F55213K

Potentiometer, wire harness, knob & sticker

Potentiometer with wire harness* F55213PC1

- Pelton & Crane

Potentiometer with wire harness* F55213PC2

- Pelton & Crane

*Tubing reference depends on Pelton delivery head series

Knob and Sticker for the SP NEWTRON

Knob and Sticker for the SP NEWTRON

TRANSFORMER 36V 2.78A+CORD DTHPU101-36C13M

Power supply for SP NEWTRON + Mini LED + BUILT-IN SOPRO Dock Needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation

Air Electric Switch E00596-OEM

For automatic operation needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation

Water Adapter M00100

Water Adapter for SP NEWTRON

Water Solenoid Kit M00103

Water Solenoid Kit for installation without water

block

Air Bleed Pack

Air Bleed Pack for SP NEWTRON

ACTEON INTEGR

F00526

M00101

Handpiece Tubing A-dec (61")	F00184
Handpiece Tubing A-dec GREY3 (61")	F01153
Handpiece Tubing A-dec (81")	F00188
Handpiece Tubing A-dec GREY3 (81")	F01154
Handpiece Tubing LED Surf6 (61")	F00616
Handpiece Tubing LED Surf6 (81")	F00629
Handpiece Tubing LED Grey3 (61")	F00632
Handpiece Tubing LED Grey3 (81")	F00633
Handpiece Tubing - Light grey silicone	F00195
Handpiece Tubing - Light grey silicone with connectors	F00197

Handpiece Tubing - Light grey silicone F00197PC1 with connectors - Pelton and Crane. Tubing reference depends on Pelton delivery head series.

Handpiece Tubing - Light grey silicone F00197PC2 with connectors - Pelton and Crane. Tubing reference depends on Pelton delivery head series.

Handpiece Tubing Slim LED G3 61"	F01205
Handpiece Tubing Slim LED G3 81"	F01206

Handpiece Tubing- Light grey silicone F00631



Mini LED™ OEM (Built In)

Mini LED Module, Light Grey Cord, Handpiece, Light Guide 7.5mm, Protective Light Shield (Fast, Pulse & Ramp) *Note - Requires a 24 VAC Transformer.

with light grey cord

F02720

Mini LED™ OEM handpiece

FO2710 ACTEON INTEGRATED





Mini LED™ SuperCharged OEM (Built In)

2,200 mWcm2 Output, (Modes: Full, Pulse & Ramp) Handpiece, Module, Light Grey Cord, Light Guide 7.5mm, Oval Light Shield) Note - Requires a 24 VAC Transformer. with light grey cord F05205

Mini LED™ SuperCharged OEM handpiece

F05218



POLYMERIZATION Modules and Cords

	OEM module	F02730
	TRANSFORMER 36V 2.78A+CORD DTHPU101-36C13M Power supply for SP NEWTRON + Mini LED + BUILT-IN SOPRO Dock Needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation	
	Air Electric Switch For automatic operation needed for retro fitte MINI LED or CAMERA OEM installation	E00596-OEM d NEWTRON,
6	Mini LED Handpiece Tubing OEM - light grey Mini LED Handpiece Tubing OEM - light grey - Pelton & Crane	F02721 F02721PC1
	ADEC ⁻ Surf 6 cord	F02326
	ADEC [*] 500 Gray 3 (medium grey) cord	F02327

*Adec is a registred trademark

INTRAORAL CAMERAS Docking Stations

BUILT-IN DOCKING STATIONS

Dock M-USB-Light Grey S 717 1036

BUILT-IN DOCKING STATION WITH MEMORY: Features Analog and USB 2.0 Outputs for Dental Unit - LIGHT GREY colored tubing - One built-in docking station with Image Memory, one 21 foot linking cable with LIGHT GREY color tubing, one S-Video Cable, one Composite


LINKING CABLES FOR BUILT-IN DOCKING



922067

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental



922070

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Light Grey Linking cable 21ft

922073

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Light Grey Linking cable 21ft

922073PC1

Pelton & Crane FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit - cable needed depends on Pelton delivery head series number.

Light Grey Linking cable 21ft

922073PC2

Pelton & Crane FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit - cable needed depends on Pelton delivery head series number.



ACTEON INTEGRATED J

INTRAORAL CAMERAS Linking cables



White Linking cable 21ft

922067

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Grey Linking cable 21ft

922070

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Light Grey Linking cable 21ft

922073

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Light Grey Linking cable 21ft

922073PC1

Pelton & Crane FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit - linking cable needed depends on Pelton delivery head.

Light Grey Linking cable 21ft

922073PC2

Pelton & Crane FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit - linking cable needed depends on Pelton delivery head.

SURF6 LINKING CABLE

922118

FOR DOCK MU USB w/swivel connector (7m/23ft) w/clip

SURF6 LINKING CABLE

922117

FOR DOCK MU USB w/swivel connector (2m/6ft) w/clip

GREY3 LINKING CABLE

922127

FOR DOCK MU USB w/swivel connector (7m/23ft) w/clip

GREY3 LINKING CABLE

922126

FOR DOCK MU USB w/swivel connector (2m/6ft) w/clip

SURF6 LINKING CABLE

922139

FOR DOCK U USB w/swivel connector (7m/23ft) w/clip

SURF6 LINKING CABLE

22122

FOR DOCK U USB w/swivel connector (2m/6ft) w/clip

GREY3 LINKING CABLE

FOR DOCK U USB w/swivel connector (7m/23ft) w/clip

GREY3 LINKING CABLE

922140

FOR DOCK U USB w/swivel connector (2m/6ft) w/clip

INTRAORAL CAMERAS Linking cables



GREY LINKING CABLE 21ft

922163

FOR BUILT-IN U USB DOCK SOPRO camera linking cable to build into the arm of the dental unit

GREY3 LINKING CABLE

922141

FOR DOCK U USB w/swivel connector (7m/23ft) w/clip

LIGHT GREY LINKING CABLE 21ft

922143

FOR BUILT-IN U USB DOCK

SOPRO camera linking cable to build into the arm of the dental unit

7M MU USB Linking Cable f/MidMark 7M MU USB Linking cable f/MidMark 922133

.... g aase

Light Grey Linking cable 21ft 7M U USB Linking cable f/MidMark 922169

123540

USB CABLE 26 ft for Dock MUSB reference S 717 1001

TRANSFORMER 36V 2.78A+CORD

DTHPU101-36C13M

Power supply for SP NEWTRON + Mini LED + BUILT-IN SOPRO Dock Needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation

8M USB Cable

Air Electric Switch

E00596-OEM

For automatic operation needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation



RELAX. WE'VE GOT YOU COVERED

ACTEON® offers a One-Year, No-Questions-Asked Warranty on all its devices.* If anything goes wrong within the first year, we will fix it for free.**

MINIMALLY
INVASIVE

SOLUTIONS

^{*} Excludes all consumable products including but not limited to: ultrasonic tips, wrenches, imaging plates, prophy powder, air-water syringes, curing light guides.

** ACTEON reserves the right to change this policy without notice or modify accordingly due to

^{**} ACTEON reserves the right to change this policy without notice or modify accordingly due to malicious or deliberate damage.