



Taking simplicity to the next level.


Maxcem Elite knows a thing or two about building success. That's because this remarkable product has been on the market for over 10 years. With exceptional performance and ongoing innovation, it will help you cement your own success — now and long into the future.

Maxcem Elite self-etch, self-adhesive resin cement for indirect restorations is simplicity refined — with higher bond strengths, predictable handling and exceptional ease of use — offering you best in class adhesion when compared to other leading self-adhesive cements.

90% Greater Bond Strength to Dentin.

Maxcem Elite™ employs the same unique adhesive monomer used in Kerr's OptiBond family of bonding agents. The optimized resin matrix and filler systems improve wettability, resulting in greater adhesion to dentin.

Bond Strength of Several Self-Adhesive Cements to Enamel and Dentin



Source: 3rd Party University Study. Data on file.

One-Peel™ Cleanup for consistent outcomes.

Maxcem Elite's thixotropic property facilitates the placement of material and easy cleanup.



Photo: Dr. Sam Halabo, DMD - San Diego, CA

Other Features.


Ease of use. No need to refrigerate, easy cleanup, and no hand mixing required.

Compatible with all indirect restorations. Anterior or posterior ceramics, PFMs and metal restorations, and CAD/CAM materials.

Efficient self-cure mechanism. Sets up completely in the absence of light.

No Color Shift for Better Esthetics!

The color change in Maxcem Elite is imperceptible, when compared to systems that contain tertiary amines. Color shift over time is better than leading self-adhesive or adhesive resin cements.



After aging for 4 weeks in 60°C water. (~1.5 years in 37°C water).

Ordering

34417 Maxcem Elite Intro Kit

Contains: 3 dual syringes (5 g): 2 clear, 1 white; 16 4:1 automix tips (regular); 8 4:1 automix tips (wide); 8 intraoral tips; 8 root-canal tips; Directions for use; Technique guide

Maxcem Elite Individual Shade Refill Packs - 4 Shades

Contains: 2 dual syringes (5 g) choice of one shade; 10 4:1 automix tips (regular); 8 4:1 automix tips (wide); 8 intraoral tips; Directions for use; Technique guide

- 33872 Clear
- 33873 White
- 33874 White Opaque
- 33875 Yellow

34418 Maxcem Elite Bulk Pack

Contains: 4 dual syringes Maxcem Elite (5 g): clear; 16 4:1 automix tips (regular); 8 4:1 automix tips (wide); 8 intraoral tips; 8 root-canal tips; Directions for use

Accessories

- 32862 50 4:1 automix tips (regular)
- 33879 8 4:1 automix tips (wide)
- 33880 8 intraoral tips
- 33881 8 root-canal tips



©2015 Kerr Corporation Lit. No. DMC1888 Rev. 1

All trademarks are property of their respective owners. RelyX™ Unicem 2 is a trademark of 3M. G-CEM LinkAce™ is a trademark of GC America. SpeedCEM manufactured by Ivoclar Vivadent.