



EXPERIENCE THE NEXT GENERATION
COMPOSITE THAT ELEVATES YOUR

ARTISTRY

Introducing


Harmonize™

The next generation composite infused
with Adaptive Response Technology

NEW

LIFELIKE RESTORATIONS NEED YOUR TOUCH— AND **ART**

Harmonize™ begins with **ART**. **Adaptive Response Technology**, a nanoparticle filler network that helps you achieve lifelike restorations with more ease and simplicity than ever. With better blending capabilities and enhanced structural integrity, the **ART** of Harmonize provides your restorations with exceptional strength and unmatched esthetics. It's science and beauty coming together at last.



ESTHETICS

HARMONIZE IS REMARKABLY DIFFICULT TO DETECT


The **ART** filler system in Harmonize diffuses and reflects light in a similar way as human enamel, for an enhanced chameleon effect!

Better blending than leading composites such as Filtek™ Supreme Ultra!

Better and more consistent polishability with ProGloss™ than other materials.

Superior gloss retention compared to leading composites (such as TPH Spectra®, GrandioSO and Tetric EvoCeram®).

Harmonize provided highest gloss value and most consistent polish²



²Study conducted at University of Alabama, 2016.

STRENGTH

EXPERIENCE INCREASED STRENGTH AND DURABILITY


Careful control of the size and shape of the **ART** filler system permits high loading - 81% - for excellent strength.

Reinforced nano-scale filler particle network has improved mechanical properties and is more reactive with resin for more efficient polymerization.

Superior strength to Filtek Supreme Ultra, Tetric EvoCeram, Estelite Sigma Quick and Admira Fusion.

Longer lasting, more likely to resist chipping and fracture over time!

Greater fracture toughness compared to most other composites¹



HANDLING

ADAPTIVE VISCOSITY OFFERS EASIER HANDLING

The **ART** filler network features high loading and a spherical shape to minimize stickiness while also incorporating a rheological modifier. This provides adaptive viscosity: it's softer while sculpting, making it easy to shape the composite. When sculpting is complete, the viscosity increases, which prevents the material from slumping.

Less pullback and fewer cavosurface voids because Harmonize doesn't stick to your composite instruments.¹


Reduced time required for trimming and finishing.²

94%

OF LEADING US DENTISTS RATED HARMONIZE EXCELLENT IN HANDLING³

In a blind evaluation, 94% of leading US dentists stated Harmonize had excellent handling characteristics including stickiness, carvability, spreadability, and viscosity.

An esthetic restoration is a feat of engineering. Harmonize™ provides the optical and handling characteristics to help you build something both beautiful and functional.



ORDER TODAY!
 800.KERR.123
 KerrHarmonize.com



ENAMEL

20 PACK
 .25 Grams per Tip

SHADES	UNIDOSE	SYRINGE
A1E	36440	36410
A2E	36441	36411
A3E	36442	36412
A3.5E	36443	36413
A4E	36444	36414
B1E	36445	36415
C2E	36450	36420
XL	36439	36409

4 Grams Each

SHADES	UNIDOSE	SYRINGE
B2E	36446	36416
B3E	36447	36417
B4E	36448	36418
C1E	36449	36419
C3E	36451	36421
C4E	36452	36422
D2E	36453	36423
D3E	36454	36424
D4E	36455	36425

10 PACK
 .25 Grams per Tip

SHADES	UNIDOSE	SYRINGE
B2E	36446	36416
B3E	36447	36417
B4E	36448	36418
C1E	36449	36419
C3E	36451	36421
C4E	36452	36422
D2E	36453	36423
D3E	36454	36424
D4E	36455	36425

4 Grams Each

DENTIN

10 PACK
 .25 Grams per Tip

SHADES	UNIDOSE	SYRINGE
A1D	36432	36402
A2D	36433	36403
A3D	36434	36404
A3.5D	36435	36405
A4D	36436	36406
B3D	36437	36407
C4D	36438	36408
XL2	36456	36426

4 Grams Each

TRANSLUCENT

10 PACK
 .25 Grams per Tip

SHADES	UNIDOSE	SYRINGE
TA-AMBER	36512	36427
TB-BLUE	36513	36428
TC-CLEAR	36514	36429
TG-GRAY	36515	36430
TSC-SUPER CLEAR		36431

4 Grams Each

Available from Oraltec

Freephone: 0508 478 222
 www.oraltec.co.nz
 orders@oraltec.co.nz



Harmonize™

The next generation composite infused with Adaptive Response Technology

Kerr™

RESTORATIVES