


Mucopren[®] Soft for Relining. Strong Bond Guaranteed!



Mucopren® Soft



Monopren® Soft is a soft relining material based on A silicone for direct and indirect application.



Comfortable processing

- ✓ Can be used chairside and is applied in just a few minutes.
- ✓ Mucopren® Soft can be easily processed with scalpel and bur.

Very comfortable for patients

- ✓ The particularly smooth, hydrophobic silicone surface offers protection against microbial contamination.
- ✓ Permanently elastic.

Durable

- ✓ Outstanding adhesion, does not detach from the prosthesis.
- ✓ High tear resistance, long service life.



Mixing tips, green, 100 tips
 Mixing tips, blue, 60 tips
 Dispensing guns

REF 17235
 REF 17217

Page 52
 Page 52
 Page 53



Mucopren® Soft

Direct reline

Indirect reline

highly recommended recommended



Mucopren® Soft

Base set:

50 mL Mucopren Soft, 50 mL Mucopren silicone sealant, 10 mL Mucopren adhesive, 7 green mixing tips, 20 blue mixing tips, 1 brush holder, 20 single-use brushes, 1 steel bur, accessories

REF 28105

100 mL Normal pack 1:1
 2 x 50 mL Mucopren Soft, 6 green mixing tips

REF 15687

Silicone sealant

50 mL Normal pack 1:1
 1 x 50 mL Mucopren silicone sealant, 10 blue mixing tips

REF 15686

Adhesive

10 mL Normal pack
 10 mL Mucopren® adhesive

REF 14203