



# Vinyl Gloves Progenics™ System

Recent studies show gloves dispensed from Bastion Progenics™ system have **92.5% less bacterial contamination** than traditionally dispensed gloves and with zero waste can **reduce costs by up to 20%**

## Progenics™ System Features

- Examination gloves
- Gloves are dispensed cuff first – reducing the risk of contamination by not touching the finger and palm areas of the glove
- Fast and easy to use
- Cost effective due to reduced waste as gloves are dispensed individually
- Unused gloves cannot be reinserted into the box further reducing cross contamination
- The acrylic wall mounted dispenser bracket keeps gloves away from contaminated surfaces and boxes are easily removed and replaced
- Brackets are supplied in kits with screw fixings, double sided tape and mounting instructions.
- Technical data sheets available
- See page 15 for the Progenics Dispenser Wall Bracket

### Progenics™ Vinyl Clear - Lightly Powdered

Code	Size	Inners Per Carton	Quantity Per Inner	Quantity Per Carton
BPV4582	S	10	100	1000
BPV4583	M	10	100	1000
BPV4584	L	10	100	1000
BPV4585	XL	10	100	1000

### Progenics™ Vinyl Blue - Lightly Powdered

Code	Size	Inners Per Carton	Quantity Per Inner	Quantity Per Carton
BPV4592	S	10	100	1000
BPV4593	M	10	100	1000
BPV4594	L	10	100	1000
BPV4595	XL	10	100	1000

### Progenics™ Vinyl Clear - Powder Free

Code	Size	Inners Per Carton	Quantity Per Inner	Quantity Per Carton
BPV4882	S	10	100	1000
BPV4883	M	10	100	1000
BPV4884	L	10	100	1000
BPV4885	XL	10	100	1000

### Progenics™ Vinyl Blue - Powder Free

Code	Size	Inners Per Carton	Quantity Per Inner	Quantity Per Carton
BPV4892	S	10	100	1000
BPV4893	M	10	100	1000
BPV4894	L	10	100	1000
BPV4895	XL	10	100	1000