



VERIDIA

Certification Pack.

Long Sleeve Apron.

67296

This pack contains
Information Sheet
TGA Certificate

Date Released
September 2020

Certification Pack 67296 | Veridia Australia



Features & Benefits

- Individually folded and wrapped
- Thumb loop style wrists keep sleeves from sliding
- Provides high-level fluid protection
- Apron-style neck to enable gown to be donned and removed quickly.
- Latex free

Applications

- Medical & Hospital
- Dental
- Aged Care & Nursing

Long Sleeve Apron, Water Resistant

67296

Fabric	Polyethylene (PE) Material
Colour	Blue
Weight	42gms
Sizing	One size fits all
Style	Apron style Over the head neck Waist ties
TGA Approved	ARTG ID: 335692
Stock Code	67296



Australian Government
Department of Health
Therapeutic Goods Administration

Public Summary

Summary for ARTG Entry:	335692	Reliance Healthcare Supplies - Gown, isolation, single use
ARTG entry for	Medical Device Included Class 1	
Sponsor	Reliance Healthcare Supplies	
Postal Address	33 Wright Road, Keilor Park, VIC, 3043 Australia	
ARTG Start Date	4/05/2020	
Product Category	Medical Device Class 1	
Status	Active	
Approval Area	Medical Devices	

Conditions

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Manufacturers

Name	Address
Reliance Healthcare Supplies	33 Wright Road Keilor Park, VIC, 3043 Australia

Products

1 . Gown, isolation, single use

Product Type	Single Device Product	Effective Date	4/05/2020
GMDN	35492 Gown, isolation, single use		
Intended Purpose	A multipurpose garment, made of natural, synthetic or a combination of materials worn by healthcare providers or visitors and patients to isolate them from the patient. The isolation gown protects the patient from the transfer of infectious agents carried by the healthcare provider or visitor; or it can protect the healthcare provider or visitor from a contagious agent which has infected the patient. Depending on the device's intended purpose, it can be fluid resistant or impervious to fluid. It can be disposable or reusable.		

Specific Conditions

No Specific Conditions included on Record

© Commonwealth of Australia. This work is copyright. You are not permitted to re-transmit, distribute or commercialise the material without obtaining prior written approval from the Commonwealth. Further details can be found at <http://www.tga.gov.au/about/website-copyright.htm>.

Public Summary