



Cavex laboratory materials are the natural extension of the dental practice and must therefore meet stringent quality requirements. Thanks to Cavex's close collaboration with dental technicians and prosthetists these materials are also closely matched to the specific requirements of the dental laboratory.

For dental use only.

Cavex LC Dental Tray

The main purpose of the Cavex LC Dental Tray, light curing tray material is the fabrication of individual and functional impression trays. The material is also indicated for producing base plates for wax-ups and for bite templates.

The sheets are anatomically shaped and vestibular enlarged for better use on moulds of dentate jaws. The design and the handling characteristics make the sheet suitable for the upper as well as the lower jaw.

The material can be cured with halogen and UV. Once cured properly the material is dimensional stable.

INDICATIONS

- L individual impression trays
- L functional impression trays
- L base plates for wax-ups
- L bite templates



FEATURES / BENEFITS

- L vestibular enlarged / for better use with dentate jaws
- L high molecular acrylate with inorganic filler / easy manipulating before - dimensional stable after curing
- L only one sheet model for upper and lower / less stock needed
- L pink opaque colour
- L curing in 3 - 5 minutes with halogen or UV
- L light curing / cures on demand
- L 3 years shelf life
- L 50 sheets per package

ORDERING INFORMATION

Article code	Article description	Contents
BP010	Cavex LC Dental Tray	50 pcs

Cavex Holland BV
P.O. Box 852
2003 RW Haarlem

T+31 23 530 77 00
F+31 23 535 64 82
dental@cavex.nl
www.cavex.nl