



PALAVENEER

TOOTH-COLOURED DENTURE ACRYLICS

PalaVeneer Dentine

Tooth coloured acrylic for deposition of PalaVeneers veneer shells

PalaVeneer Dentine is a self-curing, tooth colored PMMA prosthetic acrylic, particularly suited to fix PalaVeneer veneer shells.

PalaVeneer Dentine distinguishes itself through excellent processing properties.

Together with PalaVeneer veneer shells, the technician obtains a perfect, Premium-identical layering in the anterior region.

PalaVeneer Dentine was modified with the proven NanoPearls® material for maximum abrasion stability and fracture-resistance.

Benefits

- Optimal processing properties
- Good adhesive bond
- Simple processing
- Rapid swelling
- Colour stable
- Fluorescent material
- Adapted to the mechanical properties of PalaVeneer 6 and PalaVeneer 8 material
- Good polishability

Indication

- Fixing and backing PalaVeneer veneer shells on metal frameworks
- Fixing and backing milled acrylic teeth on metal frameworks
- To restore/bond facets
- To produce long term temporaries (facet technique)
- To aesthetically finish pre-fabricated teeth

System components

- PalaVeneer 16 shades A1 – D4
- Palabond
- Signum opaque 16 shades A1 – D4
- Signum metal bond
- Pala Lab Putty



PalaVeneer Dentine refills

Delivery form	Art. Code
Liquid 25ml DE/GB	6605 7636
Liquid 25ml IT/ES	6607 0747
Liquid 25ml FR/NL	6607 0749
Liquid 80ml DE/GB	6605 7637
Liquid 80ml IT/ES	6607 0748
Liquid 80ml FR/NL	6607 0750
Powder 35g	
A1	6605 7540
A2	6605 7565
A3	6605 7566
A3.5	6605 7567
A4	6605 7568
B1	6605 7569
B2	6605 7570
B3	6605 7621
B4	6605 7622
C1	6605 7623
C2	6605 7624
C3	6605 7625
C4	6605 7626
D2	6605 7627
D3	6605 7628
D4	6605 7629
I	6605 7632

For more Information or to place orders call:

PH: 0508 478 222

orders@oraltec.co.nz

www.oraltec.co.nz

