

BIOPLAST® XTREME Kit*

Two different degrees of hardness in one pressure moulding blank thanks to the integrated reinforcement in the delicate anterior region. Absorbs hard and flat impacts in contact sports. V-shaped harder insert for variable size of front tooth protection.

Kits come with 1 BIOPLAST® XTREME 5.0 x 125 mm, 1 bag of 15 g CETRON® cleansing powder, 1 BIOPLAST® insulating agent, retainer box.

- BIOPLAST® XTREME Kit, red #3283
- BIOPLAST® XTREME Kit, blue #3284
- BIOPLAST® XTREME Kit, black #3285

* only for MINISTAR® from 4/92 and BIOSTAR® from 4/99



BIOPLAST® XTREME PRO Kit*

For professional sports where hard and flat impacts often occur. Additional hard-elastic layer for integrated reinforcement in the anterior region.

Kits come with 1 BIOPLAST® XTREME 4.0 x 125 mm, 1 BIOPLAST® clear 2.0 x 125 mm, 1 DURAN® 0.75 x 125 mm, 1 bag of 15 g CETRON® cleansing powder, 1 BIOPLAST® insulating agent, retainer box.

- BIOPLAST® XTREME Kit, red #3292
- BIOPLAST® XTREME Kit, blue #3293
- BIOPLAST® XTREME Kit, black #3294

* only for MINISTAR® from 4/92 and BIOSTAR® from 4/99



BIOPLAST® XTREME*

BIOPLAST® XTREME material is ideally suited for simple and safe preparation of individual mouth guard. Hard V-shaped reinforcement layer ensures additional protection in the anterior region, allowing precise positioning in the desired extension. Soft layer can be used for lamination with BIOPLAST® or for insertion of the counter bite.

- 4.0 x 125 mm red #3295
- 4.0 x 125 mm blue #3296
- 4.0 x 125 mm black #3297
- 4.0 x 125 mm clear #3298
- 5.0 x 125 mm red #3286
- 5.0 x 125 mm blue #3287
- 5.0 x 125 mm black #3288
- 5.0 x 125 mm clear #3289

10 pieces * only for MINISTAR® from 4/92 and BIOSTAR® from 4/99





BIOPLAST® XTREME Deco*

The proven XTREME material for the fabrication of individual mouth guard is now available in many individual designs. A new manufacturing technology enables photo-realistic, multi-colour designing. Material does **not** bond to acrylic.

○ 5.0 × 125 mm	<i>chinese</i>	#3312
○ 5.0 × 125 mm	<i>dragon</i>	#3313
○ 5.0 × 125 mm	<i>kiss</i>	#3314
○ 5.0 × 125 mm	<i>aqua</i>	#3315
○ 5.0 × 125 mm	<i>stone</i>	#3316
○ 5.0 × 125 mm	<i>lion</i>	#3317

5 pieces

CE

* only for MINISTAR® from 4/92 and BIOSTAR® from 4/99

BIOPLAST® Color

Soft, pliable rubber like material, does **not** bond to acrylic.

Application

Coloured mouth guards, positioners

● 3.0 × 125 mm	<i>blue-transparent</i>	#3356
● 3.0 × 125 mm	<i>red-transparent</i>	#3355
○ 3.0 × 125 mm	<i>white</i>	#3445
● 3.0 × 125 mm	<i>black</i>	#3448
● 3.0 × 125 mm	<i>neon red</i>	#3352
● 3.0 × 125 mm	<i>neon green</i>	#3353
● 3.0 × 125 mm	<i>azure</i>	#3354
○ 3.0 × 125 mm	<i>Color set</i>	#3271

10 pieces

CE



BIOPLAST® Multicolor

Soft, pliable rubber like material, does **not** bond to acrylic.

Application

Coloured mouth guards, positioners

● 3.0 × 125 mm	<i>Germany</i>	#3390
● 3.0 × 125 mm	<i>France</i>	#3391
● 3.0 × 125 mm	<i>Italy</i>	#3392
● 3.0 × 125 mm	<i>Spain</i>	#3393
● 3.0 × 125 mm	<i>Sweden</i>	#3394
● 3.0 × 125 mm	<i>Brazil</i>	#3395
● 3.0 × 125 mm	<i>Transylvania</i>	#3396
● 3.0 × 125 mm	<i>Atlantis</i>	#3397
● 3.0 × 125 mm	<i>White dots</i>	#3398
● 3.0 × 125 mm	<i>Green dots</i>	#3399
○ 3.0 × 125 mm	<i>Multicolor-Set</i>	#3400

10 pieces

CE





BIOPLAST®

Soft pliable rubber like material, does **not** bond to acrylic, clear-transparent.

Application

Duplicating moulds, mouth guard, bruxism splints, positioner, bleaching splints, OSAMU-Retainer®

- 1.0 × 125 mm #3270
- 1.5 × 125 mm #3183
- 2.0 × 125 mm #3185
- 3.0 × 125 mm #3187
- 4.0 × 125 mm** #3188
- 5.0 × 125 mm*/** #3374

10/100 pieces * 10 pieces
**only for MINISTAR® from 4/92 and BIOSTAR® from 4/99

CE



BIOPLAST® bleach

Soft pliable rubber like material, best suited for fabrication of bleaching splints. Material does **not** bond to acrylic .

- 1.0 × 125 mm

#3190

10/100 pieces

CE



NEW: BRUX CHECKER® white



BRUX CHECKER® acc. to Prof. Sato

Foil with one-side colour coating for analysis and therapy planning as well as for control of prosthetic restorations.

- 0.1 × 125 mm, red
- 0.1 × 125 mm, white

#3209

#3210

10 pieces

CE

Manual "Diagnosis of occlusion in bruxism" (Kanagawa Dental College, Japan) can be found at www.scheu-dental.com/en/downloads/other